

Surrogate Summary

SDG No.: Q1984

Client: Nobis Group

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030124.D	PIBLK-PS030124.D	2,4-DCAA	1	500	419	84		32	138
		2,4-DCAA	2	500	418	84		32	138
I.BLK-PS030235.D	PIBLK-PS030235.D	2,4-DCAA	1	500	451	90		32	138
		2,4-DCAA	2	500	462	92		32	138
Q1982-08MS	TP-9MS	2,4-DCAA	1	500	493	99		27	122
		2,4-DCAA	2	500	458	92		27	122
Q1982-08MSD	TP-9MSD	2,4-DCAA	1	500	529	106		27	122
		2,4-DCAA	2	500	499	100		27	122
Q1984-01	OU4-PCS-TC-33-050725	2,4-DCAA	1	500	383	77		27	122
		2,4-DCAA	2	500	347	69		27	122
Q1984-03	OU4-PCS-TC-34-050725	2,4-DCAA	1	500	373	75		27	122
		2,4-DCAA	2	500	334	67		27	122
I.BLK-PS030247.D	PIBLK-PS030247.D	2,4-DCAA	1	500	519	104		32	138
		2,4-DCAA	2	500	531	106		32	138
Q1984-05	OU4-PCS-TC-35-050725	2,4-DCAA	1	500	323	65		27	122
		2,4-DCAA	2	500	291	58		27	122
Q1984-07	OU4-TS-24-050725	2,4-DCAA	1	500	196	39		27	122
		2,4-DCAA	2	500	183	37		27	122
Q1984-09	OU4-TS-25-050725	2,4-DCAA	1	500	220	44		27	122
		2,4-DCAA	2	500	207	41		27	122
Q1984-11	OU4-TS-26-050725	2,4-DCAA	1	500	249	50		27	122
		2,4-DCAA	2	500	228	46		27	122
Q1984-13	OU4-TS-27-050725	2,4-DCAA	1	500	291	58		27	122
		2,4-DCAA	2	500	287	57		27	122
Q1984-15	OU4-TS-28-050725	2,4-DCAA	1	500	339	68		27	122
		2,4-DCAA	2	500	326	65		27	122
I.BLK-PS030255.D	PIBLK-PS030255.D	2,4-DCAA	1	500	508	102		32	138
		2,4-DCAA	2	500	527	105		32	138
I.BLK-PS030275.D	PIBLK-PS030275.D	2,4-DCAA	1	500	471	94		32	138
		2,4-DCAA	2	500	469	94		32	138
PB167996BL	PB167996BL	2,4-DCAA	1	500	473	95		27	122
		2,4-DCAA	2	500	450	90		27	122
PB167996BS	PB167996BS	2,4-DCAA	1	500	554	111		27	122
		2,4-DCAA	2	500	539	108		27	122
I.BLK-PS030287.D	PIBLK-PS030287.D	2,4-DCAA	1	500	463	93		32	138
		2,4-DCAA	2	500	462	92		32	138