

Surrogate Summary

SDG No.: Q2001

Client: ENTACT

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030124.D	PIBLK-PS030124.D	2,4-DCAA	1	500	419	84		70 (39)	130 (175)
		2,4-DCAA	2	500	418	84		70 (39)	130 (175)
I.BLK-PS030189.D	PIBLK-PS030189.D	2,4-DCAA	1	500	482	96		70 (39)	130 (175)
		2,4-DCAA	2	500	491	98		70 (39)	130 (175)
PB167988BL	PB167988BL	2,4-DCAA	1	500	480	96		70 (39)	130 (175)
		2,4-DCAA	2	500	460	92		70 (39)	130 (175)
PB167963TB	PB167963TB	2,4-DCAA	1	500	831	166	*	70 (39)	130 (175)
		2,4-DCAA	2	500	497	99		70 (39)	130 (175)
Q2001-03	WC-A4-03-C	2,4-DCAA	1	500	483	97		70 (39)	130 (175)
		2,4-DCAA	2	500	481	96		70 (39)	130 (175)
I.BLK-PS030201.D	PIBLK-PS030201.D	2,4-DCAA	1	500	485	97		70 (39)	130 (175)
		2,4-DCAA	2	500	493	99		70 (39)	130 (175)
Q2001-03MS	WC-A4-03-CMS	2,4-DCAA	1	500	383	77		70 (39)	130 (175)
		2,4-DCAA	2	500	388	78		70 (39)	130 (175)
Q2001-03MSD	WC-A4-03-CMSD	2,4-DCAA	1	500	386	77		70 (39)	130 (175)
		2,4-DCAA	2	500	390	78		70 (39)	130 (175)
Q2001-07	WC-A1-05-C	2,4-DCAA	1	500	466	93		70 (39)	130 (175)
		2,4-DCAA	2	500	458	92		70 (39)	130 (175)
Q2001-11	WC-A1-06-C	2,4-DCAA	1	500	477	95		70 (39)	130 (175)
		2,4-DCAA	2	500	477	95		70 (39)	130 (175)
Q2001-15	WC-A1-07-C	2,4-DCAA	1	500	461	92		70 (39)	130 (175)
		2,4-DCAA	2	500	482	96		70 (39)	130 (175)
I.BLK-PS030213.D	PIBLK-PS030213.D	2,4-DCAA	1	500	487	97		70 (39)	130 (175)
		2,4-DCAA	2	500	495	99		70 (39)	130 (175)
I.BLK-PS030223.D	PIBLK-PS030223.D	2,4-DCAA	1	500	491	98		70 (39)	130 (175)
		2,4-DCAA	2	500	503	101		70 (39)	130 (175)
PB167988BS	PB167988BS	2,4-DCAA	1	500	578	116		70 (39)	130 (175)
		2,4-DCAA	2	500	579	116		70 (39)	130 (175)
I.BLK-PS030235.D	PIBLK-PS030235.D	2,4-DCAA	1	500	451	90		70 (39)	130 (175)
		2,4-DCAA	2	500	462	92		70 (39)	130 (175)