

Surrogate Summary

SDG No.: Q2010

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030124.D	PIBLK-PS030124.D	2,4-DCAA	1	500	419	84		39	175
		2,4-DCAA	2	500	418	84		39	175
I.BLK-PS030179.D	PIBLK-PS030179.D	2,4-DCAA	1	500	469	94		39	175
		2,4-DCAA	2	500	473	95		39	175
Q2007-07MS	OR-636-COMP-11MS	2,4-DCAA	1	500	328	66		10	141
		2,4-DCAA	2	500	294	59		10	141
Q2007-07MSD	OR-636-COMP-11MSD	2,4-DCAA	1	500	340	68		10	141
		2,4-DCAA	2	500	309	62		10	141
I.BLK-PS030189.D	PIBLK-PS030189.D	2,4-DCAA	1	500	482	96		39	175
		2,4-DCAA	2	500	491	98		39	175
I.BLK-PS030213.D	PIBLK-PS030213.D	2,4-DCAA	1	500	487	97		39	175
		2,4-DCAA	2	500	495	99		39	175
PB167976BL	PB167976BL	2,4-DCAA	1	500	514	103		10	141
		2,4-DCAA	2	500	492	98		10	141
PB167976BS	PB167976BS	2,4-DCAA	1	500	532	106		10	141
		2,4-DCAA	2	500	528	106		10	141
I.BLK-PS030220.D	PIBLK-PS030220.D	2,4-DCAA	1	500	493	99		39	175
		2,4-DCAA	2	500	502	100		39	175
I.BLK-PS030235.D	PIBLK-PS030235.D	2,4-DCAA	1	500	451	90		39	175
		2,4-DCAA	2	500	462	92		39	175
Q2010-01	TP-10	2,4-DCAA	1	500	57.8	12		10	141
		2,4-DCAA	2	500	55.1	11		10	141
Q2010-02	TP-13	2,4-DCAA	1	500	306	61		10	141
		2,4-DCAA	2	500	319	64		10	141
Q2010-03	TP-14	2,4-DCAA	1	500	355	71		10	141
		2,4-DCAA	2	500	336	67		10	141
Q2010-04	TP-17	2,4-DCAA	1	500	360	72		10	141
		2,4-DCAA	2	500	343	69		10	141
I.BLK-PS030247.D	PIBLK-PS030247.D	2,4-DCAA	1	500	519	104		39	175
		2,4-DCAA	2	500	531	106		39	175