

Surrogate Summary

SDG No.: Q2033

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PP071388.D	PIBLK-PP071388.D	Tetrachloro-m-xylene	1	20	17.1	85		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	17.6	88		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	17.2	86		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	17.6	88		70 (60)	130 (140)
I.BLK-PP072080.D	PIBLK-PP072080.D	Tetrachloro-m-xylene	1	20	16.7	83		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	18.4	92		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	19.4	97		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	20.9	104		70 (60)	130 (140)
PB168005BL	PB168005BL	Tetrachloro-m-xylene	1	20	18.9	94		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	19.6	98		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	21.0	105		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	21.9	109		30 (10)	150 (173)
PB168005BS	PB168005BS	Tetrachloro-m-xylene	1	20	19.4	97		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	20.7	103		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	21.2	106		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	23.0	115		30 (10)	150 (173)
PB168005BSD	PB168005BSD	Tetrachloro-m-xylene	1	20	19.2	96		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	20.5	103		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	21.0	105		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	23.4	117		30 (10)	150 (173)
Q2033-01	TW-WTS-08	Tetrachloro-m-xylene	1	20	18.2	91		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	13.3	66		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	21.8	109		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	15.4	77		30 (10)	150 (173)
I.BLK-PP072095.D	PIBLK-PP072095.D	Tetrachloro-m-xylene	1	20	16.9	85		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	17.9	90		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	19.2	96		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	21.0	105		70 (60)	130 (140)