

Instrument ID: BNA_P

Daily Analysis Runlog For Sequence/QC Batch ID # BP052325

Review By	Rahul	Review On	5/27/2025 10:49:36 AM
Supervise By	Jagrut	Supervise On	5/27/2025 2:39:24 PM
SubDirectory	BP052325	HP Acquire Method	BNA_P
		HP Processing Method	BP051325

STD. NAME	STD REF.#
Tune/Reschk	SP6757
Initial Calibration Stds	SP6784,SP6785,SP6786,SP6787,SP6788,SP6790,SP6789,SP6791
CCC	SP6787
Internal Standard/PEM	S12666,10ul/1000ul sample
ICV/I.BLK	SP6770
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	DFTPP	DFTPP	BP024787.D	23 May 2025 10:59		RC/JU	Ok
2	SSTDCCC040	SSTDCCC040	BP024788.D	23 May 2025 11:40		RC/JU	Ok,M
3	PB168131BL	PB168131BL	BP024789.D	23 May 2025 12:20		RC/JU	Ok
4	PB168131BS	PB168131BS	BP024790.D	23 May 2025 13:01		RC/JU	Ok,M
5	Q2093-14DL	WATER-COMPDL	BP024791.D	23 May 2025 13:45		RC/JU	Ok
6	Q2097-18	RT-3888	BP024792.D	23 May 2025 14:26		RC/JU	Ok
7	Q2097-12	RT3419	BP024793.D	23 May 2025 15:06		RC/JU	Ok
8	Q2097-16	72-12013	BP024794.D	23 May 2025 15:47		RC/JU	Ok
9	Q2094-01	COMP-1	BP024795.D	23 May 2025 16:28		RC/JU	Ok
10	Q2094-02	COMP-2	BP024796.D	23 May 2025 17:09		RC/JU	Ok
11	Q2094-03	COMP-3	BP024797.D	23 May 2025 17:49		RC/JU	Ok
12	Q2094-04	COMP-4	BP024798.D	23 May 2025 18:30		RC/JU	Ok
13	Q2094-05	COMP-5	BP024799.D	23 May 2025 19:11		RC/JU	Ok
14	Q2073-01	GDW1	BP024800.D	23 May 2025 19:51		RC/JU	Ok
15	Q2073-02	GDW2	BP024801.D	23 May 2025 20:32		RC/JU	Ok,M
16	Q2101-04	TP-1-MHE	BP024802.D	23 May 2025 21:13		RC/JU	Ok
17	Q2095-08	WCS-TP2	BP024803.D	23 May 2025 21:54		RC/JU	Ok
18	Q2106-01	TW-WTS-09	BP024804.D	23 May 2025 22:34		RC/JU	Ok

Instrument ID: BNA_P

Daily Analysis Runlog For Sequence/QC Batch ID # BP052325

Review By	Rahul	Review On	5/27/2025 10:49:36 AM		
Supervise By	Jagrut	Supervise On	5/27/2025 2:39:24 PM		
SubDirectory	BP052325	HP Acquire Method	BNA_P	HP Processing Method	BP051325
STD. NAME	STD REF.#				
Tune/Reschk	SP6757				
Initial Calibration Stds	SP6784,SP6785,SP6786,SP6787,SP6788,SP6790,SP6789,SP6791				
CCC	SP6787				
Internal Standard/PEM	S12666,10ul/1000ul sample				
ICV/I.BLK	SP6770				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

M : Manual Integration