

Surrogate Summary

SDG No.: Q2074

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030124.D	PIBLK-PS030124.D	2,4-DCAA	1	500	419	84		39	175
		2,4-DCAA	2	500	418	84		39	175
I.BLK-PS030313.D	PIBLK-PS030313.D	2,4-DCAA	1	500	492	98		39	175
		2,4-DCAA	2	500	474	95		39	175
PB168097BL	PB168097BL	2,4-DCAA	1	500	417	83		10	141
		2,4-DCAA	2	500	357	71		10	141
PB168097BS	PB168097BS	2,4-DCAA	1	500	535	107		10	141
		2,4-DCAA	2	500	506	101		10	141
I.BLK-PS030323.D	PIBLK-PS030323.D	2,4-DCAA	1	500	497	99		39	175
		2,4-DCAA	2	500	492	98		39	175
Q2071-01MS	L1-WC-7MS	2,4-DCAA	1	500	517	103		10	141
		2,4-DCAA	2	500	400	80		10	141
Q2071-01MSD	L1-WC-7MSD	2,4-DCAA	1	500	518	104		10	141
		2,4-DCAA	2	500	403	81		10	141
I.BLK-PS030333.D	PIBLK-PS030333.D	2,4-DCAA	1	500	501	100		39	175
		2,4-DCAA	2	500	496	99		39	175
I.BLK-PS030358.D	PIBLK-PS030358.D	2,4-DCAA	1	500	513	103		39	175
		2,4-DCAA	2	500	501	100		39	175
Q2074-01	TP-12	2,4-DCAA	1	500	317	63		10	141
		2,4-DCAA	2	500	290	58		10	141
Q2074-02	TP-7	2,4-DCAA	1	500	323	65		10	141
		2,4-DCAA	2	500	271	54		10	141
Q2074-03	TP-15	2,4-DCAA	1	500	335	67		10	141
		2,4-DCAA	2	500	317	63		10	141
Q2074-04	TP-20	2,4-DCAA	1	500	345	69		10	141
		2,4-DCAA	2	500	305	61		10	141
Q2074-05	TP-38	2,4-DCAA	1	500	374	75		10	141
		2,4-DCAA	2	500	343	69		10	141
Q2074-06	TP-19	2,4-DCAA	1	500	352	70		10	141
		2,4-DCAA	2	500	310	62		10	141
Q2074-07	TP-40	2,4-DCAA	1	500	280	56		10	141
		2,4-DCAA	2	500	239	48		10	141
Q2074-08	TP-18	2,4-DCAA	1	500	324	65		10	141
		2,4-DCAA	2	500	297	59		10	141
I.BLK-PS030368.D	PIBLK-PS030368.D	2,4-DCAA	1	500	514	103		39	175
		2,4-DCAA	2	500	513	103		39	175