

Surrogate Summary

SDG No.: Q2078

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PP072163.D	PIBLK-PP072163.D	Tetrachloro-m-xylene	1	20	16.3	82		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	16.7	84		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	16.6	83		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	17.0	85		70 (60)	130 (140)
I.BLK-PP072243.D	PIBLK-PP072243.D	Tetrachloro-m-xylene	1	20	17.6	88		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	18.2	91		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	18.4	92		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	20.7	104		70 (60)	130 (140)
PB168099BL	PB168099BL	Tetrachloro-m-xylene	1	20	19.4	97		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	20.1	101		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	20.3	101		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	21.2	106		30 (32)	150 (175)
PB168099BS	PB168099BS	Tetrachloro-m-xylene	1	20	19.7	99		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	19.9	99		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	19.1	96		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	21.7	109		30 (32)	150 (175)
Q2078-02	WC-A4-04-C	Tetrachloro-m-xylene	1	20	16.6	83		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	12.0	60		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	16.3	82		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	14.1	70		30 (32)	150 (175)
Q2078-06	WC-A4-05-C	Tetrachloro-m-xylene	1	20	16.6	83		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	10.9	54		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	17.6	88		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	14.3	71		30 (32)	150 (175)
Q2078-10	WC-A1-06A-C	Tetrachloro-m-xylene	1	20	14.5	73		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	12.1	61		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	20.2	101		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	14.3	72		30 (32)	150 (175)
Q2078-14	WC-A1-07A-C	Tetrachloro-m-xylene	1	20	13.2	66		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	10.7	54		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	16.2	81		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	13.1	65		30 (32)	150 (175)
Q2078-18	WC-A4-06-C	Tetrachloro-m-xylene	1	20	12.3	61		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	11.3	57		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	16.2	81		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	13.7	68		30 (32)	150 (175)
I.BLK-PP072258.D	PIBLK-PP072258.D	Tetrachloro-m-xylene	1	20	17.8	89		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	18.5	92		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	19.0	95		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	20.9	105		70 (60)	130 (140)
Q2084-01MS	OR-640-COMP-66MS	Tetrachloro-m-xylene	1	20	18.5	92		30 (32)	150 (144)

() = LABORATORY INHOUSE LIMIT

Surrogate Summary

SDG No.: Q2078

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
Q2084-01MS	OR-640-COMP-66MS	Decachlorobiphenyl	1	20	18.2	91		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	19.8	99		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	20.1	100		30 (32)	150 (175)
Q2084-01MSD	OR-640-COMP-66MSD	Tetrachloro-m-xylene	1	20	18.6	93		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	19.3	96		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	21.2	106		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	21.2	106		30 (32)	150 (175)
I.BLK-PP072273.D	PIBLK-PP072273.D	Tetrachloro-m-xylene	1	20	17.8	89		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	18.2	91		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	18.8	94		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	20.2	101		70 (60)	130 (140)