

Surrogate Summary

SDG No.: Q2100

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030124.D	PIBLK-PS030124.D	2,4-DCAA	1	500	419	84		39	175
		2,4-DCAA	2	500	418	84		39	175
I.BLK-PS030379.D	PIBLK-PS030379.D	2,4-DCAA	1	500	507	101		39	175
		2,4-DCAA	2	500	493	99		39	175
PB168137BL	PB168137BL	2,4-DCAA	1	500	449	90		10	141
		2,4-DCAA	2	500	406	81		10	141
PB168137BS	PB168137BS	2,4-DCAA	1	500	510	102		10	141
		2,4-DCAA	2	500	435	87		10	141
I.BLK-PS030391.D	PIBLK-PS030391.D	2,4-DCAA	1	500	531	106		39	175
		2,4-DCAA	2	500	526	105		39	175
Q2100-01	TP-16	2,4-DCAA	1	500	359	72		10	141
		2,4-DCAA	2	500	322	64		10	141
Q2100-02	TP-22	2,4-DCAA	1	500	299	60		10	141
		2,4-DCAA	2	500	279	56		10	141
Q2100-03	TP-21	2,4-DCAA	1	500	252	50		10	141
		2,4-DCAA	2	500	234	47		10	141
I.BLK-PS030403.D	PIBLK-PS030403.D	2,4-DCAA	1	500	522	104		39	175
		2,4-DCAA	2	500	531	106		39	175
Q2100-04	TP-45	2,4-DCAA	1	500	280	56		10	141
		2,4-DCAA	2	500	254	51		10	141
I.BLK-PS030408.D	PIBLK-PS030408.D	2,4-DCAA	1	500	531	106		39	175
		2,4-DCAA	2	500	535	107		39	175
I.BLK-PS030426.D	PIBLK-PS030426.D	2,4-DCAA	1	500	488	98		39	175
		2,4-DCAA	2	500	482	96		39	175
Q2100-03MS	TP-21MS	2,4-DCAA	1	500	282	56		10	141
		2,4-DCAA	2	500	249	50		10	141
Q2100-03MSD	TP-21MSD	2,4-DCAA	1	500	293	59		10	141
		2,4-DCAA	2	500	262	52		10	141
I.BLK-PS030433.D	PIBLK-PS030433.D	2,4-DCAA	1	500	499	100		39	175
		2,4-DCAA	2	500	483	97		39	175