

Surrogate Summary

SW-846

SDG No.: Q2149

Client: Aramark Uniforms

Analytical Method: 8270E

Lab Sample ID	Client ID	Parameter	Spike (PPM)	Result (PPM)	Recovery (%)	Qual	Limits (%)	
							Low	High
PB168190TB	PB168190TB	2-Fluorophenol	150	120	80		23	138
		Phenol-d6	150	121	81		10	134
		Nitrobenzene-d5	100	75.3	75		67	132
		2-Fluorobiphenyl	100	71.7	72		52	132
		2,4,6-Tribromophenol	150	124	83		44	137
		Terphenyl-d14	100	69.4	69		42	152
PB168219BL	PB168219BL	2-Fluorophenol	150	120	80		23	138
		Phenol-d6	150	120	80		10	134
		Nitrobenzene-d5	100	74.7	75		67	132
		2-Fluorobiphenyl	100	71.6	72		52	132
		2,4,6-Tribromophenol	150	124	83		44	137
		Terphenyl-d14	100	69.7	70		42	152
PB168219BS	PB168219BS	2-Fluorophenol	150	122	81		23	138
		Phenol-d6	150	123	82		10	134
		Nitrobenzene-d5	100	75.3	75		67	132
		2-Fluorobiphenyl	100	71.0	71		52	132
		2,4,6-Tribromophenol	150	126	84		44	137
		Terphenyl-d14	100	76.9	77		42	152
Q2136-05MS	OR-646-COMP-52MS	2-Fluorophenol	150	115	77		23	138
		Phenol-d6	150	108	72		10	134
		Nitrobenzene-d5	100	79.5	79		67	132
		2-Fluorobiphenyl	100	77.9	78		52	132
		2,4,6-Tribromophenol	150	126	84		44	137
		Terphenyl-d14	100	76.9	77		42	152
Q2136-05MSD	OR-646-COMP-52MSD	2-Fluorophenol	150	106	71		23	138
		Phenol-d6	150	99.1	66		10	134
		Nitrobenzene-d5	100	73.2	73		67	132
		2-Fluorobiphenyl	100	71.4	71		52	132
		2,4,6-Tribromophenol	150	114	76		44	137
		Terphenyl-d14	100	71.5	72		42	152
Q2149-01	FILTER-CAKE	2-Fluorophenol	150	102	68		23	138
		Phenol-d6	150	91.2	61		10	134
		Nitrobenzene-d5	100	74.0	74		67	132
		2-Fluorobiphenyl	100	73.4	73		52	132
		2,4,6-Tribromophenol	150	97.1	65		44	137
		Terphenyl-d14	100	45.2	45		42	152
Q2149-01DL	FILTER-CAKEDL	2-Fluorophenol	150	129	86		23	138
		Phenol-d6	150	118	78		10	134
		Nitrobenzene-d5	100	93.9	94		67	132
		2-Fluorobiphenyl	100	95.0	95		52	132
		2,4,6-Tribromophenol	150	122	81		44	137
		Terphenyl-d14	100	58.5	59		42	152