

Surrogate Summary

SDG No.: Q2176

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030475.D	PIBLK-PS030475.D	2,4-DCAA	1	500	414	83		61	136
		2,4-DCAA	2	500	445	89		61	136
I.BLK-PS030496.D	PIBLK-PS030496.D	2,4-DCAA	1	500	445	89		61	136
		2,4-DCAA	2	500	459	92		61	136
Q2160-05MS	TP05-MHH-WCMS	2,4-DCAA	1	500	494	99		10	141
		2,4-DCAA	2	500	480	96		10	141
Q2160-05MSD	TP05-MHH-WCMSD	2,4-DCAA	1	500	496	99		10	141
		2,4-DCAA	2	500	499	100		10	141
I.BLK-PS030505.D	PIBLK-PS030505.D	2,4-DCAA	1	500	442	88		61	136
		2,4-DCAA	2	500	460	92		61	136
PB168254BL	PB168254BL	2,4-DCAA	1	500	435	87		10	141
		2,4-DCAA	2	500	472	94		10	141
PB168254BS	PB168254BS	2,4-DCAA	1	500	508	102		10	141
		2,4-DCAA	2	500	502	100		10	141
I.BLK-PS030515.D	PIBLK-PS030515.D	2,4-DCAA	1	500	447	89		61	136
		2,4-DCAA	2	500	460	92		61	136
Q2176-01	TP-46	2,4-DCAA	1	500	320	64		10	141
		2,4-DCAA	2	500	339	68		10	141
Q2176-02	TP-56	2,4-DCAA	1	500	300	60		10	141
		2,4-DCAA	2	500	290	58		10	141
Q2176-03	TP-25	2,4-DCAA	1	500	317	63		10	141
		2,4-DCAA	2	500	340	68		10	141
Q2176-05	TP-28	2,4-DCAA	1	500	346	69		10	141
		2,4-DCAA	2	500	306	61		10	141
Q2176-06	TP-27	2,4-DCAA	1	500	247	49		10	141
		2,4-DCAA	2	500	254	51		10	141
Q2176-07	TP-31	2,4-DCAA	1	500	333	67		10	141
		2,4-DCAA	2	500	346	69		10	141
Q2176-08	TP-65	2,4-DCAA	1	500	316	63		10	141
		2,4-DCAA	2	500	360	72		10	141
I.BLK-PS030525.D	PIBLK-PS030525.D	2,4-DCAA	1	500	447	89		61	136
		2,4-DCAA	2	500	462	92		61	136
I.BLK-PS030528.D	PIBLK-PS030528.D	2,4-DCAA	1	500	448	90		61	136
		2,4-DCAA	2	500	466	93		61	136
Q2176-04	TP-26	2,4-DCAA	1	500	71.0	14		10	141
		2,4-DCAA	2	500	71.8	14		10	141
I.BLK-PS030531.D	PIBLK-PS030531.D	2,4-DCAA	1	500	451	90		61	136
		2,4-DCAA	2	500	457	91		61	136