

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE051425AL

Review By	yogesh	Review On	5/14/2025 11:46:43 AM
Supervise By	mohammad	Supervise On	5/16/2025 1:55:25 AM
SubDirectory	FE051425AL	HP Acquire Method	HP Processing Method FE051425AL

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178
CCC Internal Standard/PEM ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24176 PP24174,PP24179

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	MECL2	MECL2	FE053806.D	14 May 2025 08:48		YP\AJ	Ok
2	I.BLK	I.BLK	FE053807.D	14 May 2025 09:18		YP\AJ	Ok
3	100 PPM ALIPHATIC HC	100 PPM ALIPHATIC HC	FE053808.D	14 May 2025 09:48		YP\AJ	Ok
4	50 PPM ALIPHATIC HC	50 PPM ALIPHATIC HC	FE053809.D	14 May 2025 10:18		YP\AJ	Ok
5	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053810.D	14 May 2025 10:48		YP\AJ	Ok
6	10 PPM ALIPHATIC HC	10 PPM ALIPHATIC HC	FE053811.D	14 May 2025 11:18		YP\AJ	Ok
7	5 PPM ALIPHATIC HC	5 PPM ALIPHATIC HC	FE053812.D	14 May 2025 11:48		YP\AJ	Ok
8	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053813.D	14 May 2025 12:18		YP\AJ	Ok
9	I.BLK	I.BLK	FE053814.D	14 May 2025 12:48		YP\AJ	Ok
10	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053815.D	14 May 2025 13:18		YP\AJ	Ok
11	Q2019-01	MH-K	FE053816.D	14 May 2025 14:44		YP\AJ	Ok
12	Q2019-02	MH-K-EPH	FE053817.D	14 May 2025 15:13		YP\AJ	Ok
13	Q2020-01	TP-A	FE053818.D	14 May 2025 15:44		YP\AJ	Ok
14	Q2020-02	TP-A-EPH	FE053819.D	14 May 2025 16:14		YP\AJ	Ok
15	Q2022-02	COMP-1	FE053820.D	14 May 2025 16:44		YP\AJ	Ok
16	Q2022-03	COMP-1-EPH	FE053821.D	14 May 2025 17:14		YP\AJ	Ok
17	Q2022-04	COMP-2	FE053822.D	14 May 2025 18:04		YP\AJ	Ok
18	Q2022-05	COMP-2-EPH	FE053823.D	14 May 2025 18:34		YP\AJ	Ok

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE051425AL

Review By	yogesh	Review On	5/14/2025 11:46:43 AM
Supervise By	mohammad	Supervise On	5/16/2025 1:55:25 AM
SubDirectory	FE051425AL	HP Acquire Method	HP Processing Method FE051425AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC	PP24176		
Internal Standard/PEM	PP24174,PP24179		
ICV/I.BLK			
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

19	Q2024-01	PL-02-051325	FE053824.D	14 May 2025 19:04		YP\AJ	Ok,M
20	Q2024-02	PL-02-051325-E2	FE053825.D	14 May 2025 19:34		YP\AJ	Ok
21	Q2014-01	MOO-25-0134	FE053826.D	14 May 2025 20:04	need 50x dilution	YP\AJ	Dilution
22	Q2014-02	MOO-25-0134-E2	FE053827.D	14 May 2025 20:34	need 5x & 200x dilution	YP\AJ	Dilution
23	Q2014-03	MOO-25-0145	FE053828.D	14 May 2025 21:05	need 100x dilution	YP\AJ	Dilution
24	Q2014-04	MOO-25-0145-E2	FE053829.D	14 May 2025 21:35	need 100x dilution	YP\AJ	Dilution
25	Q2014-05	MOO-25-0148	FE053830.D	14 May 2025 22:05	need 1000x dilution	YP\AJ	Dilution
26	Q2014-06	MOO-25-0148-E2	FE053831.D	14 May 2025 22:35	need 500x dilution	YP\AJ	Dilution
27	I.BLK	I.BLK	FE053832.D	14 May 2025 23:35		YP\AJ	Ok
28	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053833.D	15 May 2025 00:05		YP\AJ	Ok

M : Manual Integration

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE060325AL

Review By	yogesh	Review On	6/3/2025 10:21:26 AM
Supervise By	mohammad	Supervise On	6/5/2025 1:34:56 AM
SubDirectory	FE060325AL	HP Acquire Method	HP Processing Method FE051425AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC	PP24176		
Internal Standard/PEM	PP24174,PP24179		
ICV/I.BLK			
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	MECL2	MECL2	FE054140.D	03 Jun 2025 06:43		YP\AJ	Ok
2	I.BLK	I.BLK	FE054141.D	03 Jun 2025 07:13		YP\AJ	Ok
3	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE054142.D	03 Jun 2025 07:43		YP\AJ	Ok
4	Q2152-01DL	OK-02-05292025DL	FE054143.D	03 Jun 2025 08:13		YP\AJ	Ok
5	Q2152-02DL	OK-02-05292025-E2DL	FE054144.D	03 Jun 2025 08:43		YP\AJ	Ok
6	Q2153-01DL	TR-04-0592025DL	FE054145.D	03 Jun 2025 09:14		YP\AJ	Ok
7	Q2153-02DL	TR-04-0592025-E2DL	FE054146.D	03 Jun 2025 09:43		YP\AJ	Ok
8	I.BLK	I.BLK	FE054147.D	03 Jun 2025 10:44		YP\AJ	Ok
9	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE054148.D	03 Jun 2025 11:14		YP\AJ	Ok
10	PB168239BL	PB168239BL	FE054149.D	03 Jun 2025 13:54		YP\AJ	Ok
11	PB168239BS	PB168239BS	FE054150.D	03 Jun 2025 14:24		YP\AJ	Ok,M
12	PB168239BSD	PB168239BSD	FE054151.D	03 Jun 2025 14:55		YP\AJ	Ok,M
13	Q2147-03	1	FE054152.D	03 Jun 2025 15:25		YP\AJ	Ok
14	Q2147-04	2	FE054153.D	03 Jun 2025 15:55		YP\AJ	Ok
15	Q2147-05	3	FE054154.D	03 Jun 2025 16:25		YP\AJ	Ok
16	Q2147-05MS	3MS	FE054155.D	03 Jun 2025 16:56	FE054154.D	YP\AJ	Ok,M
17	Q2147-05MSD	3MSD	FE054156.D	03 Jun 2025 17:26	FE054154.D!FE054155.D	YP\AJ	Ok,M
18	Q2147-06	4	FE054157.D	03 Jun 2025 17:56		YP\AJ	Ok

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE060325AL

Review By	yogesh	Review On	6/3/2025 10:21:26 AM
Supervise By	mohammad	Supervise On	6/5/2025 1:34:56 AM
SubDirectory	FE060325AL	HP Acquire Method	HP Processing Method FE051425AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC	PP24176		
Internal Standard/PEM	PP24174,PP24179		
ICV/I.BLK			
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

19	Q2147-07	5	FE054158.D	03 Jun 2025 18:26		YPIAJ	Ok
20	Q2174-01	DP1	FE054159.D	03 Jun 2025 18:56		YPIAJ	Ok
21	Q2174-02	D1	FE054160.D	03 Jun 2025 19:26		YPIAJ	Ok
22	Q2174-03	D2	FE054161.D	03 Jun 2025 19:57		YPIAJ	Ok,M
23	Q2179-01	D3	FE054162.D	03 Jun 2025 20:27		YPIAJ	Ok
24	Q2179-02	D4	FE054163.D	03 Jun 2025 20:57		YPIAJ	Ok
25	Q2179-03	D5	FE054164.D	03 Jun 2025 21:27		YPIAJ	Ok
26	Q2179-04	D6	FE054165.D	03 Jun 2025 21:57		YPIAJ	Ok
27	I.BLK	I.BLK	FE054166.D	03 Jun 2025 22:57		YPIAJ	Ok
28	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE054167.D	03 Jun 2025 23:27		YPIAJ	Ok,M

M : Manual Integration