

### Surrogate Summary

SDG No.: Q2235

Client: ENTACT

Analytical Method: 8081B

| Lab Sample ID    | Client ID        | Parameter            | Column | Spike | Result | Rec | Qual | Limits  |           |
|------------------|------------------|----------------------|--------|-------|--------|-----|------|---------|-----------|
|                  |                  |                      |        |       |        |     |      | Low     | High      |
| I.BLK-PD088583.D | PIBLK-PD088583.D | Decachlorobiphenyl   | 1      | 20    | 18.4   | 92  |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 16.5   | 83  |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 18.0   | 90  |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 16.9   | 84  |      | 30 (61) | 150 (148) |
| I.BLK-PD088820.D | PIBLK-PD088820.D | Decachlorobiphenyl   | 1      | 20    | 24.1   | 120 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 23.7   | 119 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 23.0   | 115 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 21.4   | 107 |      | 30 (61) | 150 (148) |
| PB168330BL       | PB168330BL       | Decachlorobiphenyl   | 1      | 20    | 25.8   | 129 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 24.1   | 121 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 24.4   | 122 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 21.8   | 109 |      | 30 (61) | 150 (148) |
| PB168330BS       | PB168330BS       | Decachlorobiphenyl   | 1      | 20    | 23.8   | 119 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 22.0   | 110 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 22.9   | 115 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 20.4   | 102 |      | 30 (61) | 150 (148) |
| PB168311TB       | PB168311TB       | Decachlorobiphenyl   | 1      | 20    | 25.4   | 127 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 23.7   | 118 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 24.3   | 121 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 21.5   | 107 |      | 30 (61) | 150 (148) |
| Q2198-04MS       | B-207-SB02MS     | Decachlorobiphenyl   | 1      | 20    | 25.9   | 129 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 24.2   | 121 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 24.8   | 124 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 21.1   | 105 |      | 30 (61) | 150 (148) |
| Q2198-04MSD      | B-207-SB02MSD    | Decachlorobiphenyl   | 1      | 20    | 26.1   | 131 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 24.3   | 122 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 24.9   | 124 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 21.3   | 106 |      | 30 (61) | 150 (148) |
| I.BLK-PD088840.D | PIBLK-PD088840.D | Decachlorobiphenyl   | 1      | 20    | 24.7   | 123 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 23.3   | 117 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 23.7   | 118 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 20.8   | 104 |      | 30 (61) | 150 (148) |
| I.BLK-PD088853.D | PIBLK-PD088853.D | Decachlorobiphenyl   | 1      | 20    | 25.2   | 126 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 24.7   | 123 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 24.6   | 123 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 22.5   | 113 |      | 30 (61) | 150 (148) |
| Q2235-03         | WC-A2-08-C       | Decachlorobiphenyl   | 1      | 20    | 23.9   | 120 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 1      | 20    | 21.5   | 108 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 21.6   | 108 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 19.6   | 98  |      | 30 (61) | 150 (148) |
| I.BLK-PD088865.D | PIBLK-PD088865.D | Decachlorobiphenyl   | 1      | 20    | 25.0   | 125 |      | 30 (57) | 150 (171) |

() = LABORATORY INHOUSE LIMIT

### Surrogate Summary

SDG No.: Q2235

Client: ENTACT

Analytical Method: 8081B

| Lab Sample ID    | Client ID        | Parameter            | Column | Spike | Result | Rec | Qual | Limits  |           |
|------------------|------------------|----------------------|--------|-------|--------|-----|------|---------|-----------|
|                  |                  |                      |        |       |        |     |      | Low     | High      |
| I.BLK-PD088865.D | PIBLK-PD088865.D | Tetrachloro-m-xylene | 1      | 20    | 24.2   | 121 |      | 30 (61) | 150 (148) |
|                  |                  | Decachlorobiphenyl   | 2      | 20    | 24.3   | 121 |      | 30 (57) | 150 (171) |
|                  |                  | Tetrachloro-m-xylene | 2      | 20    | 22.3   | 112 |      | 30 (61) | 150 (148) |