

Surrogate Summary

SDG No.: Q2236

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PO111057.D	PIBLK-PO111057.D	Tetrachloro-m-xylene	1	20	16.4	82		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	17.5	87		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	15.4	77		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	17.9	90		70 (60)	130 (140)
I.BLK-PO111516.D	PIBLK-PO111516.D	Tetrachloro-m-xylene	1	20	20.9	104		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	19.6	98		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	19.0	95		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	18.3	91		70 (60)	130 (140)
Q2236-02	WC-A4-05A-C	Tetrachloro-m-xylene	1	20	16.7	83		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	12.3	62		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	19.3	97		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	12.1	61		30 (32)	150 (175)
Q2236-06	WC-A2-04-C	Tetrachloro-m-xylene	1	20	14.8	74		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	13.6	68		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	15.5	77		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	13.9	70		30 (32)	150 (175)
Q2236-10	WC-A2-05-C	Tetrachloro-m-xylene	1	20	14.9	75		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	9.63	48		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	15.0	75		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	9.72	49		30 (32)	150 (175)
Q2236-14	WC-A2-06-C	Tetrachloro-m-xylene	1	20	20.4	102		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	16.4	82		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	21.2	106		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	16.0	80		30 (32)	150 (175)
Q2236-18	WC-A2-07-C	Tetrachloro-m-xylene	1	20	18.2	91		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	16.4	82		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	17.0	85		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	16.5	83		30 (32)	150 (175)
I.BLK-PO111532.D	PIBLK-PO111532.D	Tetrachloro-m-xylene	1	20	20.1	100		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	17.5	88		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	18.7	94		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	17.5	87		70 (60)	130 (140)
I.BLK-PP072163.D	PIBLK-PP072163.D	Tetrachloro-m-xylene	1	20	16.3	82		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	16.7	84		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	16.6	83		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	17.0	85		70 (60)	130 (140)
I.BLK-PP072687.D	PIBLK-PP072687.D	Tetrachloro-m-xylene	1	20	17.6	88		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	18.9	94		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	19.3	97		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	21.1	106		70 (60)	130 (140)
PB168319BL	PB168319BL	Tetrachloro-m-xylene	1	20	19.7	99		30 (32)	150 (144)

() = LABORATORY INHOUSE LIMIT

Surrogate Summary

SDG No.: Q2236

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
PB168319BL	PB168319BL	Decachlorobiphenyl	1	20	20.7	104		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	21.0	105		30 (32)	150 (144)
PB168319BS	PB168319BS	Decachlorobiphenyl	2	20	23.4	117		30 (32)	150 (175)
		Tetrachloro-m-xylene	1	20	20.2	101		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	21.6	108		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	20.5	103		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	25.0	125		30 (32)	150 (175)
		Tetrachloro-m-xylene	1	20	15.9	80		30 (32)	150 (144)
Q2223-01MS	HR-01-06042025MS	Decachlorobiphenyl	1	20	18.4	92		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	17.5	88		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	21.9	109		30 (32)	150 (175)
		Tetrachloro-m-xylene	1	20	14.0	70		30 (32)	150 (144)
Q2223-01MSD	HR-01-06042025MSD	Decachlorobiphenyl	1	20	15.9	79		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	15.4	77		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	19.0	95		30 (32)	150 (175)
		Tetrachloro-m-xylene	1	20	17.0	85		70 (60)	130 (140)
I.BLK-PP072702.D	PIBLK-PP072702.D	Decachlorobiphenyl	1	20	15.7	78		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	18.8	94		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	19.9	99		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20					