

Surrogate Summary

SDG No.: Q2301

Client: ENTACT

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030680.D	PIBLK-PS030680.D	2,4-DCAA	1	500	396	79		70 (61)	130 (136)
		2,4-DCAA	2	500	472	94		70 (61)	130 (136)
I.BLK-PS030697.D	PIBLK-PS030697.D	2,4-DCAA	1	500	446	89		70 (61)	130 (136)
		2,4-DCAA	2	500	477	95		70 (61)	130 (136)
PB168495BL	PB168495BL	2,4-DCAA	1	500	427	85		70 (61)	130 (136)
		2,4-DCAA	2	500	470	94		70 (61)	130 (136)
PB168495BS	PB168495BS	2,4-DCAA	1	500	501	100		70 (61)	130 (136)
		2,4-DCAA	2	500	500	100		70 (61)	130 (136)
PB168428TB	PB168428TB	2,4-DCAA	1	500	448	90		70 (61)	130 (136)
		2,4-DCAA	2	500	451	90		70 (61)	130 (136)
Q2301-03	WC-URBAN-FILL-C	2,4-DCAA	1	500	440	88		70 (61)	130 (136)
		2,4-DCAA	2	500	448	90		70 (61)	130 (136)
Q2301-03MS	WC-URBAN-FILL-CMS	2,4-DCAA	1	500	484	97		70 (61)	130 (136)
		2,4-DCAA	2	500	465	93		70 (61)	130 (136)
Q2301-03MSD	WC-URBAN-FILL-CMSD	2,4-DCAA	1	500	490	98		70 (61)	130 (136)
		2,4-DCAA	2	500	476	95		70 (61)	130 (136)
I.BLK-PS030705.D	PIBLK-PS030705.D	2,4-DCAA	1	500	440	88		70 (61)	130 (136)
		2,4-DCAA	2	500	478	96		70 (61)	130 (136)