

Surrogate Summary

SDG No.: Q2413

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030738.D	PIBLK-PS030738.D	2,4-DCAA	1	500	414	83		61	136
		2,4-DCAA	2	500	487	97		61	136
I.BLK-PS030837.D	PIBLK-PS030837.D	2,4-DCAA	1	500	530	106		61	136
		2,4-DCAA	2	500	577	115		61	136
PB168634BS	PB168634BS	2,4-DCAA	1	500	515	103		10	141
		2,4-DCAA	2	500	507	101		10	141
Q2405-01MS	MH-M/NMS	2,4-DCAA	1	500	309	62		10	141
		2,4-DCAA	2	500	307	61		10	141
Q2405-01MSD	MH-M/NMSD	2,4-DCAA	1	500	309	62		10	141
		2,4-DCAA	2	500	308	62		10	141
Q2413-01	TP-35	2,4-DCAA	1	500	216	43		10	141
		2,4-DCAA	2	500	232	46		10	141
Q2413-02	TP-34	2,4-DCAA	1	500	289	58		10	141
		2,4-DCAA	2	500	263	53		10	141
Q2413-03	TP-77	2,4-DCAA	1	500	415	83		10	141
		2,4-DCAA	2	500	610	122		10	141
Q2413-04	TP-74	2,4-DCAA	1	500	335	67		10	141
		2,4-DCAA	2	500	550	110		10	141
Q2413-05	TP-73	2,4-DCAA	1	500	253	51		10	141
		2,4-DCAA	2	500	273	55		10	141
I.BLK-PS030849.D	PIBLK-PS030849.D	2,4-DCAA	1	500	532	106		61	136
		2,4-DCAA	2	500	573	115		61	136
Q2413-06	TP-72	2,4-DCAA	1	500	238	48		10	141
		2,4-DCAA	2	500	211	42		10	141
I.BLK-PS030858.D	PIBLK-PS030858.D	2,4-DCAA	1	500	549	110		61	136
		2,4-DCAA	2	500	591	118		61	136
I.BLK-PS030861.D	PIBLK-PS030861.D	2,4-DCAA	1	500	464	93		61	136
		2,4-DCAA	2	500	484	97		61	136
PB168634BL	PB168634BL	2,4-DCAA	1	500	446	89		10	141
		2,4-DCAA	2	500	563	113		10	141
I.BLK-PS030865.D	PIBLK-PS030865.D	2,4-DCAA	1	500	464	93		61	136
		2,4-DCAA	2	500	486	97		61	136