

Surrogate Summary

SDG No.: Q2436

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS030738.D	PIBLK-PS030738.D	2,4-DCAA	1	500	414	83		61	136
		2,4-DCAA	2	500	487	97		61	136
I.BLK-PS030880.D	PIBLK-PS030880.D	2,4-DCAA	1	500	500	100		61	136
		2,4-DCAA	2	500	524	105		61	136
PB168675BL	PB168675BL	2,4-DCAA	1	500	320	64		10	141
		2,4-DCAA	2	500	407	81		10	141
PB168675BS	PB168675BS	2,4-DCAA	1	500	511	102		10	141
		2,4-DCAA	2	500	494	99		10	141
Q2436-01	TP-70	2,4-DCAA	1	500	289	58		10	141
		2,4-DCAA	2	500	324	65		10	141
Q2436-02	TP-69	2,4-DCAA	1	500	226	45		10	141
		2,4-DCAA	2	500	250	50		10	141
Q2436-03	TP-85	2,4-DCAA	1	500	228	46		10	141
		2,4-DCAA	2	500	252	50		10	141
Q2436-04	TP-86	2,4-DCAA	1	500	283	57		10	141
		2,4-DCAA	2	500	326	65		10	141
Q2436-05	TP-84	2,4-DCAA	1	500	366	73		10	141
		2,4-DCAA	2	500	400	80		10	141
Q2436-06	TP-83	2,4-DCAA	1	500	338	68		10	141
		2,4-DCAA	2	500	459	92		10	141
Q2436-07	TP-87	2,4-DCAA	1	500	349	70		10	141
		2,4-DCAA	2	500	405	81		10	141
Q2436-08	TP-100	2,4-DCAA	1	500	267	53		10	141
		2,4-DCAA	2	500	308	62		10	141
I.BLK-PS030892.D	PIBLK-PS030892.D	2,4-DCAA	1	500	498	100		61	136
		2,4-DCAA	2	500	517	103		61	136
Q2436-09	TP-99	2,4-DCAA	1	500	152	30		10	141
		2,4-DCAA	2	500	176	35		10	141
Q2436-10	TP-82	2,4-DCAA	1	500	215	43		10	141
		2,4-DCAA	2	500	246	49		10	141
Q2436-10MS	TP-82MS	2,4-DCAA	1	500	291	58		10	141
		2,4-DCAA	2	500	328	66		10	141
Q2436-10MSD	TP-82MSD	2,4-DCAA	1	500	317	63		10	141
		2,4-DCAA	2	500	321	64		10	141
I.BLK-PS030898.D	PIBLK-PS030898.D	2,4-DCAA	1	500	507	101		61	136
		2,4-DCAA	2	500	528	106		61	136