

Surrogate Summary

SDG No.: Q2458

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS030738.D	PIBLK-PS030738.D	2,4-DCAA	1	500	414	83		61	136
		2,4-DCAA	2	500	487	97		61	136
I.BLK-PS030916.D	PIBLK-PS030916.D	2,4-DCAA	1	500	421	84		61	136
		2,4-DCAA	2	500	481	96		61	136
Q2458-02	TP-55	2,4-DCAA	1	500	60.0	12		10	141
		2,4-DCAA	2	500	60.6	12		10	141
Q2458-03	TP-68	2,4-DCAA	1	500	262	52		10	141
		2,4-DCAA	2	500	250	50		10	141
Q2458-04	TP-67	2,4-DCAA	1	500	350	70		10	141
		2,4-DCAA	2	500	352	70		10	141
Q2458-05	TP-66	2,4-DCAA	1	500	216	43		10	141
		2,4-DCAA	2	500	260	52		10	141
Q2458-06	TP-60	2,4-DCAA	1	500	166	33		10	141
		2,4-DCAA	2	500	174	35		10	141
Q2458-07	TP-62	2,4-DCAA	1	500	69.9	14		10	141
		2,4-DCAA	2	500	75.5	15		10	141
Q2458-08	TP-63	2,4-DCAA	1	500	102	20		10	141
		2,4-DCAA	2	500	121	24		10	141
Q2458-09	TP-59	2,4-DCAA	1	500	213	43		10	141
		2,4-DCAA	2	500	223	45		10	141
I.BLK-PS030927.D	PIBLK-PS030927.D	2,4-DCAA	1	500	468	94		61	136
		2,4-DCAA	2	500	502	100		61	136
I.BLK-PS030943.D	PIBLK-PS030943.D	2,4-DCAA	1	500	481	96		61	136
		2,4-DCAA	2	500	497	99		61	136
Q2458-01	TP-76	2,4-DCAA	1	500	278	56		10	141
		2,4-DCAA	2	500	275	55		10	141
PB168706BL	PB168706BL	2,4-DCAA	1	500	452	90		10	141
		2,4-DCAA	2	500	489	98		10	141
Q2458-01MS	TP-76MS	2,4-DCAA	1	500	341	68		10	141
		2,4-DCAA	2	500	375	75		10	141
Q2458-01MSD	TP-76MSD	2,4-DCAA	1	500	347	69		10	141
		2,4-DCAA	2	500	378	76		10	141
PB168717BL	PB168717BL	2,4-DCAA	1	500	448	90		61	136
		2,4-DCAA	2	500	479	96		61	136
PB168717BS	PB168717BS	2,4-DCAA	1	500	534	107		61	136
		2,4-DCAA	2	500	509	102		61	136
PB168717BSD	PB168717BSD	2,4-DCAA	1	500	539	108		61	136
		2,4-DCAA	2	500	510	102		61	136
Q2458-10	FB-06272025	2,4-DCAA	1	500	437	87		61	136
		2,4-DCAA	2	500	448	90		61	136
I.BLK-PS030955.D	PIBLK-PS030955.D	2,4-DCAA	1	500	455	91		61	136

Surrogate Summary

SDG No.: Q2458

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS030955.D	PIBLK-PS030955.D	2,4-DCAA	2	500	469	94		61	136
I.BLK-PS030958.D	PIBLK-PS030958.D	2,4-DCAA	1	500	456	91		61	136
		2,4-DCAA	2	500	462	92		61	136
PB168706BS	PB168706BS	2,4-DCAA	1	500	538	108		10	141
		2,4-DCAA	2	500	503	101		10	141
I.BLK-PS030966.D	PIBLK-PS030966.D	2,4-DCAA	1	500	449	90		61	136
		2,4-DCAA	2	500	467	93		61	136