

SOP ID : M1311-TCLP-16
SDG No : N/A
Weigh By : JP
Balance ID : WC SC-7
pH Meter ID : WC PH METER-1
Extraction By : JP
Filter By : JP
Pipette ID : WC
Tumbler ID : T-1 / T-2
TCLP Filter ID : 115525

Start Prep Date : 07/02/2025 Time : 15:00
End Prep Date : 07/03/2025 Time : 08:20
Combination Ratio : 20
ZHE Cleaning Batch : 18A
Initial Room Temperature: 24 °C
Final Room Temperature: 22 °C
TCLP Technician Signature : *SP*
Supervisor By : 12

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
TCLP-FLUID-1	N/A	WP112795
HCL-TCLP, 1N	N/A	WP112797
HNO3-TCLP, 1N	N/A	WP112799
pH Strips	N/A	W1931, W1934, W3171, W3172
pH Strips	W1940, W1941, W1942	W3166, W1938, W1939,
1 Liter Amber	N/A	90924-08
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MPB4041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. TUMBLER T-1 / T-2 checked, 30 rpm. Particle size reduction is not required. q2493-04 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
07/03/25 11:00	<i>SP</i> HGA Room	SLB. RJ 1641
	Preparation Group	Analysis Group <i>Detm'y</i>