

SOP ID : M1311-TCLP-16

SDG No : N/A

Weight By : N/A

Balance ID : N/A

pH Meter ID : WC PH METER-1

Extraction By : JP

Filter By : JP

Pipette ID : N/A

Tumbler ID : N/A

TCLP Filter ID : 115525

Start Prep Date : N/A Time : N/A

End Prep Date : N/A Time : N/A

Combination Ratio : N/A

ZHE Cleaning Batch : N/A

Initial Room Temperature: N/A

Final Room Temperature: N/A

TCLP Technician Signature : *JP*

Supervisor By : *12*

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
N/A	N/A	N/A
N/A	N/A	N/A
HNO3-TCLP, 1N	N/A	WP112799
pH Strips	N/A	W1931, W1934, W3171, W3172
pH Strips	N/A	W3166, W1938, W1939,
N/A	N/A	N/A
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. q2481-01 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
07/03/15 11:00	<i>JP</i> <i>100 Room</i>	<i>SUB</i> <i>104019</i>
	Preparation Group	Analysis Group