

Surrogate Summary

SDG No.: Q2484

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS030738.D	PIBLK-PS030738.D	2,4-DCAA	1	500	414	83		61	136
		2,4-DCAA	2	500	487	97		61	136
I.BLK-PS030931.D	PIBLK-PS030931.D	2,4-DCAA	1	500	470	94		61	136
		2,4-DCAA	2	500	497	99		61	136
Q2484-01	TP-58	2,4-DCAA	1	500	263	53		10	141
		2,4-DCAA	2	500	274	55		10	141
Q2484-02	TP-57	2,4-DCAA	1	500	304	61		10	141
		2,4-DCAA	2	500	285	57		10	141
Q2484-04	TP-107	2,4-DCAA	1	500	385	77		10	141
		2,4-DCAA	2	500	437	87		10	141
Q2484-05	TP-106	2,4-DCAA	1	500	514	103		10	141
		2,4-DCAA	2	500	607	121		10	141
Q2484-06	TP-104	2,4-DCAA	1	500	363	73		10	141
		2,4-DCAA	2	500	395	79		10	141
I.BLK-PS030940.D	PIBLK-PS030940.D	2,4-DCAA	1	500	472	94		61	136
		2,4-DCAA	2	500	494	99		61	136
I.BLK-PS030943.D	PIBLK-PS030943.D	2,4-DCAA	1	500	481	96		61	136
		2,4-DCAA	2	500	497	99		61	136
PB168706BL	PB168706BL	2,4-DCAA	1	500	452	90		10	141
		2,4-DCAA	2	500	489	98		10	141
Q2458-01MS	TP-76MS	2,4-DCAA	1	500	341	68		10	141
		2,4-DCAA	2	500	375	75		10	141
Q2458-01MSD	TP-76MSD	2,4-DCAA	1	500	347	69		10	141
		2,4-DCAA	2	500	378	76		10	141
I.BLK-PS030955.D	PIBLK-PS030955.D	2,4-DCAA	1	500	455	91		61	136
		2,4-DCAA	2	500	469	94		61	136
I.BLK-PS030958.D	PIBLK-PS030958.D	2,4-DCAA	1	500	456	91		61	136
		2,4-DCAA	2	500	462	92		61	136
PB168706BS	PB168706BS	2,4-DCAA	1	500	538	108		10	141
		2,4-DCAA	2	500	503	101		10	141
I.BLK-PS030966.D	PIBLK-PS030966.D	2,4-DCAA	1	500	449	90		61	136
		2,4-DCAA	2	500	467	93		61	136
I.BLK-PS030979.D	PIBLK-PS030979.D	2,4-DCAA	1	500	474	95		61	136
		2,4-DCAA	2	500	473	95		61	136
Q2484-03	TP-64	2,4-DCAA	1	500	96.8	19		10	141
		2,4-DCAA	2	500	96.2	19		10	141
I.BLK-PS030989.D	PIBLK-PS030989.D	2,4-DCAA	1	500	476	95		61	136
		2,4-DCAA	2	500	476	95		61	136