

**Instrument ID:** FILTER/GRAVIMETRIC

**Daily Analysis Runlog For Sequence/QCBatch ID # LB136381**

|                  |                  |              |                      |
|------------------|------------------|--------------|----------------------|
| Review By        | rubina           | Review On    | 7/7/2025 12:21:49 PM |
| Supervise By     | Iwona            | Supervise On | 7/7/2025 12:22:58 PM |
| SubDirectory     | LB136381         | Test         | Paint Filter         |
| <b>STD. NAME</b> | <b>STD REF.#</b> |              |                      |
| ICAL Standard    | N/A              |              |                      |
| ICV Standard     | N/A              |              |                      |
| CCV Standard     | N/A              |              |                      |
| ICSA Standard    | N/A              |              |                      |
| CRI Standard     | N/A              |              |                      |
| LCS Standard     | N/A              |              |                      |
| Chk Standard     | N/A              |              |                      |

| Sr# | SampleId    | ClientID       | QcType | Date           | Comment | Operator | Status |
|-----|-------------|----------------|--------|----------------|---------|----------|--------|
| 1   | Q2487-09    | G4(0-6)        | SAM    | 07/07/25 08:15 |         | rubina   | OK     |
| 2   | Q2487-09DUP | G4(0-6)DUP     | DUP    | 07/07/25 08:23 |         | rubina   | OK     |
| 3   | Q2487-10    | G4(6-12)       | SAM    | 07/07/25 08:31 |         | rubina   | OK     |
| 4   | Q2487-11    | G3(0-6)        | SAM    | 07/07/25 08:38 |         | rubina   | OK     |
| 5   | Q2487-12    | G3(6-12)       | SAM    | 07/07/25 08:45 |         | rubina   | OK     |
| 6   | Q2487-13    | G2(0-6)        | SAM    | 07/07/25 08:52 |         | rubina   | OK     |
| 7   | Q2487-14    | G2(6-12)       | SAM    | 07/07/25 09:00 |         | rubina   | OK     |
| 8   | Q2487-15    | G1(0-6)        | SAM    | 07/07/25 09:07 |         | rubina   | OK     |
| 9   | Q2487-16    | G1(6-12)       | SAM    | 07/07/25 09:15 |         | rubina   | OK     |
| 10  | Q2491-01    | EO-1-070225    | SAM    | 07/07/25 09:22 |         | rubina   | OK     |
| 11  | Q2493-01    | WC-11          | SAM    | 07/07/25 09:30 |         | rubina   | OK     |
| 12  | Q2507-01    | SU-04-7.3-2025 | SAM    | 07/07/25 09:38 |         | rubina   | OK     |
| 13  | Q2513-01    | HR-2-070325    | SAM    | 07/07/25 09:45 |         | rubina   | OK     |
| 14  | Q2513-03    | HR-3-070325    | SAM    | 07/07/25 09:53 |         | rubina   | OK     |

**Instrument ID:** WC SC-3

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136383**

|                  |                  |              |                     |
|------------------|------------------|--------------|---------------------|
| Review By        | jignesh          | Review On    | 7/8/2025 2:58:29 PM |
| Supervise By     | Iwona            | Supervise On | 7/8/2025 3:06:16 PM |
| SubDirectory     | LB136383         | Test         | TVS                 |
| <b>STD. NAME</b> | <b>STD REF.#</b> |              |                     |
| ICAL Standard    | N/A              |              |                     |
| ICV Standard     | N/A              |              |                     |
| CCV Standard     | N/A              |              |                     |
| ICSA Standard    | N/A              |              |                     |
| CRI Standard     | N/A              |              |                     |
| LCS Standard     | N/A              |              |                     |
| Chk Standard     | N/A              |              |                     |

| Sr# | SampleId    | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|-------------|------------|--------|----------------|---------|----------|--------|
| 1   | LB136383BL  | LB136383BL | MB     | 07/07/25 16:00 |         | jignesh  | OK     |
| 2   | Q2487-09    | G4(0-6)    | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |
| 3   | Q2487-09DUP | G4(0-6)DUP | DUP    | 07/07/25 16:00 |         | jignesh  | OK     |
| 4   | Q2487-10    | G4(6-12)   | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |
| 5   | Q2487-11    | G3(0-6)    | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |
| 6   | Q2487-12    | G3(6-12)   | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |
| 7   | Q2487-13    | G2(0-6)    | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |
| 8   | Q2487-14    | G2(6-12)   | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |
| 9   | Q2487-15    | G1(0-6)    | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |
| 10  | Q2487-16    | G1(6-12)   | SAM    | 07/07/25 16:00 |         | jignesh  | OK     |

**Instrument ID:** KONELAB

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136387**

|                  |                                                                |              |                     |
|------------------|----------------------------------------------------------------|--------------|---------------------|
| Review By        | rubina                                                         | Review On    | 7/8/2025 1:13:21 PM |
| Supervise By     | Iwona                                                          | Supervise On | 7/8/2025 1:14:40 PM |
| SubDirectory     | LB136387                                                       | Test         | Cyanide             |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                               |              |                     |
| ICAL Standard    | WP113823,WP113824,WP113825,WP113826,WP113827,WP113828,WP113829 |              |                     |
| ICV Standard     | W3012                                                          |              |                     |
| CCV Standard     | WP113824                                                       |              |                     |
| ICSA Standard    | N/A                                                            |              |                     |
| CRI Standard     | N/A                                                            |              |                     |
| LCS Standard     | WP112995                                                       |              |                     |
| Chk Standard     | WP112643,WP112900,WP113831                                     |              |                     |

| Sr# | SampleId    | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|-------------|------------|--------|----------------|---------|----------|--------|
| 1   | 0.0PPBCN    | 0.0PPBCN   | CAL1   | 07/07/25 10:39 |         | rubina   | OK     |
| 2   | 5.0PPBCN    | 5.0PPBCN   | CAL2   | 07/07/25 10:39 |         | rubina   | OK     |
| 3   | 10PPBCN     | 10PPBCN    | CAL3   | 07/07/25 10:39 |         | rubina   | OK     |
| 4   | 50PPBCN     | 50PPBCN    | CAL4   | 07/07/25 10:39 |         | rubina   | OK     |
| 5   | 100PPBCN    | 100PPBCN   | CAL5   | 07/07/25 10:39 |         | rubina   | OK     |
| 6   | 250PPBCN    | 250PPBCN   | CAL6   | 07/07/25 10:39 |         | rubina   | OK     |
| 7   | 500PPBCN    | 500PPBCN   | CAL7   | 07/07/25 10:39 |         | rubina   | OK     |
| 8   | ICV1        | ICV1       | ICV    | 07/07/25 14:28 |         | rubina   | OK     |
| 9   | ICB1        | ICB1       | ICB    | 07/07/25 14:28 |         | rubina   | OK     |
| 10  | CCV1        | CCV1       | CCV    | 07/07/25 14:28 |         | rubina   | OK     |
| 11  | CCB1        | CCB1       | CCB    | 07/07/25 14:28 |         | rubina   | OK     |
| 12  | PB168749BL  | PB168749BL | MB     | 07/07/25 14:28 |         | rubina   | OK     |
| 13  | PB168749BS  | PB168749BS | LCS    | 07/07/25 14:35 |         | rubina   | OK     |
| 14  | Q2487-09    | G4(0-6)    | SAM    | 07/07/25 14:35 |         | rubina   | OK     |
| 15  | Q2487-09DUP | G4(0-6)DUP | DUP    | 07/07/25 14:35 |         | rubina   | OK     |
| 16  | Q2487-09MS  | G4(0-6)MS  | MS     | 07/07/25 14:35 |         | rubina   | OK     |
| 17  | Q2487-09MSD | G4(0-6)MSD | MSD    | 07/07/25 14:35 |         | rubina   | OK     |
| 18  | Q2487-10    | G4(6-12)   | SAM    | 07/07/25 14:35 |         | rubina   | OK     |

Instrument ID: KONELAB

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136387**

|                  |                                                                |              |                     |
|------------------|----------------------------------------------------------------|--------------|---------------------|
| Review By        | rubina                                                         | Review On    | 7/8/2025 1:13:21 PM |
| Supervise By     | Iwona                                                          | Supervise On | 7/8/2025 1:14:40 PM |
| SubDirectory     | LB136387                                                       | Test         | Cyanide             |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                               |              |                     |
| ICAL Standard    | WP113823,WP113824,WP113825,WP113826,WP113827,WP113828,WP113829 |              |                     |
| ICV Standard     | W3012                                                          |              |                     |
| CCV Standard     | WP113824                                                       |              |                     |
| ICSA Standard    | N/A                                                            |              |                     |
| CRI Standard     | N/A                                                            |              |                     |
| LCS Standard     | WP112995                                                       |              |                     |
| Chk Standard     | WP112643,WP112900,WP113831                                     |              |                     |

|    |          |          |     |                |  |        |    |
|----|----------|----------|-----|----------------|--|--------|----|
| 19 | Q2487-11 | G3(0-6)  | SAM | 07/07/25 14:35 |  | rubina | OK |
| 20 | Q2487-12 | G3(6-12) | SAM | 07/07/25 14:35 |  | rubina | OK |
| 21 | Q2487-13 | G2(0-6)  | SAM | 07/07/25 14:35 |  | rubina | OK |
| 22 | CCV2     | CCV2     | CCV | 07/07/25 14:35 |  | rubina | OK |
| 23 | CCB2     | CCB2     | CCB | 07/07/25 14:35 |  | rubina | OK |
| 24 | Q2487-14 | G2(6-12) | SAM | 07/07/25 14:41 |  | rubina | OK |
| 25 | Q2487-15 | G1(0-6)  | SAM | 07/07/25 14:41 |  | rubina | OK |
| 26 | Q2487-16 | G1(6-12) | SAM | 07/07/25 14:41 |  | rubina | OK |
| 27 | Q2515-01 | WC-1     | SAM | 07/07/25 14:41 |  | rubina | OK |
| 28 | CCV3     | CCV3     | CCV | 07/07/25 14:41 |  | rubina | OK |
| 29 | CCB3     | CCB3     | CCB | 07/07/25 14:42 |  | rubina | OK |

**Instrument ID:** WC SC-3

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136393**

|                  |                  |              |                     |
|------------------|------------------|--------------|---------------------|
| Review By        | jignesh          | Review On    | 7/8/2025 2:43:53 PM |
| Supervise By     | Iwona            | Supervise On | 7/8/2025 3:06:01 PM |
| SubDirectory     | LB136393         | Test         | TS                  |
| <b>STD. NAME</b> | <b>STD REF.#</b> |              |                     |
| ICAL Standard    | N/A              |              |                     |
| ICV Standard     | N/A              |              |                     |
| CCV Standard     | N/A              |              |                     |
| ICSA Standard    | N/A              |              |                     |
| CRI Standard     | N/A              |              |                     |
| LCS Standard     | N/A              |              |                     |
| Chk Standard     | N/A              |              |                     |

| Sr# | SampleId    | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|-------------|------------|--------|----------------|---------|----------|--------|
| 1   | LB136393BL  | LB136393BL | MB     | 07/07/25 11:00 |         | jignesh  | OK     |
| 2   | Q2487-09    | G4(0-6)    | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |
| 3   | Q2487-09DUP | G4(0-6)DUP | DUP    | 07/07/25 11:00 |         | jignesh  | OK     |
| 4   | Q2487-10    | G4(6-12)   | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |
| 5   | Q2487-11    | G3(0-6)    | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |
| 6   | Q2487-12    | G3(6-12)   | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |
| 7   | Q2487-13    | G2(0-6)    | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |
| 8   | Q2487-14    | G2(6-12)   | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |
| 9   | Q2487-15    | G1(0-6)    | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |
| 10  | Q2487-16    | G1(6-12)   | SAM    | 07/07/25 11:00 |         | jignesh  | OK     |

**Instrument ID:** SPECTROPHOTOMETER-1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136397**

|                  |                                     |              |                     |
|------------------|-------------------------------------|--------------|---------------------|
| Review By        | rubina                              | Review On    | 7/8/2025 4:24:26 PM |
| Supervise By     | Iwona                               | Supervise On | 7/9/2025 9:00:41 AM |
| SubDirectory     | LB136397                            | Test         | Hexavalent Chromium |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                     |
| ICAL Standard    | N/A                                 |              |                     |
| ICV Standard     | N/A                                 |              |                     |
| CCV Standard     | N/A                                 |              |                     |
| ICSA Standard    | N/A                                 |              |                     |
| CRI Standard     | N/A                                 |              |                     |
| LCS Standard     | N/A                                 |              |                     |
| Chk Standard     | WP113840,WP112831,WP112830,WP113087 |              |                     |

| Sr# | SampleId   | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|------------|------------|--------|----------------|---------|----------|--------|
| 1   | CAL1       | CAL1       | CAL    | 07/08/25 15:40 |         | rubina   | OK     |
| 2   | CAL2       | CAL2       | CAL    | 07/08/25 15:41 |         | rubina   | OK     |
| 3   | CAL3       | CAL3       | CAL    | 07/08/25 15:42 |         | rubina   | OK     |
| 4   | CAL4       | CAL4       | CAL    | 07/08/25 15:43 |         | rubina   | OK     |
| 5   | CAL5       | CAL5       | CAL    | 07/08/25 15:44 |         | rubina   | OK     |
| 6   | CAL6       | CAL6       | CAL    | 07/08/25 15:45 |         | rubina   | OK     |
| 7   | CAL7       | CAL7       | CAL    | 07/08/25 15:46 |         | rubina   | OK     |
| 8   | ICV        | ICV        | ICV    | 07/08/25 15:47 |         | rubina   | OK     |
| 9   | ICB        | ICB        | ICB    | 07/08/25 15:48 |         | rubina   | OK     |
| 10  | CCV1       | CCV1       | CCV    | 07/08/25 15:49 |         | rubina   | OK     |
| 11  | CCB1       | CCB1       | CCB    | 07/08/25 15:50 |         | rubina   | OK     |
| 12  | RL Check   | RL Check   | RL     | 07/08/25 15:51 |         | rubina   | OK     |
| 13  | PB168732BL | PB168732BL | MB     | 07/08/25 15:52 |         | rubina   | OK     |
| 14  | PB168732BS | PB168732BS | LCS    | 07/08/25 15:53 |         | rubina   | OK     |
| 15  | Q2487-09   | G4(0-6)    | SAM    | 07/08/25 15:54 |         | rubina   | OK     |
| 16  | Q2487-10   | G4(6-12)   | SAM    | 07/08/25 15:55 |         | rubina   | OK     |
| 17  | Q2487-11   | G3(0-6)    | SAM    | 07/08/25 15:56 |         | rubina   | OK     |
| 18  | Q2487-12   | G3(6-12)   | SAM    | 07/08/25 15:57 |         | rubina   | OK     |

Instrument ID: SPECTROPHOTOMETER-1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136397**

|                  |                                     |              |                     |
|------------------|-------------------------------------|--------------|---------------------|
| Review By        | rubina                              | Review On    | 7/8/2025 4:24:26 PM |
| Supervise By     | Iwona                               | Supervise On | 7/9/2025 9:00:41 AM |
| SubDirectory     | LB136397                            | Test         | Hexavalent Chromium |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                     |
| ICAL Standard    | N/A                                 |              |                     |
| ICV Standard     | N/A                                 |              |                     |
| CCV Standard     | N/A                                 |              |                     |
| ICSA Standard    | N/A                                 |              |                     |
| CRI Standard     | N/A                                 |              |                     |
| LCS Standard     | N/A                                 |              |                     |
| Chk Standard     | WP113840,WP112831,WP112830,WP113087 |              |                     |

|    |                 |          |     |                |  |        |    |
|----|-----------------|----------|-----|----------------|--|--------|----|
| 19 | Q2487-13        | G2(0-6)  | SAM | 07/08/25 15:58 |  | rubina | OK |
| 20 | Q2487-14        | G2(6-12) | SAM | 07/08/25 15:59 |  | rubina | OK |
| 21 | Q2487-15        | G1(0-6)  | SAM | 07/08/25 16:00 |  | rubina | OK |
| 22 | Q2487-16        | G1(6-12) | SAM | 07/08/25 16:01 |  | rubina | OK |
| 23 | CCV2            | CCV2     | CCV | 07/08/25 16:02 |  | rubina | OK |
| 24 | CCB2            | CCB2     | CCB | 07/08/25 16:03 |  | rubina | OK |
| 25 | Q2493-01        | WC-11    | SAM | 07/08/25 16:04 |  | rubina | OK |
| 26 | Q2517-01        | TP-14    | SAM | 07/08/25 16:05 |  | rubina | OK |
| 27 | Q2517-01DUP     | TP-14DUP | DUP | 07/08/25 16:06 |  | rubina | OK |
| 28 | Q2517-01MSPre   | TP-14MS  | MS  | 07/08/25 16:07 |  | rubina | OK |
| 29 | Q2517-01MS2Ins  | TP-14MS  | MS  | 07/08/25 16:08 |  | rubina | OK |
| 30 | Q2517-01MS3Post | TP-14MS  | MS  | 07/08/25 16:09 |  | rubina | OK |
| 31 | CCV3            | CCV3     | CCV | 07/08/25 16:10 |  | rubina | OK |
| 32 | CCB3            | CCB3     | CCB | 07/08/25 16:11 |  | rubina | OK |

**Instrument ID:** SPECTROPHOTOMETER-2

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136399**

|                  |                                                                                            |              |                     |
|------------------|--------------------------------------------------------------------------------------------|--------------|---------------------|
| Review By        | Iwona                                                                                      | Review On    | 7/8/2025 4:08:29 PM |
| Supervise By     | jignesh                                                                                    | Supervise On | 7/8/2025 4:09:31 PM |
| SubDirectory     | LB136399                                                                                   | Test         | COD                 |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                                           |              |                     |
| ICAL Standard    | N/A                                                                                        |              |                     |
| ICV Standard     | N/A                                                                                        |              |                     |
| CCV Standard     | N/A                                                                                        |              |                     |
| ICSA Standard    | N/A                                                                                        |              |                     |
| CRI Standard     | N/A                                                                                        |              |                     |
| LCS Standard     | N/A                                                                                        |              |                     |
| Chk Standard     | WP113238,WP113237,WP113235,WP113234,WP113233,WP113240,WP113236,W3129,WP113843,WP113844,WP1 |              |                     |

| Sr# | SampleId    | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|-------------|------------|--------|----------------|---------|----------|--------|
| 1   | CAL1        | CAL1       | CAL    | 05/28/25 13:10 |         | Iwona    | OK     |
| 2   | CAL2        | CAL2       | CAL    | 05/28/25 13:10 |         | Iwona    | OK     |
| 3   | CAL3        | CAL3       | CAL    | 05/28/25 13:11 |         | Iwona    | OK     |
| 4   | CAL4        | CAL4       | CAL    | 05/28/25 13:11 |         | Iwona    | OK     |
| 5   | CAL5        | CAL5       | CAL    | 05/28/25 13:12 |         | Iwona    | OK     |
| 6   | CAL6        | CAL6       | CAL    | 05/28/25 13:12 |         | Iwona    | OK     |
| 7   | ICV         | ICV        | ICV    | 05/28/25 13:13 |         | Iwona    | OK     |
| 8   | ICB         | ICB        | ICB    | 05/28/25 13:13 |         | Iwona    | OK     |
| 9   | CCV1        | CCV1       | CCV    | 07/08/25 15:00 |         | Iwona    | OK     |
| 10  | CCB1        | CCB1       | CCB    | 07/08/25 15:00 |         | Iwona    | OK     |
| 11  | RL Check    | RL Check   | RL     | 07/08/25 15:01 |         | Iwona    | OK     |
| 12  | LB136399BL  | LB136399BL | MB     | 07/08/25 15:01 |         | Iwona    | OK     |
| 13  | LB136399BS  | LB136399BS | LCS    | 07/08/25 15:02 |         | Iwona    | OK     |
| 14  | Q2487-09    | G4(0-6)    | SAM    | 07/08/25 15:02 |         | Iwona    | OK     |
| 15  | Q2487-09DUP | G4(0-6)DUP | DUP    | 07/08/25 15:03 |         | Iwona    | OK     |
| 16  | Q2487-09MS  | G4(0-6)MS  | MS     | 07/08/25 15:03 |         | Iwona    | OK     |
| 17  | Q2487-09MSD | G4(0-6)MSD | MSD    | 07/08/25 15:04 |         | Iwona    | OK     |
| 18  | Q2487-10    | G4(6-12)   | SAM    | 07/08/25 15:04 |         | Iwona    | OK     |

**Instrument ID:** SPECTROPHOTOMETER-2

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136399**

|                  |                                                                                            |              |                     |
|------------------|--------------------------------------------------------------------------------------------|--------------|---------------------|
| Review By        | Iwona                                                                                      | Review On    | 7/8/2025 4:08:29 PM |
| Supervise By     | jignesh                                                                                    | Supervise On | 7/8/2025 4:09:31 PM |
| SubDirectory     | LB136399                                                                                   | Test         | COD                 |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                                           |              |                     |
| ICAL Standard    | N/A                                                                                        |              |                     |
| ICV Standard     | N/A                                                                                        |              |                     |
| CCV Standard     | N/A                                                                                        |              |                     |
| ICSA Standard    | N/A                                                                                        |              |                     |
| CRI Standard     | N/A                                                                                        |              |                     |
| LCS Standard     | N/A                                                                                        |              |                     |
| Chk Standard     | WP113238,WP113237,WP113235,WP113234,WP113233,WP113240,WP113236,W3129,WP113843,WP113844,WP1 |              |                     |

|    |          |          |     |                |  |       |    |
|----|----------|----------|-----|----------------|--|-------|----|
| 19 | Q2487-11 | G3(0-6)  | SAM | 07/08/25 15:05 |  | Iwona | OK |
| 20 | Q2487-12 | G3(6-12) | SAM | 07/08/25 15:05 |  | Iwona | OK |
| 21 | Q2487-13 | G2(0-6)  | SAM | 07/08/25 15:06 |  | Iwona | OK |
| 22 | CCV2     | CCV2     | CCV | 07/08/25 15:06 |  | Iwona | OK |
| 23 | CCB2     | CCB2     | CCB | 07/08/25 15:07 |  | Iwona | OK |
| 24 | Q2487-14 | G2(6-12) | SAM | 07/08/25 15:07 |  | Iwona | OK |
| 25 | Q2487-15 | G1(0-6)  | SAM | 07/08/25 15:08 |  | Iwona | OK |
| 26 | Q2487-16 | G1(6-12) | SAM | 07/08/25 15:08 |  | Iwona | OK |
| 27 | CCV3     | CCV3     | CCV | 07/08/25 15:09 |  | Iwona | OK |
| 28 | CCB3     | CCB3     | CCB | 07/08/25 15:09 |  | Iwona | OK |

**Instrument ID:** KONELAB

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136405**

|                  |                                     |              |                       |
|------------------|-------------------------------------|--------------|-----------------------|
| Review By        | rubina                              | Review On    | 7/10/2025 12:27:55 PM |
| Supervise By     | Iwona                               | Supervise On | 7/10/2025 12:57:25 PM |
| SubDirectory     | LB136405                            | Test         | Ammonia               |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                       |
| ICAL Standard    | WP113849                            |              |                       |
| ICV Standard     | WP113850                            |              |                       |
| CCV Standard     | WP113850                            |              |                       |
| ICSA Standard    | N/A                                 |              |                       |
| CRI Standard     | N/A                                 |              |                       |
| LCS Standard     | WP113449                            |              |                       |
| Chk Standard     | WP113852,WP111745,WP111385,WP111660 |              |                       |

| Sr# | SampleId    | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|-------------|------------|--------|----------------|---------|----------|--------|
| 1   | 0.0PPM      | 0.0PPM     | CAL1   | 07/09/25 09:20 |         | rubina   | OK     |
| 2   | 0.1PPM      | 0.1PPM     | CAL2   | 07/09/25 09:20 |         | rubina   | OK     |
| 3   | 0.2PPM      | 0.2PPM     | CAL3   | 07/09/25 09:20 |         | rubina   | OK     |
| 4   | 0.4PPM      | 0.4PPM     | CAL4   | 07/09/25 09:21 |         | rubina   | OK     |
| 5   | 1.0PPM      | 1.0PPM     | CAL5   | 07/09/25 09:21 |         | rubina   | OK     |
| 6   | 1.3PPM      | 1.3PPM     | CAL6   | 07/09/25 09:21 |         | rubina   | OK     |
| 7   | 2.0PPM      | 2.0PPM     | CAL7   | 07/09/25 09:21 |         | rubina   | OK     |
| 8   | ICV1        | ICV1       | ICV    | 07/09/25 10:03 |         | rubina   | OK     |
| 9   | ICB1        | ICB1       | ICB    | 07/09/25 10:03 |         | rubina   | OK     |
| 10  | CCV1        | CCV1       | CCV    | 07/09/25 10:03 |         | rubina   | OK     |
| 11  | CCB1        | CCB1       | CCB    | 07/09/25 10:03 |         | rubina   | OK     |
| 12  | RL          | RL         | SAM    | 07/09/25 10:03 |         | rubina   | OK     |
| 13  | PB168757BL  | PB168757BL | MB     | 07/09/25 10:14 |         | rubina   | OK     |
| 14  | PB168757BS  | PB168757BS | LCS    | 07/09/25 10:14 |         | rubina   | OK     |
| 15  | Q2487-09    | G4(0-6)    | SAM    | 07/09/25 10:14 |         | rubina   | OK     |
| 16  | Q2487-09DUP | G4(0-6)DUP | DUP    | 07/09/25 10:14 |         | rubina   | OK     |
| 17  | Q2487-09MS  | G4(0-6)MS  | MS     | 07/09/25 10:14 |         | rubina   | OK     |
| 18  | Q2487-09MSD | G4(0-6)MSD | MSD    | 07/09/25 10:14 |         | rubina   | OK     |

Instrument ID: KONELAB

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136405**

|                  |                                     |              |                       |
|------------------|-------------------------------------|--------------|-----------------------|
| Review By        | rubina                              | Review On    | 7/10/2025 12:27:55 PM |
| Supervise By     | Iwona                               | Supervise On | 7/10/2025 12:57:25 PM |
| SubDirectory     | LB136405                            | Test         | Ammonia               |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                       |
| ICAL Standard    | WP113849                            |              |                       |
| ICV Standard     | WP113850                            |              |                       |
| CCV Standard     | WP113850                            |              |                       |
| ICSA Standard    | N/A                                 |              |                       |
| CRI Standard     | N/A                                 |              |                       |
| LCS Standard     | WP113449                            |              |                       |
| Chk Standard     | WP113852,WP111745,WP111385,WP111660 |              |                       |

|    |             |                 |     |                |      |        |          |
|----|-------------|-----------------|-----|----------------|------|--------|----------|
| 19 | Q2487-10    | G4(6-12)        | SAM | 07/09/25 10:14 |      | rubina | OK       |
| 20 | Q2487-11    | G3(0-6)         | SAM | 07/09/25 10:14 |      | rubina | OK       |
| 21 | Q2487-12    | G3(6-12)        | SAM | 07/09/25 10:25 |      | rubina | OK       |
| 22 | CCV2        | CCV2            | CCV | 07/09/25 10:25 |      | rubina | OK       |
| 23 | CCB2        | CCB2            | CCB | 07/09/25 10:25 |      | rubina | OK       |
| 24 | Q2487-13    | G2(0-6)         | SAM | 07/09/25 10:25 |      | rubina | OK       |
| 25 | Q2487-14    | G2(6-12)        | SAM | 07/09/25 10:25 |      | rubina | OK       |
| 26 | Q2487-15    | G1(0-6)         | SAM | 07/09/25 10:25 |      | rubina | OK       |
| 27 | Q2487-16    | G1(6-12)        | SAM | 07/09/25 10:25 |      | rubina | OK       |
| 28 | PB168758BL  | PB168758BL      | MB  | 07/09/25 10:25 |      | rubina | OK       |
| 29 | PB168758BS  | PB168758BS      | LCS | 07/09/25 10:36 |      | rubina | OK       |
| 30 | Q2512-01    | WATER TREATMENT | SAM | 07/09/25 10:36 | High | rubina | Dilution |
| 31 | Q2512-01DUP | WATER TREATMENT | DUP | 07/09/25 10:36 | High | rubina | Dilution |
| 32 | Q2512-01MS  | WATER TREATMENT | MS  | 07/09/25 10:36 |      | rubina | OK       |
| 33 | Q2512-01MSD | WATER TREATMENT | MSD | 07/09/25 10:36 |      | rubina | OK       |
| 34 | CCV3        | CCV3            | CCV | 07/09/25 10:36 |      | rubina | OK       |
| 35 | CCB3        | CCB3            | CCB | 07/09/25 10:36 |      | rubina | OK       |
| 36 | Q2525-01    | EFFLUENT-COMPOS | SAM | 07/09/25 10:43 |      | rubina | OK       |
| 37 | CCV4        | CCV4            | CCV | 07/09/25 10:43 |      | rubina | OK       |
| 38 | CCB4        | CCB4            | CCB | 07/09/25 10:43 |      | rubina | OK       |

Instrument ID: KONELAB

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136405**

|                  |                                     |              |                       |
|------------------|-------------------------------------|--------------|-----------------------|
| Review By        | rubina                              | Review On    | 7/10/2025 12:27:55 PM |
| Supervise By     | Iwona                               | Supervise On | 7/10/2025 12:57:25 PM |
| SubDirectory     | LB136405                            | Test         | Ammonia               |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                       |
| ICAL Standard    | WP113849                            |              |                       |
| ICV Standard     | WP113850                            |              |                       |
| CCV Standard     | WP113850                            |              |                       |
| ICSA Standard    | N/A                                 |              |                       |
| CRI Standard     | N/A                                 |              |                       |
| LCS Standard     | WP113449                            |              |                       |
| Chk Standard     | WP113852,WP111745,WP111385,WP111660 |              |                       |

|    |               |                 |     |                |           |        |          |
|----|---------------|-----------------|-----|----------------|-----------|--------|----------|
| 39 | Q2512-01DL    | WATER TREATMENT | SAM | 07/09/25 11:04 | Report 2X | rubina | Confirms |
| 40 | Q2512-01DUPDL | WATER TREATMENT | DUP | 07/09/25 11:04 | Report 2X | rubina | Confirms |
| 41 | CCV5          | CCV5            | CCV | 07/09/25 11:04 |           | rubina | OK       |
| 42 | CCB5          | CCB5            | CCB | 07/09/25 11:04 |           | rubina | OK       |

**Instrument ID:** KONELAB

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136408**

|                  |                                     |              |                     |
|------------------|-------------------------------------|--------------|---------------------|
| Review By        | rubina                              | Review On    | 7/9/2025 3:52:38 PM |
| Supervise By     | Iwona                               | Supervise On | 7/9/2025 4:17:54 PM |
| SubDirectory     | LB136408                            | Test         | TKN                 |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                     |
| ICAL Standard    | WP113832                            |              |                     |
| ICV Standard     | WP113834                            |              |                     |
| CCV Standard     | WP113833                            |              |                     |
| ICSA Standard    | N/A                                 |              |                     |
| CRI Standard     | N/A                                 |              |                     |
| LCS Standard     | WP113835                            |              |                     |
| Chk Standard     | WP113852,WP111745,WP111385,WP111660 |              |                     |

| Sr# | SampleId    | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|-------------|------------|--------|----------------|---------|----------|--------|
| 1   | 0.0PPM      | 0.0PPM     | CAL1   | 07/09/25 11:35 |         | rubina   | OK     |
| 2   | 0.5PPM      | 0.5PPM     | CAL2   | 07/09/25 11:35 |         | rubina   | OK     |
| 3   | 1.0PPM      | 1.0PPM     | CAL3   | 07/09/25 11:35 |         | rubina   | OK     |
| 4   | 2.5PPM      | 2.5PPM     | CAL4   | 07/09/25 11:35 |         | rubina   | OK     |
| 5   | 5.0PPM      | 5.0PPM     | CAL5   | 07/09/25 11:35 |         | rubina   | OK     |
| 6   | 6.7PPM      | 6.7PPM     | CAL6   | 07/09/25 11:35 |         | rubina   | OK     |
| 7   | 10.0PPM     | 10.0PPM    | CAL7   | 07/09/25 11:35 |         | rubina   | OK     |
| 8   | ICV1        | ICV1       | ICV    | 07/09/25 12:18 |         | rubina   | OK     |
| 9   | ICB1        | ICB1       | ICB    | 07/09/25 12:18 |         | rubina   | OK     |
| 10  | CCV1        | CCV1       | CCV    | 07/09/25 12:18 |         | rubina   | OK     |
| 11  | CCB1        | CCB1       | CCB    | 07/09/25 12:18 |         | rubina   | OK     |
| 12  | RL          | RL         | LOQ    | 07/09/25 12:18 |         | rubina   | OK     |
| 13  | PB168731BL  | PB168731BL | MB     | 07/09/25 12:18 |         | rubina   | OK     |
| 14  | PB168731BS  | PB168731BS | LCS    | 07/09/25 12:29 |         | rubina   | OK     |
| 15  | Q2487-09    | G4(0-6)    | SAM    | 07/09/25 12:29 |         | rubina   | OK     |
| 16  | Q2487-09DUP | G4(0-6)DUP | DUP    | 07/09/25 12:29 |         | rubina   | OK     |
| 17  | Q2487-09MS  | G4(0-6)MS  | MS     | 07/09/25 12:29 |         | rubina   | OK     |
| 18  | Q2487-09MSD | G4(0-6)MSD | MSD    | 07/09/25 12:29 |         | rubina   | OK     |

Instrument ID: KONELAB

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136408**

|                  |                                     |              |                     |
|------------------|-------------------------------------|--------------|---------------------|
| Review By        | rubina                              | Review On    | 7/9/2025 3:52:38 PM |
| Supervise By     | Iwona                               | Supervise On | 7/9/2025 4:17:54 PM |
| SubDirectory     | LB136408                            | Test         | TKN                 |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                     |
| ICAL Standard    | WP113832                            |              |                     |
| ICV Standard     | WP113834                            |              |                     |
| CCV Standard     | WP113833                            |              |                     |
| ICSA Standard    | N/A                                 |              |                     |
| CRI Standard     | N/A                                 |              |                     |
| LCS Standard     | WP113835                            |              |                     |
| Chk Standard     | WP113852,WP111745,WP111385,WP111660 |              |                     |

|    |             |            |     |                |             |        |          |
|----|-------------|------------|-----|----------------|-------------|--------|----------|
| 19 | Q2487-10    | G4(6-12)   | SAM | 07/09/25 12:29 | TKN is high | rubina | Dilution |
| 20 | Q2487-11    | G3(0-6)    | SAM | 07/09/25 12:29 |             | rubina | OK       |
| 21 | Q2487-12    | G3(6-12)   | SAM | 07/09/25 12:29 | TKN is high | rubina | Dilution |
| 22 | CCV2        | CCV2       | CCV | 07/09/25 12:39 |             | rubina | OK       |
| 23 | CCB2        | CCB2       | CCB | 07/09/25 12:39 |             | rubina | OK       |
| 24 | Q2487-13    | G2(0-6)    | SAM | 07/09/25 12:40 |             | rubina | OK       |
| 25 | Q2487-14    | G2(6-12)   | SAM | 07/09/25 12:40 | TKN is high | rubina | Dilution |
| 26 | Q2487-15    | G1(0-6)    | SAM | 07/09/25 12:40 |             | rubina | OK       |
| 27 | Q2487-16    | G1(6-12)   | SAM | 07/09/25 12:40 |             | rubina | OK       |
| 28 | PB168699BL  | PB168699BL | MB  | 07/09/25 12:40 |             | rubina | OK       |
| 29 | PB168699BS  | PB168699BS | LCS | 07/09/25 12:50 |             | rubina | OK       |
| 30 | Q2470-01    | SW-1       | SAM | 07/09/25 12:50 |             | rubina | OK       |
| 31 | Q2470-01DUP | SW-1DUP    | DUP | 07/09/25 12:50 |             | rubina | OK       |
| 32 | Q2470-01MS  | SW-1MS     | MS  | 07/09/25 12:50 |             | rubina | OK       |
| 33 | Q2470-01MSD | SW-1MSD    | MSD | 07/09/25 12:50 |             | rubina | OK       |
| 34 | CCV3        | CCV3       | CCV | 07/09/25 12:50 |             | rubina | OK       |
| 35 | CCB3        | CCB3       | CCB | 07/09/25 12:50 |             | rubina | OK       |
| 36 | Q2471-01    | SW-1       | SAM | 07/09/25 12:50 |             | rubina | OK       |
| 37 | CCV4        | CCV4       | CCV | 07/09/25 12:57 |             | rubina | OK       |
| 38 | CCB4        | CCB4       | CCB | 07/09/25 12:57 |             | rubina | OK       |

Instrument ID: KONELAB

**Daily Analysis Runlog For Sequence/QCBatch ID # LB136408**

|                  |                                     |              |                     |
|------------------|-------------------------------------|--------------|---------------------|
| Review By        | rubina                              | Review On    | 7/9/2025 3:52:38 PM |
| Supervise By     | Iwona                               | Supervise On | 7/9/2025 4:17:54 PM |
| SubDirectory     | LB136408                            | Test         | TKN                 |
| <b>STD. NAME</b> | <b>STD REF.#</b>                    |              |                     |
| ICAL Standard    | WP113832                            |              |                     |
| ICV Standard     | WP113834                            |              |                     |
| CCV Standard     | WP113833                            |              |                     |
| ICSA Standard    | N/A                                 |              |                     |
| CRI Standard     | N/A                                 |              |                     |
| LCS Standard     | WP113835                            |              |                     |
| Chk Standard     | WP113852,WP111745,WP111385,WP111660 |              |                     |

|    |            |            |     |                |             |        |          |
|----|------------|------------|-----|----------------|-------------|--------|----------|
| 39 | Q2487-10DL | G4(6-12)DL | SAM | 07/09/25 13:32 | 10X For TKN | rubina | Confirms |
| 40 | Q2487-12DL | G3(6-12)DL | SAM | 07/09/25 13:32 | 5X For TKN  | rubina | Confirms |
| 41 | Q2487-14DL | G2(6-12)DL | SAM | 07/09/25 13:32 | 2X For TKN  | rubina | Confirms |
| 42 | CCV5       | CCV5       | CCV | 07/09/25 13:32 |             | rubina | OK       |
| 43 | CCB5       | CCB5       | CCB | 07/09/25 13:32 |             | rubina | OK       |

**Instrument ID:** WC SC-3

**Daily Analysis Runlog For Sequence/QC Batch ID # LB136422**

|                  |                                                        |              |                      |
|------------------|--------------------------------------------------------|--------------|----------------------|
| Review By        | jignesh                                                | Review On    | 7/10/2025 1:15:14 PM |
| Supervise By     | Iwona                                                  | Supervise On | 7/11/2025 1:04:54 PM |
| SubDirectory     | LB136422                                               | Test         | Oil and Grease       |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                       |              |                      |
| ICAL Standard    | N/A                                                    |              |                      |
| ICV Standard     | N/A                                                    |              |                      |
| CCV Standard     | N/A                                                    |              |                      |
| ICSA Standard    | N/A                                                    |              |                      |
| CRI Standard     | N/A                                                    |              |                      |
| LCS Standard     | N/A                                                    |              |                      |
| Chk Standard     | W3204,NA,EP2624,WP112782,NA,W3186,WP112785,NA,WP112786 |              |                      |

| Sr# | SampleId    | ClientID   | QcType | Date           | Comment | Operator | Status |
|-----|-------------|------------|--------|----------------|---------|----------|--------|
| 1   | LB136422BL  | LB136422BL | MB     | 07/10/25 14:30 |         | jignesh  | OK     |
| 2   | LB136422BS  | LB136422BS | LCS    | 07/10/25 14:30 |         | jignesh  | OK     |
| 3   | Q2487-09    | G4(0-6)    | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |
| 4   | Q2487-09DUP | G4(0-6)DUP | DUP    | 07/10/25 14:30 |         | jignesh  | OK     |
| 5   | Q2487-09MS  | G4(0-6)MS  | MS     | 07/10/25 14:30 |         | jignesh  | OK     |
| 6   | Q2487-09MSD | G4(0-6)MSD | MSD    | 07/10/25 14:30 |         | jignesh  | OK     |
| 7   | Q2487-10    | G4(6-12)   | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |
| 8   | Q2487-11    | G3(0-6)    | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |
| 9   | Q2487-12    | G3(6-12)   | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |
| 10  | Q2487-13    | G2(0-6)    | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |
| 11  | Q2487-14    | G2(6-12)   | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |
| 12  | Q2487-15    | G1(0-6)    | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |
| 13  | Q2487-16    | G1(6-12)   | SAM    | 07/10/25 14:30 |         | jignesh  | OK     |

**Instrument ID:** GRAVIMETRIC

**Daily Analysis Runlog For Sequence/QCBatch ID # LB136506**

|                  |                  |              |                      |
|------------------|------------------|--------------|----------------------|
| Review By        | Ayul             | Review On    | 7/16/2025 2:58:02 PM |
| Supervise By     | amarnath         | Supervise On | 7/16/2025 2:58:55 PM |
| SubDirectory     | LB136506         | Test         | Trivalent Chromium   |
| <b>STD. NAME</b> | <b>STD REF.#</b> |              |                      |
| ICAL Standard    | N/A              |              |                      |
| ICV Standard     | N/A              |              |                      |
| CCV Standard     | N/A              |              |                      |
| ICSA Standard    | N/A              |              |                      |
| CRI Standard     | N/A              |              |                      |
| LCS Standard     | N/A              |              |                      |
| Chk Standard     | N/A              |              |                      |

| Sr# | SampleId | ClientID | QcType | Date           | Comment | Operator | Status |
|-----|----------|----------|--------|----------------|---------|----------|--------|
| 1   | Q2487-15 | G1(0-6)  | SAM    | 07/09/25 17:14 |         | Ayul     | OK     |
| 2   | Q2487-09 | G4(0-6)  | SAM    | 07/09/25 18:12 |         | Ayul     | OK     |
| 3   | Q2487-10 | G4(6-12) | SAM    | 07/09/25 18:17 |         | Ayul     | OK     |
| 4   | Q2487-11 | G3(0-6)  | SAM    | 07/09/25 18:21 |         | Ayul     | OK     |
| 5   | Q2487-12 | G3(6-12) | SAM    | 07/09/25 18:26 |         | Ayul     | OK     |
| 6   | Q2487-13 | G2(0-6)  | SAM    | 07/09/25 18:31 |         | Ayul     | OK     |
| 7   | Q2487-16 | G1(6-12) | SAM    | 07/09/25 19:03 |         | Ayul     | OK     |
| 8   | Q2487-14 | G2(6-12) | SAM    | 07/09/25 21:03 |         | Ayul     | OK     |