

SOP ID :	M1311-TCLP-16	
SDG No :	N/A	Start Prep Date : 07/02/2025 Time : 15:00
Weigh By :	JP	End Prep Date : 07/03/2025 Time : 08:20
Balance ID :	WC SC-7	Combination Ratio : 20
pH Meter ID :	WC PH METER-1	ZHE Cleaning Batch : N/A
Extraction By :	JP	Initial Room Temperature: 24 °C
Filter By :	JP	Final Room Temperature: 22 °C
Pipette ID :	WC	TCLP Technician Signature : <i>JP</i>
Tumbler ID :	ZHE-1	Supervisor By : <i>12</i>
TCLP Filter ID :	50223706	

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
TCLP-FLUID-1	N/A	WP112795
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
40ml VOA Vials	430992	N/A

Extraction Conformance/Non-Conformance Comments:

TUMBLER ZHE-1 checked, 30 rpm. ALL ZHE samples are extracted and given as vial A & B. Leak checked after 10 minutes of tumbling.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
07/03/25 11:00	<i>JP</i> 1 Cell Room	<i>JP</i> 1 VDC Lab
	Preparation Group	Analysis Group

Sample ID	ClientID	ZHE Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB168690TB	LEB690	11	N/A	500	N/A	N/A	N/A	4.94	N/A	ZHE-2
Q2478-04	WC-1	01	25.02	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-17	G4(1.5)	02	25.03	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-18	G4(10)	03	25.02	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-19	G3(9)	04	25.03	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-20	G3(3)	05	25.02	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-21	G2(2.5)	06	25.03	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-22	G2(9)	07	25.02	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-23	G1(10)	08	25.03	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2487-24	G1(4.5)	09	25.02	500	N/A	N/A	N/A	N/A	N/A	ZHE-1
Q2493-04	WC-11	10	25.03	500	N/A	N/A	N/A	N/A	N/A	ZHE-1

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB168690TB	LEB690	N/A	N/A	N/A	N/A	N/A	N/A
Q2478-04	WC-1	N/A	N/A	N/A	N/A	100	N/A
Q2487-17	G4(1.5)	N/A	N/A	N/A	N/A	100	N/A
Q2487-18	G4(10)	N/A	N/A	N/A	N/A	100	N/A
Q2487-19	G3(9)	N/A	N/A	N/A	N/A	100	N/A
Q2487-20	G3(3)	N/A	N/A	N/A	N/A	100	N/A
Q2487-21	G2(2.5)	N/A	N/A	N/A	N/A	100	N/A
Q2487-22	G2(9)	N/A	N/A	N/A	N/A	100	N/A
Q2487-23	G1(10)	N/A	N/A	N/A	N/A	100	N/A
Q2487-24	G1(4.5)	N/A	N/A	N/A	N/A	100	N/A
Q2493-04	WC-11	N/A	N/A	N/A	N/A	100	N/A