

SOP ID : M1311-TCLP-16

SDG No : N/A

Weigh By : JP

Balance ID : WC SC-7

pH Meter ID : WC PH METER-1

Extraction By : JP

Filter By : JP

Pipette ID : WC

Tumbler ID : T-1

TCLP Filter ID : 115525

Start Prep Date : 07/03/2025 Time : 15:00

End Prep Date : 07/04/2025 Time : 08:20

Combination Ratio : 20

ZHE Cleaning Batch : ~~N/A~~ 10

Initial Room Temperature: 24 °C

Final Room Temperature: 22 °C

TCLP Technician Signature : 


Supervisor By : 

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
TCLP-FLUID-1	N/A	WP112795
HCL-TCLP,1N	N/A	WP112797
HNO3-TCLP,1N	N/A	WP112799
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	W1940,W1941,W1942	W3166,W1938,W1939,
1 Liter Amber	N/A	90924-08
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. TUMBLER T-1 checked,30 rpm. q2489-08 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
07/10/25 11:00	 Tech Room	SKG RJ 164T
	Preparation Group	Analysis Group 

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB168728TB	LEB728	09	N/A	2000	N/A	N/A	N/A	4.93	1.0	T-1
Q2489-01	G4(0-6)	01	100.02	2000	N/A	N/A	N/A	6.0	1.0	T-1
Q2489-02	G4(6-12)	02	100.03	2000	N/A	N/A	N/A	6.2	1.5	T-1
Q2489-03	G3(0-6)	03	100.02	2000	N/A	N/A	N/A	6.0	1.0	T-1
Q2489-04	G3(6-12)	04	100.03	2000	N/A	N/A	N/A	6.2	1.5	T-1
Q2489-05	G2(0-6)	05	100.04	2000	N/A	N/A	N/A	7.6	1.0	T-1
Q2489-06	G2(6-12)	06	100.03	2000	N/A	N/A	N/A	6.0	1.5	T-1
Q2489-07	G1(0-6)	07	100.02	2000	N/A	N/A	N/A	7.6	1.0	T-1
Q2489-08	G1(6-12)	08	100.01	2000	N/A	N/A	N/A	6.0	1.5	T-1

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB168728TB	LEB728	N/A	N/A	N/A	N/A	N/A	N/A
Q2489-01	G4(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-02	G4(6-12)	N/A	N/A	N/A	N/A	100	N/A
Q2489-03	G3(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-04	G3(6-12)	N/A	N/A	N/A	N/A	100	N/A
Q2489-05	G2(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-06	G2(6-12)	N/A	N/A	N/A	N/A	100	N/A
Q2489-07	G1(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-08	G1(6-12)	N/A	N/A	N/A	N/A	100	N/A

Hot Block ID : WC S-1 / WC S-2

Thermometer ID : FLASHPOINT

SampleID	ClientID	Sample Weight (g)	Volume DI Water (mL)	PH after 5 min stir	PH after 10 min stir	Extraction Fluid 1 or 2	pH Extraction Fluid
PB168728TB	LEB728	N/A	N/A	N/A	N/A	#1	4.93
Q2489-01	G4(0-6)	5.02	96.5	7.6	2.5	#1	4.93
Q2489-02	G4(6-12)	5.01	96.5	8.0	3.0	#1	4.93
Q2489-03	G3(0-6)	5.02	96.5	8.0	3.0	#1	4.93
Q2489-04	G3(6-12)	5.03	96.5	8.0	3.0	#1	4.93
Q2489-05	G2(0-6)	5.04	96.5	9.0	3.5	#1	4.93
Q2489-06	G2(6-12)	5.03	96.5	8.0	3.0	#1	4.93
Q2489-07	G1(0-6)	5.02	96.5	9.2	3.5	#1	4.93
Q2489-08	G1(6-12)	5.01	96.5	8.6	3.0	#1	4.93