

SOP ID : M1312-SPLP-10
SDG No : N/A
Weigh By : JP
Balance ID : WC SC-7
pH Meter ID : WC PH METER-1
Extraction By : JP
Filter By : JP
Pipette ID : WC
Tumbler ID : T-2
TCLP Filter ID : 115525

Start Prep Date : 07/08/2025 Time : 14:00
End Prep Date : 07/09/2025 Time : 07:15
Combination Ratio : 20
ZHE Cleaning Batch : N/A
Initial Room Temperature: 23 °C
Final Room Temperature: 22 °C
TCLP Technician Signature : *[Signature]*
Supervisor By : *[Signature]*

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
SPLP FLUID	WP112802	N/A
N/A	N/A	N/A
HNO3-TCLP,1N	N/A	WP112799
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	N/A	W3166,W1938,W1939,
1 Liter Amber	N/A	90924-08
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Tumbler T-2 checked,30 rpm. Matrix spikes are added after filtration and before preservation. Particle size reduction is not required.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
07/10/25 11:00	<i>[Signature]</i> HCA Room	<i>[Signature]</i> RS 1 EAT
	Preparation Group	Analysis Group <i>[Signature]</i>

SPLP EXTRACTION LOGPAGE

PB168755

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB168755TB	LEB755	19	N/A	2000	N/A	N/A	N/A	4.24	1.0	T-2
Q2489-01	G4(0-6)	11	100.02	2000	N/A	N/A	N/A	7.0	1.0	T-2
Q2489-02	G4(6-12)	12	100.03	2000	N/A	N/A	N/A	6.0	1.5	T-2
Q2489-03	G3(0-6)	13	100.02	2000	N/A	N/A	N/A	7.0	1.5	T-2
Q2489-04	G3(6-12)	14	100.03	2000	N/A	N/A	N/A	5.6	1.0	T-2
Q2489-05	G2(0-6)	15	100.04	2000	N/A	N/A	N/A	7.0	1.0	T-2
Q2489-06	G2(6-12)	16	100.02	2000	N/A	N/A	N/A	7.2	1.5	T-2
Q2489-07	G1(0-6)	17	100.01	2000	N/A	N/A	N/A	7.0	1.0	T-2
Q2489-08	G1(6-12)	18	100.02	2000	N/A	N/A	N/A	7.2	1.5	T-2

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB168755TB	LEB755	N/A	N/A	N/A	N/A	N/A	N/A
Q2489-01	G4(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-02	G4(6-12)	N/A	N/A	N/A	N/A	100	N/A
Q2489-03	G3(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-04	G3(6-12)	N/A	N/A	N/A	N/A	100	N/A
Q2489-05	G2(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-06	G2(6-12)	N/A	N/A	N/A	N/A	100	N/A
Q2489-07	G1(0-6)	N/A	N/A	N/A	N/A	100	N/A
Q2489-08	G1(6-12)	N/A	N/A	N/A	N/A	100	N/A

Hot Block ID : WC S-1 / WC S-2

Thermometer ID : FLASHPOINT

SampleID	ClientID	Sample Weight (g)	Volume DI Water (mL)	PH after 5 min stir	PH after 10 min stir	Extraction Fluid 1 or 2	pH Extraction Fluid
PB168755TB	LEB755	N/A	N/A	N/A	N/A	#1	4.24
Q2489-01	G4(0-6)	N/A	N/A	N/A	N/A	#1	4.24
Q2489-02	G4(6-12)	N/A	N/A	N/A	N/A	#1	4.24
Q2489-03	G3(0-6)	N/A	N/A	N/A	N/A	#1	4.24
Q2489-04	G3(6-12)	N/A	N/A	N/A	N/A	#1	4.24
Q2489-05	G2(0-6)	N/A	N/A	N/A	N/A	#1	4.24
Q2489-06	G2(6-12)	N/A	N/A	N/A	N/A	#1	4.24
Q2489-07	G1(0-6)	N/A	N/A	N/A	N/A	#1	4.24
Q2489-08	G1(6-12)	N/A	N/A	N/A	N/A	#1	4.24