

Surrogate Summary

SDG No.: Q2529

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS031005.D	PIBLK-PS031005.D	2,4-DCAA	1	500	356	71		61	136
		2,4-DCAA	2	500	496	99		61	136
I.BLK-PS031021.D	PIBLK-PS031021.D	2,4-DCAA	1	500	459	92		61	136
		2,4-DCAA	2	500	490	98		61	136
PB168811BL	PB168811BL	2,4-DCAA	1	500	401	80		10	141
		2,4-DCAA	2	500	474	95		10	141
PB168811BS	PB168811BS	2,4-DCAA	1	500	517	103		10	141
		2,4-DCAA	2	500	456	91		10	141
I.BLK-PS031027.D	PIBLK-PS031027.D	2,4-DCAA	1	500	428	86		61	136
		2,4-DCAA	2	500	483	97		61	136
I.BLK-PS031131.D	PIBLK-PS031131.D	2,4-DCAA	1	500	472	94		61	136
		2,4-DCAA	2	500	505	101		61	136
Q2529-01	TP-91	2,4-DCAA	1	500	121	24		10	141
		2,4-DCAA	2	500	128	26		10	141
Q2529-02	TP-80	2,4-DCAA	1	500	63.2	13		10	141
		2,4-DCAA	2	500	261	52		10	141
Q2529-03	TP-79	2,4-DCAA	1	500	111	22		10	141
		2,4-DCAA	2	500	247	49		10	141
Q2529-04	TP-95	2,4-DCAA	1	500	139	28		10	141
		2,4-DCAA	2	500	311	62		10	141
Q2529-05	TP-98	2,4-DCAA	1	500	116	23		10	141
		2,4-DCAA	2	500	216	43		10	141
Q2529-06	TP-102	2,4-DCAA	1	500	149	30		10	141
		2,4-DCAA	2	500	188	38		10	141
Q2529-07	TP-101	2,4-DCAA	1	500	164	33		10	141
		2,4-DCAA	2	500	268	54		10	141
Q2529-08	TP-89	2,4-DCAA	1	500	115	23		10	141
		2,4-DCAA	2	500	185	37		10	141
Q2529-09	TP-33	2,4-DCAA	1	500	107	21		10	141
		2,4-DCAA	2	500	59.6	12		10	141
Q2529-10	TP-30	2,4-DCAA	1	500	118	24		10	141
		2,4-DCAA	2	500	190	38		10	141
I.BLK-PS031143.D	PIBLK-PS031143.D	2,4-DCAA	1	500	473	95		61	136
		2,4-DCAA	2	500	506	101		61	136
Q2529-10MS	TP-30MS	2,4-DCAA	1	500	155	31		10	141
		2,4-DCAA	2	500	220	44		10	141
Q2529-10MSD	TP-30MSD	2,4-DCAA	1	500	180	36		10	141
		2,4-DCAA	2	500	257	51		10	141
I.BLK-PS031153.D	PIBLK-PS031153.D	2,4-DCAA	1	500	478	96		61	136
		2,4-DCAA	2	500	512	102		61	136