

Surrogate Summary

SDG No.: **Q2543**

Client: **CDM Smith**

Analytical Method: **8151A**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS031005.D	PIBLK-PS031005.D	2,4-DCAA	1	500	356	71		61	136
		2,4-DCAA	2	500	496	99		61	136
I.BLK-PS031021.D	PIBLK-PS031021.D	2,4-DCAA	1	500	459	92		61	136
		2,4-DCAA	2	500	490	98		61	136
PB168811BL	PB168811BL	2,4-DCAA	1	500	401	80		10	141
		2,4-DCAA	2	500	474	95		10	141
PB168811BS	PB168811BS	2,4-DCAA	1	500	517	103		10	141
		2,4-DCAA	2	500	456	91		10	141
I.BLK-PS031027.D	PIBLK-PS031027.D	2,4-DCAA	1	500	428	86		61	136
		2,4-DCAA	2	500	483	97		61	136
I.BLK-PS031143.D	PIBLK-PS031143.D	2,4-DCAA	1	500	473	95		61	136
		2,4-DCAA	2	500	506	101		61	136
Q2529-10MS	TP-30MS	2,4-DCAA	1	500	155	31		10	141
		2,4-DCAA	2	500	220	44		10	141
Q2529-10MSD	TP-30MSD	2,4-DCAA	1	500	180	36		10	141
		2,4-DCAA	2	500	257	51		10	141
Q2543-01	TP-41	2,4-DCAA	1	500	244	49		10	141
		2,4-DCAA	2	500	449	90		10	141
Q2543-02	TP-49	2,4-DCAA	1	500	268	54		10	141
		2,4-DCAA	2	500	377	75		10	141
Q2543-03	TP-23	2,4-DCAA	1	500	326	65		10	141
		2,4-DCAA	2	500	414	83		10	141
Q2543-04	TP-23-99	2,4-DCAA	1	500	326	65		10	141
		2,4-DCAA	2	500	401	80		10	141
I.BLK-PS031153.D	PIBLK-PS031153.D	2,4-DCAA	1	500	478	96		61	136
		2,4-DCAA	2	500	512	102		61	136