

Surrogate Summary

SDG No.: Q2579

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PO112087.D	PIBLK-PO112087.D	Tetrachloro-m-xyl	1	20	19.4	97		70 (60)	130 (140)
		Decachlorobiphen	1	20	18.8	94		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	20.0	100		70 (60)	130 (140)
		Decachlorobiphen	2	20	19.3	97		70 (60)	130 (140)
I.BLK-PO112263.D	PIBLK-PO112263.D	Tetrachloro-m-xyl	1	20	19.9	100		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.4	97		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	19.3	97		70 (60)	130 (140)
		Decachlorobiphen	2	20	19.3	97		70 (60)	130 (140)
Q2590-01MS	NB-7-071125MS	Tetrachloro-m-xyl	1	20	22.3	112		30 (32)	150 (144)
		Decachlorobiphen	1	20	13.3	66		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	21.3	107		30 (32)	150 (144)
		Decachlorobiphen	2	20	15.3	76		30 (32)	150 (175)
Q2590-01MSD	NB-7-071125MSD	Tetrachloro-m-xyl	1	20	22.2	111		30 (32)	150 (144)
		Decachlorobiphen	1	20	14.4	72		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	21.3	107		30 (32)	150 (144)
		Decachlorobiphen	2	20	16.8	84		30 (32)	150 (175)
I.BLK-PO112274.D	PIBLK-PO112274.D	Tetrachloro-m-xyl	1	20	22.1	110		70 (60)	130 (140)
		Decachlorobiphen	1	20	17.0	85		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	21.1	106		70 (60)	130 (140)
		Decachlorobiphen	2	20	18.7	93		70 (60)	130 (140)
I.BLK-PP073553.D	PIBLK-PP073553.D	Tetrachloro-m-xyl	1	20	16.2	81		70 (60)	130 (140)
		Decachlorobiphen	1	20	17.3	87		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	16.2	81		70 (60)	130 (140)
		Decachlorobiphen	2	20	17.1	85		70 (60)	130 (140)
I.BLK-PP073772.D	PIBLK-PP073772.D	Tetrachloro-m-xyl	1	20	20.0	100		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.6	98		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	19.0	95		70 (60)	130 (140)
		Decachlorobiphen	2	20	21.3	106		70 (60)	130 (140)
PB168826BL	PB168826BL	Tetrachloro-m-xyl	1	20	20.0	100		30 (32)	150 (144)
		Decachlorobiphen	1	20	20.0	100		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	18.7	94		30 (32)	150 (144)
		Decachlorobiphen	2	20	21.3	106		30 (32)	150 (175)
PB168826BS	PB168826BS	Tetrachloro-m-xyl	1	20	21.2	106		30 (32)	150 (144)
		Decachlorobiphen	1	20	21.0	105		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	18.9	95		30 (32)	150 (144)
		Decachlorobiphen	2	20	22.8	114		30 (32)	150 (175)
Q2579-02	WC-A2-13-C	Tetrachloro-m-xyl	1	20	15.3	77		30 (32)	150 (144)
		Decachlorobiphen	1	20	8.73	44		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	17.4	87		30 (32)	150 (144)
		Decachlorobiphen	2	20	10.9	55		30 (32)	150 (175)
Q2579-06	WC-A2-14-C	Tetrachloro-m-xyl	1	20	13.6	68		30 (32)	150 (144)

() = LABORATORY INHOUSE LIMIT

Surrogate Summary

SDG No.: Q2579

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
Q2579-06	WC-A2-14-C	Decachlorobiphen	1	20	9.66	48		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	16.3	82		30 (32)	150 (144)
		Decachlorobiphen	2	20	11.7	58		30 (32)	150 (175)
I.BLK-PP073787.D	PIBLK-PP073787.D	Tetrachloro-m-xyl	1	20	19.1	95		70 (60)	130 (140)
		Decachlorobiphen	1	20	18.9	95		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	18.1	90		70 (60)	130 (140)
		Decachlorobiphen	2	20	20.3	102		70 (60)	130 (140)