

Surrogate Summary

SDG No.: Q2579

Client: ENTACT

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS031005.D	PIBLK-PS031005.D	2,4-DCAA	1	500	356	71		70 (61)	130 (136)
		2,4-DCAA	2	500	496	99		70 (61)	130 (136)
I.BLK-PS031030.D	PIBLK-PS031030.D	2,4-DCAA	1	500	445	89		70 (61)	130 (136)
		2,4-DCAA	2	500	482	96		70 (61)	130 (136)
PB168844BL	PB168844BL	2,4-DCAA	1	500	423	85		70 (61)	130 (136)
		2,4-DCAA	2	500	499	100		70 (61)	130 (136)
PB168844BS	PB168844BS	2,4-DCAA	1	500	435	87		70 (61)	130 (136)
		2,4-DCAA	2	500	441	88		70 (61)	130 (136)
PB168803TB	PB168803TB	2,4-DCAA	1	500	462	92		70 (61)	130 (136)
		2,4-DCAA	2	500	472	94		70 (61)	130 (136)
Q2489-01MS	G4(0-6)MS	2,4-DCAA	1	500	609	122		70 (61)	130 (136)
		2,4-DCAA	2	500	529	106		70 (61)	130 (136)
Q2489-01MSD	G4(0-6)MSD	2,4-DCAA	1	500	636	127		70 (61)	130 (136)
		2,4-DCAA	2	500	516	103		70 (61)	130 (136)
I.BLK-PS031046.D	PIBLK-PS031046.D	2,4-DCAA	1	500	444	89		70 (61)	130 (136)
		2,4-DCAA	2	500	486	97		70 (61)	130 (136)
I.BLK-PS031061.D	PIBLK-PS031061.D	2,4-DCAA	1	500	450	90		70 (61)	130 (136)
		2,4-DCAA	2	500	496	99		70 (61)	130 (136)
Q2579-03	WC-A2-13-C	2,4-DCAA	1	500	437	87		70 (61)	130 (136)
		2,4-DCAA	2	500	489	98		70 (61)	130 (136)
Q2579-07	WC-A2-14-C	2,4-DCAA	1	500	404	81		70 (61)	130 (136)
		2,4-DCAA	2	500	444	89		70 (61)	130 (136)
I.BLK-PS031071.D	PIBLK-PS031071.D	2,4-DCAA	1	500	468	94		70 (61)	130 (136)
		2,4-DCAA	2	500	498	100		70 (61)	130 (136)