

SOP ID :	M1311-TCLP-16		
SDG No :	N/A	Start Prep Date :	N/A Time : N/A
Weigh By :	N/A	End Prep Date :	N/A Time : N/A
Balance ID :	N/A	Combination Ratio :	N/A
pH Meter ID :	WC PH METER-1	ZHE Cleaning Batch :	N/A
Extraction By :	JP	Initial Room Temperature:	N/A
Filter By :	JP	Final Room Temperature:	N/A
Pipette ID :	N/A	TCLP Technician Signature :	JP
Tumbler ID :	N/A	Supervisor By :	12
TCLP Filter ID :	115525		

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
N/A	N/A	N/A
N/A	N/A	N/A
HNO3-TCLP,1N	N/A	WP112799
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	W1940,W1941,W1942	W3166,W1938,W1939,
1 Liter Amber	N/A	90924-08
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. q2646-03 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
01/21/25 11:30	JP 1st Room	SPS. RJ / EAT
	Preparation Group	Analysis Group 1st Room

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB168926TB	LEB926	N/A	N/A	N/A	N/A	N/A	N/A	4.93	1.5	N/A
Q2646-03	FRAC TANK	N/A	N/A	N/A	N/A	N/A	N/A	4.5	1.0	N/A

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB168926TB	LEB926	N/A	N/A	N/A	N/A	N/A	N/A
Q2646-03	FRAC TANK	N/A	N/A	N/A	N/A	<0.5	N/A