

Surrogate Summary

SDG No.: Q2667

Client: Tetra Tech, EMI

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS031156.D	PIBLK-PS031156.D	2,4-DCAA	1	500	397	79		61	136
		2,4-DCAA	2	500	504	101		61	136
I.BLK-PS031232.D	PIBLK-PS031232.D	2,4-DCAA	1	500	412	82		61	136
		2,4-DCAA	2	500	491	98		61	136
Q2641-02MS	P001-CONCRETE001-01MS	2,4-DCAA	1	500	452	90		61	136
		2,4-DCAA	2	500	499	100		61	136
Q2641-02MSD	P001-CONCRETE001-01MSD	2,4-DCAA	1	500	418	84		61	136
		2,4-DCAA	2	500	463	93		61	136
PB169001BL	PB169001BL	2,4-DCAA	1	500	407	81		61	136
		2,4-DCAA	2	500	473	95		61	136
PB169001BS	PB169001BS	2,4-DCAA	1	500	488	98		61	136
		2,4-DCAA	2	500	501	100		61	136
PB168953TB	PB168953TB	2,4-DCAA	1	500	579	116		61	136
		2,4-DCAA	2	500	500	100		61	136
I.BLK-PS031242.D	PIBLK-PS031242.D	2,4-DCAA	1	500	422	84		61	136
		2,4-DCAA	2	500	498	100		61	136
Q2667-01	C0AP2	2,4-DCAA	1	500	429	86		61	136
		2,4-DCAA	2	500	487	97		61	136
Q2667-02	C0AP3	2,4-DCAA	1	500	403	81		61	136
		2,4-DCAA	2	500	483	97		61	136
I.BLK-PS031251.D	PIBLK-PS031251.D	2,4-DCAA	1	500	374	75		61	136
		2,4-DCAA	2	500	470	94		61	136