

**Laboratory Control Sample/Laboratory Control Sample Duplicate Summary**

**SW-846**

|                 |               |                           |                   |
|-----------------|---------------|---------------------------|-------------------|
| <b>SDG No.:</b> | <u>Q2732</u>  | <b>Analytical Method:</b> | <u>8151A</u>      |
| <b>Client:</b>  | <u>ENTACT</u> | <b>Datafile :</b>         | <u>PS031456.D</u> |

| Lab Sample ID            | Parameter        | Spike | Result | Units | Rec | RPD | Qual | RPD  | Limits  |           | RPD |
|--------------------------|------------------|-------|--------|-------|-----|-----|------|------|---------|-----------|-----|
|                          |                  |       |        |       |     |     |      | Qual | Low     | High      |     |
| PB169278BS<br>(Column 1) | 2,4-D            | 5     | 5.30   | ug/L  | 106 |     |      |      | 70 (83) | 130 (130) |     |
|                          | 2,4,5-TP(Silvex) | 5     | 5.50   | ug/L  | 110 |     |      |      | 70 (78) | 130 (127) |     |
| PB169278BS<br>(Column 2) | 2,4-D            | 5     | 4.80   | ug/L  | 96  |     |      |      | 70 (83) | 130 (130) |     |
|                          | 2,4,5-TP(Silvex) | 5     | 5.00   | ug/L  | 100 |     |      |      | 70 (78) | 130 (127) |     |