

### Surrogate Summary

SDG No.: Q2732

Client: ENTACT

Analytical Method: 8081B

| Lab Sample ID    | Client ID        | Parameter         | Column | Spike | Result | Recovery(%) | Qual | Limits(%) |           |
|------------------|------------------|-------------------|--------|-------|--------|-------------|------|-----------|-----------|
|                  |                  |                   |        |       |        |             |      | Low       | High      |
| I.BLK-PD089685.D | PIBLK-PD089685.D | Decachlorobiphen  | 1      | 20    | 20.9   | 104         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 19.6   | 98          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 20.9   | 104         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 20.4   | 102         |      | 30 (61)   | 150 (148) |
| I.BLK-PD089913.D | PIBLK-PD089913.D | Decachlorobiphen  | 1      | 20    | 19.6   | 98          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 19.3   | 96          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 23.1   | 116         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 24.5   | 123         |      | 30 (61)   | 150 (148) |
| Q2732-03         | WC-A7-01-C       | Decachlorobiphen  | 1      | 20    | 18.1   | 90          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 16.7   | 83          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 16.9   | 85          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 19.9   | 99          |      | 30 (61)   | 150 (148) |
| Q2832-02MS       | TG-S01MS         | Decachlorobiphen  | 1      | 20    | 16.6   | 83          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 18.9   | 94          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 17.4   | 87          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 22.7   | 113         |      | 30 (61)   | 150 (148) |
| Q2832-02MSD      | TG-S01MSD        | Decachlorobiphen  | 1      | 20    | 16.6   | 83          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 18.9   | 94          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 18.2   | 91          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 22.6   | 113         |      | 30 (61)   | 150 (148) |
| I.BLK-PD089925.D | PIBLK-PD089925.D | Decachlorobiphen  | 1      | 20    | 21.0   | 105         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 20.1   | 100         |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 23.9   | 119         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 25.4   | 127         |      | 30 (61)   | 150 (148) |
| I.BLK-PD089934.D | PIBLK-PD089934.D | Decachlorobiphen  | 1      | 20    | 21.5   | 107         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 20.1   | 100         |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 23.5   | 118         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 24.9   | 125         |      | 30 (61)   | 150 (148) |
| PB169263BL       | PB169263BL       | Decachlorobiphen  | 1      | 20    | 22.3   | 112         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 20.4   | 102         |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 26.0   | 130         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 26.3   | 131         |      | 30 (61)   | 150 (148) |
| PB169263BS       | PB169263BS       | Decachlorobiphen  | 1      | 20    | 19.7   | 98          |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 17.6   | 88          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 22.7   | 114         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 22.5   | 112         |      | 30 (61)   | 150 (148) |
| PB169212TB       | PB169212TB       | Decachlorobiphen  | 1      | 20    | 21.3   | 106         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 1      | 20    | 19.4   | 97          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 24.9   | 125         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 25.0   | 125         |      | 30 (61)   | 150 (148) |
| I.BLK-PD089945.D | PIBLK-PD089945.D | Decachlorobiphen  | 1      | 20    | 19.3   | 97          |      | 30 (57)   | 150 (171) |

### Surrogate Summary

SDG No.: Q2732

Client: ENTACT

Analytical Method: 8081B

| Lab Sample ID    | Client ID        | Parameter         | Column | Spike | Result | Recovery(%) | Qual | Limits(%) |           |
|------------------|------------------|-------------------|--------|-------|--------|-------------|------|-----------|-----------|
|                  |                  |                   |        |       |        |             |      | Low       | High      |
| I.BLK-PD089945.D | PIBLK-PD089945.D | Tetrachloro-m-xyl | 1      | 20    | 19.3   | 97          |      | 30 (61)   | 150 (148) |
|                  |                  | Decachlorobiphen  | 2      | 20    | 22.7   | 114         |      | 30 (57)   | 150 (171) |
|                  |                  | Tetrachloro-m-xyl | 2      | 20    | 25.1   | 126         |      | 30 (61)   | 150 (148) |