

SOP ID :	M1311-TCLP-16	
SDG No :	N/A	Start Prep Date : 08/12/2025 Time : 17:00
Weigh By :	JP	End Prep Date : 08/13/2025 Time : 10:35
Balance ID :	WC SC-7	Combination Ratio : 20
pH Meter ID :	WC PH METER-1	ZHE Cleaning Batch : N/A 70
Extraction By :	JP	Initial Room Temperature: 24 °C
Filter By :	JP	Final Room Temperature: 22 °C
Pipette ID :	WC	TCLP Technician Signature : <i>JP</i>
Tumbler ID :	T-1 / T-2	Supervisor By : <i>RM</i>
TCLP Filter ID :	115525	

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
TCLP-FLUID-1	N/A	WP112804
HCL-TCLP,1N	N/A	WP112797
HNO3-TCLP,1N	N/A	WP112799
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	W1940,W1941,W1942	W3166,W1938,W1939,
1 Liter Amber	N/A	90924-08
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. TUMBLER T-1 / T-2 checked,30 rpm. Particle size reduction is not required. Q2838-08 IS USED for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
08/13/25 11:50	<i>JP</i> / <i>ICN Room</i>	<i>SPS.</i> / <i>RJ / E+L</i>
	Preparation Group	Analysis Group

/ RM

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB169212TB	LEB212	13	N/A	2000	N/A	N/A	N/A	4.93	1.5	T-2
Q2732-03	WC-A7-01-C	01	100.02	2000	N/A	N/A	N/A	11.0	1.0	T-1
Q2832-02	TG-S01	02	100.02	2000	N/A	N/A	N/A	6.0	1.5	T-1
Q2832-04	TG-S02	03	100.03	2000	N/A	N/A	N/A	6.2	1.0	T-1
Q2832-06	TG-S03	04	100.02	2000	N/A	N/A	N/A	5.8	1.5	T-1
Q2832-08	TG-S04	05	100.01	2000	N/A	N/A	N/A	5.6	1.0	T-1
Q2832-10	TG-S05	06	100.02	2000	N/A	N/A	N/A	5.0	1.0	T-1
Q2836-03	WC-A2-15-C	07	100.03	2000	N/A	N/A	N/A	11.5	1.5	T-1
Q2836-07	WC-A2-16-C	08	100.02	2000	N/A	N/A	N/A	11.5	1.0	T-1
Q2836-11	WC-A2-17-C	09	100.03	2000	N/A	N/A	N/A	11.5	1.5	T-1
Q2836-15	WC-A5-02-C	10	100.02	2000	N/A	N/A	N/A	11.0	1.5	T-1
Q2838-04	TP-11	11	100.01	2000	N/A	N/A	N/A	5.6	1.0	T-2
Q2838-08	TP-10	12	100.02	2000	N/A	N/A	N/A	5.6	1.5	T-2

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB169212TB	LEB212	N/A	N/A	N/A	N/A	N/A	N/A
Q2732-03	WC-A7-01-C	N/A	N/A	N/A	N/A	100	N/A
Q2832-02	TG-S01	N/A	N/A	N/A	N/A	100	N/A
Q2832-04	TG-S02	N/A	N/A	N/A	N/A	100	N/A
Q2832-06	TG-S03	N/A	N/A	N/A	N/A	100	N/A
Q2832-08	TG-S04	N/A	N/A	N/A	N/A	100	N/A
Q2832-10	TG-S05	N/A	N/A	N/A	N/A	100	N/A
Q2836-03	WC-A2-15-C	N/A	N/A	N/A	N/A	100	N/A
Q2836-07	WC-A2-16-C	N/A	N/A	N/A	N/A	100	N/A
Q2836-11	WC-A2-17-C	N/A	N/A	N/A	N/A	100	N/A
Q2836-15	WC-A5-02-C	N/A	N/A	N/A	N/A	100	N/A
Q2838-04	TP-11	N/A	N/A	N/A	N/A	100	N/A
Q2838-08	TP-10	N/A	N/A	N/A	N/A	100	N/A

Hot Block ID : WC S-1 / WC S-2

 Thermometer ID : FLASHPOINT

SampleID	ClientID	Sample Weight (g)	Volume DI Water (mL)	PH after 5 min stir	PH after 10 min stir	Extraction Fluid 1 or 2	pH Extraction Fluid
PB169212TB	LEB212	N/A	N/A	N/A	N/A	#1	4.93
Q2732-03	WC-A7-01-C	5.02	96.5	11.5	4.5	#1	4.93
Q2832-02	TG-S01	5.02	96.5	8.4	3.5	#1	4.93
Q2832-04	TG-S02	5.01	96.5	8.4	3.5	#1	4.93
Q2832-06	TG-S03	5.02	96.5	8.2	3.0	#1	4.93
Q2832-08	TG-S04	5.03	96.5	8.0	3.0	#1	4.93
Q2832-10	TG-S05	5.04	96.5	7.0	2.5	#1	4.93
Q2836-03	WC-A2-15-C	5.03	96.5	12.0	4.5	#1	4.93
Q2836-07	WC-A2-16-C	5.02	96.5	12.0	4.5	#1	4.93
Q2836-11	WC-A2-17-C	5.03	96.5	12.5	4.5	#1	4.93
Q2836-15	WC-A5-02-C	5.02	96.5	11.5	4.5	#1	4.93
Q2838-04	TP-11	5.03	96.5	7.6	3.0	#1	4.93
Q2838-08	TP-10	5.02	96.5	7.2	2.5	#1	4.93