



PERCENT SOLID

Supervisor: Iwona
Analyst: rubina
Date: 8/1/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 07/31/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 08/01/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB136669

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q2734-01	VNJ-253	1	1.14	10.65	11.79	10.74	90.1	
Q2735-01	BC271315-1-1	2	1.00	1.00	2.00	2.00	100.0	PILC SAMPLE
Q2735-02	BC271315-1-2	3	1.00	1.00	2.00	2.00	100.0	PILC SAMPLE
Q2735-03	BC274768-1-1	4	1.00	1.00	2.00	2.00	100.0	PILC SAMPLE
Q2735-04	BC274768-1-2	5	1.00	1.00	2.00	2.00	100.0	PILC SAMPLE
Q2736-01	PAD-073125	6	1.00	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
Q2737-01	OILY-DEBRIS	7	1.14	10.72	11.86	7.87	62.8	
Q2740-01	OR-02-07312025	8	1.12	10.54	11.66	11.09	94.6	
Q2740-02	OR-02-07312025-E2	9	1.14	10.85	11.99	11.2	92.7	
Q2741-01	TR-05-07212025	10	1.16	10.34	11.5	10.83	93.5	
Q2741-02	TR-05-07212025-E2	11	1.12	10.20	11.32	10.12	88.2	
Q2743-01	WC1	12	1.13	10.67	11.8	9.48	78.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$