

**metals**  
**- 14 -**  
**ANALYSIS RUN LOG**

**Client:** Aramark Uniforms

**Contract:** ARAM01

**Lab code:** ACE

**Sdg no.:** Q2791

**Instrument id number:** \_\_\_\_\_

**Method:** \_\_\_\_\_

**Run number:** LB136774

**Start date:** 08/11/2025

**End date:** 08/11/2025

| Lab sample id. | Client Sample Id | d/f | Time | Parameter list |
|----------------|------------------|-----|------|----------------|
| S0             | S0               | 1   | 1617 | HG             |
| S0.2           | S0.2             | 1   | 1622 | HG             |
| S2.5           | S2.5             | 1   | 1624 | HG             |
| S5             | S5               | 1   | 1632 | HG             |
| S7.5           | S7.5             | 1   | 1650 | HG             |
| S10            | S10              | 1   | 1701 | HG             |
| ICV07          | ICV07            | 1   | 1707 | HG             |
| ICB07          | ICB07            | 1   | 1709 | HG             |
| CCV19          | CCV19            | 1   | 1722 | HG             |
| CCB19          | CCB19            | 1   | 1724 | HG             |
| CRA            | CRA              | 1   | 1726 | HG             |
| PB169166BL     | PB169166BL       | 1   | 1745 | HG             |
| PB169166BS     | PB169166BS       | 1   | 1747 | HG             |
| Q2791-01       | SLUDGE           | 1   | 1757 | HG             |
| CCV20          | CCV20            | 1   | 1801 | HG             |
| CCB20          | CCB20            | 1   | 1803 | HG             |
| Q2800-02DUP    | VNJ-222DUP       | 1   | 1808 | HG             |
| Q2800-02MS     | VNJ-222MS        | 1   | 1810 | HG             |
| Q2800-02MSD    | VNJ-222MSD       | 1   | 1812 | HG             |
| Q2800-02L      | VNJ-222L         | 5   | 1815 | HG             |
| PB169166TB     | PB169166TB       | 1   | 1819 | HG             |
| CCV21          | CCV21            | 1   | 1822 | HG             |
| CCB21          | CCB21            | 1   | 1824 | HG             |

**metals**  
**- 14 -**  
**ANALYSIS RUN LOG**

**Client:** Aramark Uniforms

**Contract:** ARAM01

**Lab code:** ACE

**Sdg no.:** Q2791

**Instrument id number:** \_\_\_\_\_

**Method:** \_\_\_\_\_

**Run number:** LB136778

**Start date:** 08/11/2025

**End date:** 08/11/2025

| Lab sample id. | Client Sample Id | d/f | Time | Parameter list       |
|----------------|------------------|-----|------|----------------------|
| S0             | S0               | 1   | 1043 | Ag,As,Ba,Cd,Cr,Pb,Se |
| S1             | S1               | 1   | 1048 | Ag,As,Ba,Cd,Cr,Pb,Se |
| S2             | S2               | 1   | 1052 | Ag,As,Ba,Cd,Cr,Pb,Se |
| S3             | S3               | 1   | 1056 | Ag,As,Ba,Cd,Cr,Pb,Se |
| S4             | S4               | 1   | 1100 | Ag,As,Ba,Cd,Cr,Pb,Se |
| S5             | S5               | 1   | 1104 | Ag,As,Ba,Cd,Cr,Pb,Se |
| ICV01          | ICV01            | 1   | 1109 | Ag,As,Ba,Cd,Cr,Pb,Se |
| LLICV01        | LLICV01          | 1   | 1133 | Ag,As,Ba,Cd,Cr,Pb,Se |
| ICB01          | ICB01            | 1   | 1137 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CRI01          | CRI01            | 1   | 1141 | Ag,As,Ba,Cd,Cr,Pb,Se |
| ICSA01         | ICSA01           | 1   | 1148 | Ag,As,Ba,Cd,Cr,Pb,Se |
| ICSAB01        | ICSAB01          | 1   | 1152 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CCV01          | CCV01            | 1   | 1205 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CCB01          | CCB01            | 1   | 1209 | Ag,As,Ba,Cd,Cr,Pb,Se |
| PB169178BL     | PB169178BL       | 1   | 1222 | Ag,As,Ba,Cd,Cr,Pb,Se |
| PB169178BS     | PB169178BS       | 1   | 1226 | Ag,As,Ba,Cd,Cr,Pb,Se |
| PB169149TB     | PB169149TB       | 1   | 1242 | Ag,As,Ba,Cd,Cr,Pb,Se |
| Q2791-01       | SLUDGE           | 1   | 1256 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CCV02          | CCV02            | 1   | 1300 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CCB02          | CCB02            | 1   | 1304 | Ag,As,Ba,Cd,Cr,Pb,Se |
| Q2800-02DUP    | VNJ-222DUP       | 1   | 1317 | Ag,As,Ba,Cd,Cr,Pb,Se |
| Q2800-02L      | VNJ-222L         | 5   | 1322 | Ag,As,Ba,Cd,Cr,Pb,Se |
| Q2800-02MS     | VNJ-222MS        | 1   | 1326 | Ag,As,Ba,Cd,Cr,Pb,Se |
| Q2800-02MSD    | VNJ-222MSD       | 1   | 1330 | Ag,As,Ba,Cd,Cr,Pb,Se |
| Q2800-02A      | VNJ-222A         | 1   | 1334 | Ag,Ba                |
| CCV03          | CCV03            | 1   | 1454 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CCB03          | CCB03            | 1   | 1458 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CCV04          | CCV04            | 1   | 1602 | Ag,As,Ba,Cd,Cr,Pb,Se |
| CCB04          | CCB04            | 1   | 1606 | Ag,As,Ba,Cd,Cr,Pb,Se |