

Surrogate Summary

SW-846

SDG No.: Q2791

Client: Aramark Uniforms

Analytical Method: 8270E

Lab Sample ID	Client ID	Parameter	Spike (PPM)	Result (PPM)	Recovery (%)	Qual	Limits (%)	
							Low	High
PB169149TB	PB169149TB	2-Fluorophenol	150	119	79		23	138
		Phenol-d6	150	106	71		10	134
		Nitrobenzene-d5	100	82.0	82		67	132
		2-Fluorobiphenyl	100	78.2	78		52	132
		2,4,6-Tribromophenol	150	142	94		44	137
		Terphenyl-d14	100	76.4	76		42	152
PB169180BL	PB169180BL	2-Fluorophenol	150	124	83		23	138
		Phenol-d6	150	109	73		10	134
		Nitrobenzene-d5	100	86.1	86		67	132
		2-Fluorobiphenyl	100	84.5	84		52	132
		2,4,6-Tribromophenol	150	145	97		44	137
		Terphenyl-d14	100	80.7	81		42	152
PB169180BS	PB169180BS	2-Fluorophenol	150	129	86		23	138
		Phenol-d6	150	120	80		10	134
		Nitrobenzene-d5	100	85.7	86		67	132
		2-Fluorobiphenyl	100	85.5	85		52	132
		2,4,6-Tribromophenol	150	156	104		44	137
		Terphenyl-d14	100	80.4	80		42	152
Q2779-04MS	TP-4MS	2-Fluorophenol	150	113	76		23	138
		Phenol-d6	150	96.1	64		10	134
		Nitrobenzene-d5	100	84.2	84		67	132
		2-Fluorobiphenyl	100	79.2	79		52	132
		2,4,6-Tribromophenol	150	145	97		44	137
		Terphenyl-d14	100	78.3	78		42	152
Q2779-04MSD	TP-4MSD	2-Fluorophenol	150	107	71		23	138
		Phenol-d6	150	91.7	61		10	134
		Nitrobenzene-d5	100	80.3	80		67	132
		2-Fluorobiphenyl	100	75.8	76		52	132
		2,4,6-Tribromophenol	150	135	90		44	137
		Terphenyl-d14	100	80.5	81		42	152
Q2791-01	SLUDGE	2-Fluorophenol	150	116	77		23	138
		Phenol-d6	150	104	70		10	134
		Nitrobenzene-d5	100	85.3	85		67	132
		2-Fluorobiphenyl	100	73.6	74		52	132
		2,4,6-Tribromophenol	150	138	92		44	137
		Terphenyl-d14	100	64.5	65		42	152