

Analysis Method: 9034

ANALYST: rubina

Parameter: Reactive Sulfide

SUPERVISOR REVIEW BY: jignesh

Run Number: LB136810

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3213
Starch Solution, 4L	W3149

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB169218BL		1	5.00	50	2.00	0.00	1.94	1.94	0.06	0.00	0.00	08/13/2025	15:00
2	Q2732-03		1	5.08	50	2.00	0.00	1.88	1.88	0.12	0.06	4.72	08/13/2025	15:03
3	Q2732-03DUP		1	5.08	50	2.00	0.00	1.88	1.88	0.12	0.06	4.72	08/13/2025	15:06
4	Q2826-02		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	08/13/2025	15:08
5	Q2827-04		1	5.02	50	2.00	0.00	1.86	1.86	0.14	0.08	6.37	08/13/2025	15:10
6	Q2827-08		1	5.04	50	2.00	0.00	1.88	1.88	0.12	0.06	4.76	08/13/2025	15:13
7	Q2831-03		1	5.04	50	2.00	0.00	1.86	1.86	0.14	0.08	6.35	08/13/2025	15:15
8	Q2832-02		1	5.01	50	2.00	0.00	1.86	1.86	0.14	0.08	6.39	08/13/2025	15:18
9	Q2832-04		1	5.07	50	2.00	0.00	1.90	1.90	0.10	0.04	3.16	08/13/2025	15:20
10	Q2832-06		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	08/13/2025	15:23
11	Q2832-08		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	08/13/2025	15:25
12	Q2832-10		1	5.06	50	2.00	0.00	1.90	1.90	0.10	0.04	3.16	08/13/2025	15:28

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume

Analysis Method: 9034

Parameter: Reactive Sulfide

Run Number: LB136810

ANALYST: rubina

SUPERVISOR REVIEW BY: jignesh

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3213
Starch Solution, 4L	W3149

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
13	Q2836-03		1	5.03	50	2.00	0.00	1.90	1.90	0.10	0.04	3.18	08/13/2025	15:30
14	Q2836-07		1	5.03	50	2.00	0.00	1.90	1.90	0.10	0.04	3.18	08/13/2025	15:33
15	Q2836-11		1	5.08	50	2.00	0.00	1.90	1.90	0.10	0.04	3.15	08/13/2025	13:35
16	Q2836-15		1	5.06	50	2.00	0.00	1.86	1.86	0.14	0.08	6.32	08/13/2025	15:38
17	Q2838-04		1	5.01	50	2.00	0.00	1.88	1.88	0.12	0.06	4.79	08/13/2025	15:41
18	Q2838-08		1	5.06	50	2.00	0.00	1.86	1.86	0.14	0.08	6.32	08/13/2025	15:44

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume