

Surrogate Summary

SDG No.: Q2837

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PP074167.D	PIBLK-PP074167.D	Tetrachloro-m-xyl	1	20	19.1	95		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.7	98		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.3	86		70 (60)	130 (140)
		Decachlorobiphen	2	20	18.7	93		70 (60)	130 (140)
I.BLK-PP074409.D	PIBLK-PP074409.D	Tetrachloro-m-xyl	1	20	20.5	103		70 (60)	130 (140)
		Decachlorobiphen	1	20	18.7	94		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.8	89		70 (60)	130 (140)
		Decachlorobiphen	2	20	17.5	87		70 (60)	130 (140)
PB169241BL	PB169241BL	Tetrachloro-m-xyl	1	20	23.5	118		30 (30)	150 (173)
		Decachlorobiphen	1	20	22.4	112		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	22.3	112		30 (30)	150 (173)
		Decachlorobiphen	2	20	21.5	108		30 (10)	150 (173)
PB169241BS	PB169241BS	Tetrachloro-m-xyl	1	20	25.0	125		30 (30)	150 (173)
		Decachlorobiphen	1	20	23.0	115		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	24.4	122		30 (30)	150 (173)
		Decachlorobiphen	2	20	22.0	110		30 (10)	150 (173)
PB169241BSD	PB169241BSD	Tetrachloro-m-xyl	1	20	25.8	129		30 (30)	150 (173)
		Decachlorobiphen	1	20	23.6	118		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	24.7	123		30 (30)	150 (173)
		Decachlorobiphen	2	20	22.3	111		30 (10)	150 (173)
Q2837-01	TW-WTS-13	Tetrachloro-m-xyl	1	20	24.0	120		30 (30)	150 (173)
		Decachlorobiphen	1	20	16.6	83		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	21.9	110		30 (30)	150 (173)
		Decachlorobiphen	2	20	15.4	77		30 (10)	150 (173)
I.BLK-PP074423.D	PIBLK-PP074423.D	Tetrachloro-m-xyl	1	20	20.4	102		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.3	97		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	18.5	92		70 (60)	130 (140)
		Decachlorobiphen	2	20	17.8	89		70 (60)	130 (140)