

RETENTION TIMES OF INITIAL CALIBRATION

Lab Name:	<u>Alliance</u>	Contract:	<u>ALLI03</u>
Lab Code:	<u>ACE</u>	SDG NO.:	<u>Q2893</u>
Instrument ID:	<u>ECD_S</u>	Calibration Date(s):	<u>09/09/2025</u> <u>09/09/2025</u>
		Calibration Times:	<u>11:32</u> <u>13:08</u>

GC Column: RTX-CLP ID: 0.32 (mm)

LAB FILE ID:	RT 200 = <u>PS031668.D</u>	RT 500 = <u>PS031669.D</u>
RT 750 = <u>PS031670.D</u>	RT 1000 = <u>PS031671.D</u>	RT 1500 = <u>PS031672.D</u>

COMPOUND	RT 200	RT 500	RT 750	RT 1000	RT 1500	MEAN RT	RT WINDOW	
							FROM	TO
2,4,5-T	9.59	9.59	9.58	9.58	9.58	9.58	9.48	9.68
2,4,5-TP(Silvex)	9.29	9.29	9.29	9.29	9.29	9.29	9.19	9.39
2,4-D	8.42	8.42	8.41	8.41	8.41	8.42	8.32	8.52
2,4-DB	10.17	10.16	10.16	10.16	10.15	10.16	10.06	10.26
2,4-DCAA	7.29	7.29	7.29	7.29	7.28	7.29	7.19	7.39
DICAMBA	7.47	7.47	7.47	7.47	7.47	7.47	7.37	7.57
DICHLORPROP	8.18	8.18	8.18	8.18	8.18	8.18	8.08	8.28
Dinoseb	11.36	11.36	11.35	11.35	11.35	11.35	11.25	11.45

RETENTION TIMES OF INITIAL CALIBRATION

Lab Name: Alliance

Contract: ALLI03

Lab Code: ACE

SDG NO.: Q2893

Instrument ID: ECD_S

Calibration Date(s): **09/09/2025** **09/09/2025**

Calibration Times: 11:32 13:08

GC Column: RTX-CLP2 ID: 0.32 (mm)

LAB FILE ID:	RT 200 =	<u>PS031668.D</u>	RT 500 =	<u>PS031669.D</u>
	RT 750 =	PS031670.D	RT 1000 =	PS031671.D
			RT 1500 =	PS031672.D

[illegible]