

Instrument ID: ECD_Q

Daily Analysis Runlog For Sequence/QC Batch ID # PQ100725

Review By	yogesh	Review On	10/7/2025 2:41:30 PM
Supervise By	mohammad	Supervise On	10/9/2025 1:41:11 AM
SubDirectory	PQ100725	HP Acquire Method	HP Processing Method PQ100725

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24970,PP24971,PP24972,PP24973,PP24974
CCC Internal Standard/PEM	PP24976,PP24977
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24975

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	HEXANE	HEXANE	PQ070946.D	07 Oct 2025 08:17		YPIAJ	Ok
2	M8011.504.1 0.5 PPB IC	M8011.504.1 0.5 PPB IC	PQ070947.D	07 Oct 2025 08:57		YPIAJ	Ok,M
3	M8011.504.1 0.25 PPB	M8011.504.1 0.25 PPB	PQ070948.D	07 Oct 2025 09:07		YPIAJ	Ok
4	M8011.504.1 0.1 PPB IC	M8011.504.1 0.1 PPB IC	PQ070949.D	07 Oct 2025 10:30		YPIAJ	Ok
5	M8011.504.1 0.05 PPB	M8011.504.1 0.05 PPB	PQ070950.D	07 Oct 2025 10:43		YPIAJ	Ok
6	M8011.504.1 0.025 PPB	M8011.504.1 0.025 PPB	PQ070951.D	07 Oct 2025 11:09		YPIAJ	Ok,M
7	M8011.504.1 0.1 PPB IC	CVPQ100725	PQ070952.D	07 Oct 2025 11:29		YPIAJ	Ok,M
8	RT CHECK	RT CHECK	PQ070953.D	07 Oct 2025 12:19		YPIAJ	Ok
9	M8011.504.1 0.1 PPB CM	M8011.504.1 0.1 PPB CM	PQ070954.D	07 Oct 2025 12:29		YPIAJ	Ok
10	Low-Level LFB-MDL Ch	Low-Level LFB-MDL Ch	PQ070955.D	07 Oct 2025 12:38		YPIAJ	Ok,M
11	PB170012BL	PB170012BL	PQ070956.D	07 Oct 2025 12:48		YPIAJ	Ok
12	PB170012BS	PB170012BS	PQ070957.D	07 Oct 2025 12:58		YPIAJ	Ok
13	PB170012BSD	PB170012BSD	PQ070958.D	07 Oct 2025 13:09		YPIAJ	Ok
14	Q3015-20	WP0925-RR-8011-WP	PQ070959.D	07 Oct 2025 13:18	Need 4X	YPIAJ	Dilution
15	Q3015-20DL	WP0925-RR-8011-WP	PQ070960.D	07 Oct 2025 13:47		YPIAJ	Ok
16	M8011.504.1 0.1 PPB CM	M8011.504.1 0.1 PPB CM	PQ070961.D	07 Oct 2025 13:57		YPIAJ	Ok
17	Q2893-07	LOD-MDL-WATER-01-C	PQ070962.D	07 Oct 2025 14:07		YPIAJ	Ok,M
18	Q2893-08	LOQ-WATER-02-QT3-2	PQ070963.D	07 Oct 2025 14:21		YPIAJ	Ok,M

Instrument ID: ECD_Q

Daily Analysis Runlog For Sequence/QC Batch ID # PQ100725

Review By	yogesh	Review On	10/7/2025 2:41:30 PM
Supervise By	mohammad	Supervise On	10/9/2025 1:41:11 AM
SubDirectory	PQ100725	HP Acquire Method	HP Processing Method PQ100725

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24970,PP24971,PP24972,PP24973,PP24974
CCC Internal Standard/PEM	PP24976,PP24977
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24975

19	M8011.504.1 0.1 PPB C	M8011.504.1 0.1 PPB C	PQ070964.D	07 Oct 2025 14:31		YP\AJ	Ok
----	-----------------------	-----------------------	------------	-------------------	--	-------	----

M : Manual Integration