

RETENTION TIMES OF INITIAL CALIBRATION

Lab Name:	<u>Alliance</u>	Contract:	<u>ZUCC01</u>
Lab Code:	<u>ACE</u>	SDG NO.:	<u>Q3032</u>
Instrument ID:	<u>ECD_S</u>	Calibration Date(s):	<u>09/09/2025</u> <u>09/09/2025</u>
		Calibration Times:	<u>11:32</u> <u>13:08</u>

GC Column: RTX-CLP ID: 0.32 (mm)

LAB FILE ID:	RT 200 = <u>PS031668.D</u>	RT 500 = <u>PS031669.D</u>
RT 750 = <u>PS031670.D</u>	RT 1000 = <u>PS031671.D</u>	RT 1500 = <u>PS031672.D</u>

COMPOUND	RT 200	RT 500	RT 750	RT 1000	RT 1500	MEAN RT	RT WINDOW FROM TO	
2,4,5-TP(Silvex)	9.29	9.29	9.29	9.29	9.29	9.29	9.19	9.39
2,4-D	8.42	8.42	8.41	8.41	8.41	8.42	8.32	8.52
2,4-DCAA	7.29	7.29	7.29	7.29	7.28	7.29	7.19	7.39



Lab Name:	<u>Alliance</u>	Contract:	<u>ZUCC01</u>	
Lab Code:	<u>ACE</u>	SDG NO.:	<u>Q3032</u>	
Instrument ID:	<u>ECD_S</u>	Calibration Date(s):	<u>09/09/2025</u>	<u>09/09/2025</u>
		Calibration Times:	<u>11:32</u>	<u>13:08</u>

LAB FILE ID:	RT 200 =	<u>PS031668.D</u>	RT 500 =	<u>PS031669.D</u>
RT 750 =	<u>PS031670.D</u>	RT 1000 =	<u>PS031671.D</u>	RT 1500 = <u>PS031672.D</u>

[illegible]