

**Instrument ID:** P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137198**

|                  |                                                 |              |                       |
|------------------|-------------------------------------------------|--------------|-----------------------|
| Review By        | Janvi                                           | Review On    | 9/16/2025 11:58:31 AM |
| Supervise By     | jaswal                                          | Supervise On | 9/16/2025 12:00:47 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                       |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                       |
| ICV Standard     | MP87047,MP87038                                 |              |                       |
| CCV Standard     | MP87042                                         |              |                       |
| ICSA Standard    | MP87040,MP87048                                 |              |                       |
| CRI Standard     | MP87038                                         |              |                       |
| LCS Standard     |                                                 |              |                       |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                       |

| Sr# | SampleId   | ClientID   | QcType   | Date           | Comment | Operator | Status |
|-----|------------|------------|----------|----------------|---------|----------|--------|
| 1   | S0         | S0         | CAL1     | 09/15/25 11:06 |         | Jaswal   | OK     |
| 2   | S1         | S1         | CAL2     | 09/15/25 11:10 |         | Jaswal   | OK     |
| 3   | S2         | S2         | CAL3     | 09/15/25 11:14 |         | Jaswal   | OK     |
| 4   | S3         | S3         | CAL4     | 09/15/25 11:18 |         | Jaswal   | OK     |
| 5   | S4         | S4         | CAL5     | 09/15/25 11:22 |         | Jaswal   | OK     |
| 6   | S5         | S5         | CAL6     | 09/15/25 11:26 |         | Jaswal   | OK     |
| 7   | ICV01      | ICV01      | ICV      | 09/15/25 11:31 |         | Jaswal   | OK     |
| 8   | LLICV01    | LLICV01    | LLICV    | 09/15/25 11:36 |         | Jaswal   | OK     |
| 9   | ICB01      | ICB01      | ICB      | 09/15/25 11:43 |         | Jaswal   | OK     |
| 10  | CRI01      | CRI01      | CRDL     | 09/15/25 11:47 |         | Jaswal   | OK     |
| 11  | ICSA01     | ICSA01     | ICSA     | 09/15/25 12:12 |         | Jaswal   | OK     |
| 12  | ICSAB01    | ICSAB01    | ICSAB    | 09/15/25 12:16 |         | Jaswal   | OK     |
| 13  | ICSADL     | ICSADL     | ICSA     | 09/15/25 12:26 |         | Jaswal   | OK     |
| 14  | ICSABDL    | ICSABDL    | ICSAB    | 09/15/25 12:30 |         | Jaswal   | OK     |
| 15  | CCV01      | CCV01      | CCV      | 09/15/25 12:35 |         | Jaswal   | OK     |
| 16  | CCB01      | CCB01      | CCB      | 09/15/25 12:39 |         | Jaswal   | OK     |
| 17  | LR-1       | LR-1       | HIGH STD | 09/15/25 12:44 |         | Jaswal   | OK     |
| 18  | PB169680BL | PB169680BL | MB       | 09/15/25 12:58 |         | Jaswal   | OK     |

Instrument ID: P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137198**

|                  |                                                 |              |                       |
|------------------|-------------------------------------------------|--------------|-----------------------|
| Review By        | Janvi                                           | Review On    | 9/16/2025 11:58:31 AM |
| Supervise By     | jaswal                                          | Supervise On | 9/16/2025 12:00:47 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                       |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                       |
| ICV Standard     | MP87047,MP87038                                 |              |                       |
| CCV Standard     | MP87042                                         |              |                       |
| ICSA Standard    | MP87040,MP87048                                 |              |                       |
| CRI Standard     | MP87038                                         |              |                       |
| LCS Standard     |                                                 |              |                       |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                       |

|    |             |            |          |                |                                           |        |    |
|----|-------------|------------|----------|----------------|-------------------------------------------|--------|----|
| 19 | PB169680BS  | PB169680BS | LCS      | 09/15/25 13:02 |                                           | Jaswal | OK |
| 20 | PB169682BL  | PB169682BL | MB       | 09/15/25 13:06 |                                           | Jaswal | OK |
| 21 | PB169682BS  | PB169682BS | LCS      | 09/15/25 13:10 |                                           | Jaswal | OK |
| 22 | Q3082-01    | WC-29      | SAM      | 09/15/25 13:14 |                                           | Jaswal | OK |
| 23 | Q3082-03    | WC-30      | SAM      | 09/15/25 13:18 |                                           | Jaswal | OK |
| 24 | Q3082-05    | WC-31      | SAM      | 09/15/25 13:22 |                                           | Jaswal | OK |
| 25 | CCV02       | CCV02      | CCV      | 09/15/25 13:42 |                                           | Jaswal | OK |
| 26 | CCB02       | CCB02      | CCB      | 09/15/25 13:47 |                                           | Jaswal | OK |
| 27 | Q3082-07    | WC-32      | SAM      | 09/15/25 13:51 |                                           | Jaswal | OK |
| 28 | Q3082-09    | WC-33      | SAM      | 09/15/25 13:55 |                                           | Jaswal | OK |
| 29 | Q3082-01DUP | WC-29DUP   | DUP      | 09/15/25 13:59 |                                           | Jaswal | OK |
| 30 | Q3082-01L   | WC-29L     | SD       | 09/15/25 14:03 |                                           | Jaswal | OK |
| 31 | Q3082-01MS  | WC-29MS    | MS       | 09/15/25 14:07 |                                           | Jaswal | OK |
| 32 | Q3082-01MSD | WC-29MSD   | MSD      | 09/15/25 14:11 |                                           | Jaswal | OK |
| 33 | Q3082-01A   | WC-29A     | PS       | 09/15/25 14:15 | 0.1 ML OF M6008 AND M6017 IN 10ML SAMPLE. | Jaswal | OK |
| 34 | LR-2        | LR-2       | HIGH STD | 09/15/25 14:40 |                                           | Jaswal | OK |
| 35 | LR-3        | LR-3       | HIGH STD | 09/15/25 14:48 |                                           | Jaswal | OK |
| 36 | CCV03       | CCV03      | CCV      | 09/15/25 15:11 |                                           | Jaswal | OK |
| 37 | CCB03       | CCB03      | CCB      | 09/15/25 15:15 |                                           | Jaswal | OK |

Instrument ID: P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137198**

|                  |                                                 |              |                       |
|------------------|-------------------------------------------------|--------------|-----------------------|
| Review By        | Janvi                                           | Review On    | 9/16/2025 11:58:31 AM |
| Supervise By     | jaswal                                          | Supervise On | 9/16/2025 12:00:47 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                       |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                       |
| ICV Standard     | MP87047,MP87038                                 |              |                       |
| CCV Standard     | MP87042                                         |              |                       |
| ICSA Standard    | MP87040,MP87048                                 |              |                       |
| CRI Standard     | MP87038                                         |              |                       |
| LCS Standard     |                                                 |              |                       |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                       |

|    |             |              |     |                |                                              |        |    |
|----|-------------|--------------|-----|----------------|----------------------------------------------|--------|----|
| 38 | PB169685BL  | PB169685BL   | MB  | 09/15/25 15:19 |                                              | Jaswal | OK |
| 39 | PB169685BS  | PB169685BS   | LCS | 09/15/25 15:23 |                                              | Jaswal | OK |
| 40 | Q3082-02    | WC-29        | SAM | 09/15/25 15:27 |                                              | Jaswal | OK |
| 41 | Q3082-04    | WC-30        | SAM | 09/15/25 15:32 |                                              | Jaswal | OK |
| 42 | Q3082-06    | WC-31        | SAM | 09/15/25 15:36 |                                              | Jaswal | OK |
| 43 | Q3082-08    | WC-32        | SAM | 09/15/25 15:41 |                                              | Jaswal | OK |
| 44 | Q3082-10    | WC-33        | SAM | 09/15/25 15:45 |                                              | Jaswal | OK |
| 45 | Q3085-01    | PAULDING-WC1 | SAM | 09/15/25 15:49 |                                              | Jaswal | OK |
| 46 | Q3091-07    | 12/9C-COMP   | SAM | 09/15/25 15:54 |                                              | Jaswal | OK |
| 47 | Q3091-14    | 12/9B-COMP   | SAM | 09/15/25 15:58 |                                              | Jaswal | OK |
| 48 | CCV04       | CCV04        | CCV | 09/15/25 16:12 |                                              | Jaswal | OK |
| 49 | CCB04       | CCB04        | CCB | 09/15/25 16:16 |                                              | Jaswal | OK |
| 50 | Q3092-02    | 291376       | SAM | 09/15/25 16:21 |                                              | Jaswal | OK |
| 51 | Q3092-04    | 279182       | SAM | 09/15/25 16:25 |                                              | Jaswal | OK |
| 52 | Q3092-04DUP | 279182DUP    | DUP | 09/15/25 16:29 |                                              | Jaswal | OK |
| 53 | Q3092-04L   | 279182L      | SD  | 09/15/25 16:34 |                                              | Jaswal | OK |
| 54 | Q3092-04MS  | 279182MS     | MS  | 09/15/25 16:38 |                                              | Jaswal | OK |
| 55 | Q3092-04MSD | 279182MSD    | MSD | 09/15/25 16:42 |                                              | Jaswal | OK |
| 56 | Q3092-04A   | 279182A      | PS  | 09/15/25 16:46 | 0.1 ML OF M6008 AND<br>M6017 IN 10ML SAMPLE. | Jaswal | OK |

Instrument ID: P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137198**

|              |        |              |                       |
|--------------|--------|--------------|-----------------------|
| Review By    | Janvi  | Review On    | 9/16/2025 11:58:31 AM |
| Supervise By | jaswal | Supervise On | 9/16/2025 12:00:47 PM |

| STD. NAME     | STD REF.#                                       |
|---------------|-------------------------------------------------|
| ICAL Standard | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |
| ICV Standard  | MP87047,MP87038                                 |
| CCV Standard  | MP87042                                         |
| ICSA Standard | MP87040,MP87048                                 |
| CRI Standard  | MP87038                                         |
| LCS Standard  |                                                 |
| Chk Standard  | MP87043,MP87044,MP87045,MP87046,MP87049         |

|    |            |            |     |                |  |        |    |
|----|------------|------------|-----|----------------|--|--------|----|
| 57 | PB169677TB | PB169677TB | MB  | 09/15/25 16:50 |  | Jaswal | OK |
| 58 | CCV05      | CCV05      | CCV | 09/15/25 16:55 |  | Jaswal | OK |
| 59 | CCB05      | CCB05      | CCB | 09/15/25 16:59 |  | Jaswal | OK |

**Instrument ID:** CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137199**

|              |        |              |                       |
|--------------|--------|--------------|-----------------------|
| Review By    | mohan  | Review On    | 9/16/2025 3:04:24 PM  |
| Supervise By | jaswal | Supervise On | 9/16/2025 10:31:06 PM |

| STD. NAME     | STD REF.#                                       |
|---------------|-------------------------------------------------|
| ICAL Standard | MP87199,MP87201,MP87202,MP87203,MP87204,MP87205 |
| ICV Standard  | MP87206                                         |
| CCV Standard  | MP87208                                         |
| ICSA Standard |                                                 |
| CRI Standard  | MP87210                                         |
| LCS Standard  |                                                 |
| Chk Standard  | MP87207,MP87209,MP87211,MP87214                 |

| Sr# | SampleId   | ClientID   | QcType   | Date           | Comment | Operator | Status |
|-----|------------|------------|----------|----------------|---------|----------|--------|
| 1   | S0         | S0         | CAL1     | 09/16/25 09:57 |         | mohan    | OK     |
| 2   | S0.2       | S0.2       | CAL2     | 09/16/25 10:00 |         | mohan    | OK     |
| 3   | S2.5       | S2.5       | CAL3     | 09/16/25 10:02 |         | mohan    | OK     |
| 4   | S5         | S5         | CAL4     | 09/16/25 10:04 |         | mohan    | OK     |
| 5   | S7.5       | S7.5       | CAL5     | 09/16/25 10:07 |         | mohan    | OK     |
| 6   | S10        | S10        | CAL6     | 09/16/25 10:12 |         | mohan    | OK     |
| 7   | ICV30      | ICV30      | ICV      | 09/16/25 10:15 |         | mohan    | OK     |
| 8   | ICB30      | ICB30      | ICB      | 09/16/25 10:17 |         | mohan    | OK     |
| 9   | CCV91      | CCV91      | CCV      | 09/16/25 10:19 |         | mohan    | OK     |
| 10  | CCB91      | CCB91      | CCB      | 09/16/25 10:22 |         | mohan    | OK     |
| 11  | CRA        | CRA        | CRDL     | 09/16/25 10:24 |         | mohan    | OK     |
| 12  | HighStd    | HighStd    | HIGH STD | 09/16/25 10:26 |         | mohan    | OK     |
| 13  | ChkStd     | ChkStd     | SAM      | 09/16/25 10:28 |         | mohan    | OK     |
| 14  | PB169703BL | PB169703BL | MB       | 09/16/25 10:31 |         | mohan    | OK     |
| 15  | PB169703BS | PB169703BS | LCS      | 09/16/25 10:33 |         | mohan    | OK     |
| 16  | Q3082-02   | WC-29      | SAM      | 09/16/25 10:36 |         | mohan    | OK     |
| 17  | Q3082-04   | WC-30      | SAM      | 09/16/25 10:39 |         | mohan    | OK     |
| 18  | Q3082-06   | WC-31      | SAM      | 09/16/25 10:41 |         | mohan    | OK     |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137199**

|                  |                                                 |              |                       |
|------------------|-------------------------------------------------|--------------|-----------------------|
| Review By        | mohan                                           | Review On    | 9/16/2025 3:04:24 PM  |
| Supervise By     | jaswal                                          | Supervise On | 9/16/2025 10:31:06 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                       |
| ICAL Standard    | MP87199,MP87201,MP87202,MP87203,MP87204,MP87205 |              |                       |
| ICV Standard     | MP87206                                         |              |                       |
| CCV Standard     | MP87208                                         |              |                       |
| ICSA Standard    |                                                 |              |                       |
| CRI Standard     | MP87210                                         |              |                       |
| LCS Standard     |                                                 |              |                       |
| Chk Standard     | MP87207,MP87209,MP87211,MP87214                 |              |                       |

|    |             |             |     |                |  |       |    |
|----|-------------|-------------|-----|----------------|--|-------|----|
| 19 | Q3082-08    | WC-32       | SAM | 09/16/25 10:43 |  | mohan | OK |
| 20 | Q3082-10    | WC-33       | SAM | 09/16/25 10:45 |  | mohan | OK |
| 21 | CCV92       | CCV92       | CCV | 09/16/25 10:48 |  | mohan | OK |
| 22 | CCB92       | CCB92       | CCB | 09/16/25 10:50 |  | mohan | OK |
| 23 | Q3091-07    | 12/9C-COMP  | SAM | 09/16/25 10:52 |  | mohan | OK |
| 24 | Q3091-14    | 12/9B-COMP  | SAM | 09/16/25 10:54 |  | mohan | OK |
| 25 | Q3092-02    | 291376      | SAM | 09/16/25 10:57 |  | mohan | OK |
| 26 | Q3092-04    | 279182      | SAM | 09/16/25 10:59 |  | mohan | OK |
| 27 | Q3092-04DUP | 279182DUP   | DUP | 09/16/25 11:01 |  | mohan | OK |
| 28 | Q3092-04MS  | 279182MS    | MS  | 09/16/25 11:04 |  | mohan | OK |
| 29 | Q3092-04MSD | 279182MSD   | MSD | 09/16/25 11:06 |  | mohan | OK |
| 30 | PB169704BL  | PB169704BL  | MB  | 09/16/25 11:08 |  | mohan | OK |
| 31 | PB169704BS  | PB169704BS  | LCS | 09/16/25 11:10 |  | mohan | OK |
| 32 | Q3088-02    | 90525-LQ    | SAM | 09/16/25 11:13 |  | mohan | OK |
| 33 | CCV93       | CCV93       | CCV | 09/16/25 11:18 |  | mohan | OK |
| 34 | CCB93       | CCB93       | CCB | 09/16/25 11:20 |  | mohan | OK |
| 35 | Q3088-02DUP | 90525-LQDUP | DUP | 09/16/25 11:22 |  | mohan | OK |
| 36 | Q3088-02MS  | 90525-LQMS  | MS  | 09/16/25 11:30 |  | mohan | OK |
| 37 | Q3088-02MSD | 90525-LQMSD | MSD | 09/16/25 11:32 |  | mohan | OK |
| 38 | Q3090-01    | 90825-B     | SAM | 09/16/25 11:34 |  | mohan | OK |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137199**

|                  |                                                 |              |                       |
|------------------|-------------------------------------------------|--------------|-----------------------|
| Review By        | mohan                                           | Review On    | 9/16/2025 3:04:24 PM  |
| Supervise By     | jaswal                                          | Supervise On | 9/16/2025 10:31:06 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                       |
| ICAL Standard    | MP87199,MP87201,MP87202,MP87203,MP87204,MP87205 |              |                       |
| ICV Standard     | MP87206                                         |              |                       |
| CCV Standard     | MP87208                                         |              |                       |
| ICSA Standard    |                                                 |              |                       |
| CRI Standard     | MP87210                                         |              |                       |
| LCS Standard     |                                                 |              |                       |
| Chk Standard     | MP87207,MP87209,MP87211,MP87214                 |              |                       |

|    |            |            |     |                |  |       |    |
|----|------------|------------|-----|----------------|--|-------|----|
| 39 | PB169677TB | PB169677TB | MB  | 09/16/25 11:37 |  | mohan | OK |
| 40 | Q3092-04L  | 279182L    | SD  | 09/16/25 11:39 |  | mohan | OK |
| 41 | Q3092-04A  | 279182A    | PS  | 09/16/25 11:41 |  | mohan | OK |
| 42 | Q3088-02L  | 90525-LQL  | SD  | 09/16/25 11:43 |  | mohan | OK |
| 43 | Q3088-02A  | 90525-LQA  | PS  | 09/16/25 11:46 |  | mohan | OK |
| 44 | CCV94      | CCV94      | CCV | 09/16/25 11:48 |  | mohan | OK |
| 45 | CCB94      | CCB94      | CCB | 09/16/25 11:50 |  | mohan | OK |