

Surrogate Summary

SDG No.: Q3103

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PO113719.D	PIBLK-PO113719.D	Tetrachloro-m-xyl	1	20	17.0	85		70 (60)	130 (140)
		Decachlorobiphen	1	20	17.3	86		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.0	85		70 (60)	130 (140)
		Decachlorobiphen	2	20	18.1	91		70 (60)	130 (140)
I.BLK-PO113788.D	PIBLK-PO113788.D	Tetrachloro-m-xyl	1	20	18.6	93		70 (60)	130 (140)
		Decachlorobiphen	1	20	21.8	109		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.9	90		70 (60)	130 (140)
		Decachlorobiphen	2	20	19.6	98		70 (60)	130 (140)
PB169713BL	PB169713BL	Tetrachloro-m-xyl	1	20	19.0	95		30 (30)	150 (173)
		Decachlorobiphen	1	20	22.5	112		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	18.3	92		30 (30)	150 (173)
		Decachlorobiphen	2	20	20.4	102		30 (10)	150 (173)
PB169713BS	PB169713BS	Tetrachloro-m-xyl	1	20	18.3	91		30 (30)	150 (173)
		Decachlorobiphen	1	20	23.1	116		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	17.3	87		30 (30)	150 (173)
		Decachlorobiphen	2	20	20.7	104		30 (10)	150 (173)
PB169713BSD	PB169713BSD	Tetrachloro-m-xyl	1	20	18.4	92		30 (30)	150 (173)
		Decachlorobiphen	1	20	23.0	115		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	17.7	88		30 (30)	150 (173)
		Decachlorobiphen	2	20	20.6	103		30 (10)	150 (173)
Q3103-01	TW-WTS-14	Tetrachloro-m-xyl	1	20	21.0	105		30 (30)	150 (173)
		Decachlorobiphen	1	20	19.0	95		30 (10)	150 (173)
		Tetrachloro-m-xyl	2	20	20.1	101		30 (30)	150 (173)
		Decachlorobiphen	2	20	17.3	87		30 (10)	150 (173)
I.BLK-PO113800.D	PIBLK-PO113800.D	Tetrachloro-m-xyl	1	20	18.4	92		70 (60)	130 (140)
		Decachlorobiphen	1	20	20.4	102		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.9	89		70 (60)	130 (140)
		Decachlorobiphen	2	20	19.5	98		70 (60)	130 (140)