

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137272**

|              |        |              |                       |
|--------------|--------|--------------|-----------------------|
| Review By    | mohan  | Review On    | 9/24/2025 10:22:42 AM |
| Supervise By | jaswal | Supervise On | 9/24/2025 10:54:45 AM |

| STD. NAME     | STD REF.#                                       |
|---------------|-------------------------------------------------|
| ICAL Standard | MP87286,MP87287,MP87288,MP87289,MP87290,MP87291 |
| ICV Standard  | MP87292                                         |
| CCV Standard  | MP87294                                         |
| ICSA Standard |                                                 |
| CRI Standard  | MP87296                                         |
| LCS Standard  |                                                 |
| Chk Standard  | MP87293,MP87295,MP87297,MP87301                 |

| Sr# | SampleId   | ClientID     | QcType   | Date           | Comment | Operator | Status |
|-----|------------|--------------|----------|----------------|---------|----------|--------|
| 1   | S0         | S0           | CAL1     | 09/23/25 08:40 |         | mohan    | OK     |
| 2   | S0.2       | S0.2         | CAL2     | 09/23/25 08:42 |         | mohan    | OK     |
| 3   | S2.5       | S2.5         | CAL3     | 09/23/25 08:44 |         | mohan    | OK     |
| 4   | S5         | S5           | CAL4     | 09/23/25 08:46 |         | mohan    | OK     |
| 5   | S7.5       | S7.5         | CAL5     | 09/23/25 08:52 |         | mohan    | OK     |
| 6   | S10        | S10          | CAL6     | 09/23/25 08:58 |         | mohan    | OK     |
| 7   | ICV40      | ICV40        | ICV      | 09/23/25 09:04 |         | mohan    | OK     |
| 8   | ICB40      | ICB40        | ICB      | 09/23/25 09:06 |         | mohan    | OK     |
| 9   | CCV20      | CCV20        | CCV      | 09/23/25 09:09 |         | mohan    | OK     |
| 10  | CCB20      | CCB20        | CCB      | 09/23/25 09:11 |         | mohan    | OK     |
| 11  | CRA        | CRA          | CRDL     | 09/23/25 09:16 |         | mohan    | OK     |
| 12  | HighStd    | HighStd      | HIGH STD | 09/23/25 09:18 |         | mohan    | OK     |
| 13  | ChkStd     | ChkStd       | SAM      | 09/23/25 09:21 |         | mohan    | OK     |
| 14  | PB169787BL | PB169787BL   | MB       | 09/23/25 09:23 |         | mohan    | OK     |
| 15  | PB169787BS | PB169787BS   | LCS      | 09/23/25 09:26 |         | mohan    | OK     |
| 16  | Q3145-01   | SP-A         | SAM      | 09/23/25 09:28 |         | mohan    | OK     |
| 17  | Q3145-03   | SP-B         | SAM      | 09/23/25 09:31 |         | mohan    | OK     |
| 18  | Q3152-01   | EO-02-091925 | SAM      | 09/23/25 09:33 |         | mohan    | OK     |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137272**

|                  |                                                 |              |                       |
|------------------|-------------------------------------------------|--------------|-----------------------|
| Review By        | mohan                                           | Review On    | 9/24/2025 10:22:42 AM |
| Supervise By     | jaswal                                          | Supervise On | 9/24/2025 10:54:45 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                       |
| ICAL Standard    | MP87286,MP87287,MP87288,MP87289,MP87290,MP87291 |              |                       |
| ICV Standard     | MP87292                                         |              |                       |
| CCV Standard     | MP87294                                         |              |                       |
| ICSA Standard    |                                                 |              |                       |
| CRI Standard     | MP87296                                         |              |                       |
| LCS Standard     |                                                 |              |                       |
| Chk Standard     | MP87293,MP87295,MP87297,MP87301                 |              |                       |

|    |             |                |     |                |  |       |    |
|----|-------------|----------------|-----|----------------|--|-------|----|
| 19 | Q3153-01    | OR-500-COMP-39 | SAM | 09/23/25 09:36 |  | mohan | OK |
| 20 | Q3153-05    | OR-500-COMP-40 | SAM | 09/23/25 09:38 |  | mohan | OK |
| 21 | CCV21       | CCV21          | CCV | 09/23/25 09:40 |  | mohan | OK |
| 22 | CCB21       | CCB21          | CCB | 09/23/25 09:43 |  | mohan | OK |
| 23 | Q3155-01    | COMP-1         | SAM | 09/23/25 09:45 |  | mohan | OK |
| 24 | Q3155-02    | COMP-2         | SAM | 09/23/25 09:48 |  | mohan | OK |
| 25 | Q3155-03    | COMP-3         | SAM | 09/23/25 09:50 |  | mohan | OK |
| 26 | Q3157-01    | 20-DEAD-1      | SAM | 09/23/25 09:52 |  | mohan | OK |
| 27 | Q3159-01    | VNJ-245        | SAM | 09/23/25 09:55 |  | mohan | OK |
| 28 | Q3159-03    | VNJ-255        | SAM | 09/23/25 10:02 |  | mohan | OK |
| 29 | Q3159-03DUP | VNJ-255DUP     | DUP | 09/23/25 10:05 |  | mohan | OK |
| 30 | Q3159-03MS  | VNJ-255MS      | MS  | 09/23/25 10:07 |  | mohan | OK |
| 31 | Q3159-03MSD | VNJ-255MSD     | MSD | 09/23/25 10:09 |  | mohan | OK |
| 32 | CCV22       | CCV22          | CCV | 09/23/25 10:12 |  | mohan | OK |
| 33 | CCB22       | CCB22          | CCB | 09/23/25 10:14 |  | mohan | OK |
| 34 | Q3159-03L   | VNJ-255L       | SD  | 09/23/25 10:17 |  | mohan | OK |
| 35 | Q3159-03A   | VNJ-255A       | PS  | 09/23/25 10:19 |  | mohan | OK |
| 36 | CCV23       | CCV23          | CCV | 09/23/25 10:22 |  | mohan | OK |
| 37 | CCB23       | CCB23          | CCB | 09/23/25 10:24 |  | mohan | OK |

**Instrument ID:** P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137295**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Janvi                                           | Review On    | 9/25/2025 3:54:42 PM |
| Supervise By     |                                                 | Supervise On |                      |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                      |
| ICV Standard     | MP87047,MP87038                                 |              |                      |
| CCV Standard     | MP87042                                         |              |                      |
| ICSA Standard    | MP87040,MP87048                                 |              |                      |
| CRI Standard     | MP87038                                         |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                      |

| Sr# | SampleId | ClientID       | QcType   | Date           | Comment | Operator | Status |
|-----|----------|----------------|----------|----------------|---------|----------|--------|
| 1   | S0       | S0             | CAL1     | 09/24/25 10:54 |         | Jaswal   | OK     |
| 2   | S1       | S1             | CAL2     | 09/24/25 10:58 |         | Jaswal   | OK     |
| 3   | S2       | S2             | CAL3     | 09/24/25 11:03 |         | Jaswal   | OK     |
| 4   | S3       | S3             | CAL4     | 09/24/25 11:07 |         | Jaswal   | OK     |
| 5   | S4       | S4             | CAL5     | 09/24/25 11:10 |         | Jaswal   | OK     |
| 6   | S5       | S5             | CAL6     | 09/24/25 11:15 |         | Jaswal   | OK     |
| 7   | ICV01    | ICV01          | ICV      | 09/24/25 11:34 |         | Jaswal   | OK     |
| 8   | LLICV01  | LLICV01        | LLICV    | 09/24/25 11:44 |         | Jaswal   | OK     |
| 9   | ICB01    | ICB01          | ICB      | 09/24/25 11:56 |         | Jaswal   | OK     |
| 10  | CRI01    | CRI01          | CRDL     | 09/24/25 12:00 |         | Jaswal   | OK     |
| 11  | ICSA01   | ICSA01         | ICSA     | 09/24/25 12:04 |         | Jaswal   | OK     |
| 12  | ICSAB01  | ICSAB01        | ICSAB    | 09/24/25 12:08 |         | Jaswal   | OK     |
| 13  | ICSADL   | ICSADL         | ICSA     | 09/24/25 12:13 |         | Jaswal   | OK     |
| 14  | ICSABDL  | ICSABDL        | ICSAB    | 09/24/25 12:17 |         | Jaswal   | OK     |
| 15  | CCV01    | CCV01          | CCV      | 09/24/25 12:21 |         | Jaswal   | OK     |
| 16  | CCB01    | CCB01          | CCB      | 09/24/25 12:26 |         | Jaswal   | OK     |
| 17  | LR-2     | LR-2           | HIGH STD | 09/24/25 12:35 |         | Jaswal   | OK     |
| 18  | Q3153-04 | OR-500-COMP-39 | SAM      | 09/24/25 12:45 |         | Jaswal   | OK     |

Instrument ID: P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137295**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Janvi                                           | Review On    | 9/25/2025 3:54:42 PM |
| Supervise By     |                                                 | Supervise On |                      |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                      |
| ICV Standard     | MP87047,MP87038                                 |              |                      |
| CCV Standard     | MP87042                                         |              |                      |
| ICSA Standard    | MP87040,MP87048                                 |              |                      |
| CRI Standard     | MP87038                                         |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                      |

|    |             |                  |          |                |                           |        |          |
|----|-------------|------------------|----------|----------------|---------------------------|--------|----------|
| 19 | LR-1        | LR-1             | HIGH STD | 09/24/25 12:49 |                           | Jaswal | OK       |
| 20 | LR-3        | LR-3             | HIGH STD | 09/24/25 12:53 |                           | Jaswal | OK       |
| 21 | Q3153-08    | OR-500-COMP-40   | SAM      | 09/24/25 12:58 |                           | Jaswal | OK       |
| 22 | Q3145-01    | SP-A             | SAM      | 09/24/25 13:03 |                           | Jaswal | OK       |
| 23 | Q3145-03    | SP-B             | SAM      | 09/24/25 13:07 |                           | Jaswal | OK       |
| 24 | Q3146-01    | MECHANIC-ST-PAIN | SAM      | 09/24/25 13:11 | Pb is high, need dilution | Jaswal | Dilution |
| 25 | Q3153-01    | OR-500-COMP-39   | SAM      | 09/24/25 13:16 |                           | Jaswal | OK       |
| 26 | Q3153-05    | OR-500-COMP-40   | SAM      | 09/24/25 13:20 |                           | Jaswal | OK       |
| 27 | CCV02       | CCV02            | CCV      | 09/24/25 13:40 |                           | Jaswal | OK       |
| 28 | CCB02       | CCB02            | CCB      | 09/24/25 13:44 |                           | Jaswal | OK       |
| 29 | Q3151-02    | 1402             | SAM      | 09/24/25 13:48 |                           | Jaswal | OK       |
| 30 | Q3146-01DL  | MECHANIC-ST-PAIN | SAM      | 09/24/25 14:01 | 5X For Pb                 | Jaswal | Confirms |
| 31 | Q3157-02    | 20-DEAD-1        | SAM      | 09/24/25 14:06 |                           | Jaswal | OK       |
| 32 | Q3157-02DUP | 20-DEAD-1DUP     | DUP      | 09/24/25 14:10 |                           | Jaswal | OK       |
| 33 | Q3146-01DL2 | MECHANIC-ST-PAIN | SAM      | 09/24/25 14:15 | Not use                   | Jaswal | Not Ok   |
| 34 | Q3146-01DL3 | MECHANIC-ST-PAIN | SAM      | 09/24/25 14:21 | Not use                   | Jaswal | Not Ok   |
| 35 | Q3157-02L   | 20-DEAD-1L       | SD       | 09/24/25 14:25 |                           | Jaswal | OK       |
| 36 | Q3157-02MS  | 20-DEAD-1MS      | MS       | 09/24/25 14:30 |                           | Jaswal | OK       |
| 37 | Q3157-02MSD | 20-DEAD-1MSD     | MSD      | 09/24/25 14:34 |                           | Jaswal | OK       |
| 38 | Q3157-02A   | 20-DEAD-1A       | PS       | 09/24/25 14:38 |                           | Jaswal | OK       |

Instrument ID: P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137295**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Janvi                                           | Review On    | 9/25/2025 3:54:42 PM |
| Supervise By     |                                                 | Supervise On |                      |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                      |
| ICV Standard     | MP87047,MP87038                                 |              |                      |
| CCV Standard     | MP87042                                         |              |                      |
| ICSA Standard    | MP87040,MP87048                                 |              |                      |
| CRI Standard     | MP87038                                         |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                      |

|    |             |                   |     |                |  |        |    |
|----|-------------|-------------------|-----|----------------|--|--------|----|
| 39 | CCV03       | CCV03             | CCV | 09/24/25 14:42 |  | Jaswal | OK |
| 40 | CCB03       | CCB03             | CCB | 09/24/25 14:46 |  | Jaswal | OK |
| 41 | Q3151-02DUP | 1402DUP           | DUP | 09/24/25 14:50 |  | Jaswal | OK |
| 42 | Q3151-02L   | 1402L             | SD  | 09/24/25 14:55 |  | Jaswal | OK |
| 43 | Q3151-02MS  | 1402MS            | MS  | 09/24/25 14:59 |  | Jaswal | OK |
| 44 | Q3151-02MSD | 1402MSD           | MSD | 09/24/25 15:03 |  | Jaswal | OK |
| 45 | Q3151-02A   | 1402A             | PS  | 09/24/25 15:07 |  | Jaswal | OK |
| 46 | Q3142-01    | RW8-SP100-2025091 | SAM | 09/24/25 15:11 |  | Jaswal | OK |
| 47 | Q3142-02    | RW8-SP303-2025091 | SAM | 09/24/25 15:15 |  | Jaswal | OK |
| 48 | PB169760BL  | PB169760BL        | MB  | 09/24/25 15:25 |  | Jaswal | OK |
| 49 | PB169760BS  | PB169760BS        | LCS | 09/24/25 15:33 |  | Jaswal | OK |
| 50 | Q3155-01    | COMP-1            | SAM | 09/24/25 15:37 |  | Jaswal | OK |
| 51 | CCV04       | CCV04             | CCV | 09/24/25 15:41 |  | Jaswal | OK |
| 52 | CCB04       | CCB04             | CCB | 09/24/25 15:46 |  | Jaswal | OK |
| 53 | Q3155-02    | COMP-2            | SAM | 09/24/25 15:50 |  | Jaswal | OK |
| 54 | Q3155-03    | COMP-3            | SAM | 09/24/25 15:54 |  | Jaswal | OK |
| 55 | Q3157-01    | 20-DEAD-1         | SAM | 09/24/25 15:58 |  | Jaswal | OK |
| 56 | Q3157-01DUP | 20-DEAD-1DUP      | DUP | 09/24/25 16:02 |  | Jaswal | OK |
| 57 | Q3157-01L   | 20-DEAD-1L        | SD  | 09/24/25 16:06 |  | Jaswal | OK |
| 58 | Q3157-01MS  | 20-DEAD-1MS       | MS  | 09/24/25 16:10 |  | Jaswal | OK |

Instrument ID: P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137295**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Janvi                                           | Review On    | 9/25/2025 3:54:42 PM |
| Supervise By     |                                                 | Supervise On |                      |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                      |
| ICV Standard     | MP87047,MP87038                                 |              |                      |
| CCV Standard     | MP87042                                         |              |                      |
| ICSA Standard    | MP87040,MP87048                                 |              |                      |
| CRI Standard     | MP87038                                         |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                      |

|    |             |              |     |                |  |        |    |
|----|-------------|--------------|-----|----------------|--|--------|----|
| 59 | Q3157-01MSD | 20-DEAD-1MSD | MSD | 09/24/25 16:14 |  | Jaswal | OK |
| 60 | Q3157-01A   | 20-DEAD-1A   | PS  | 09/24/25 16:18 |  | Jaswal | OK |
| 61 | PB169761BS  | PB169761BS   | LCS | 09/24/25 16:22 |  | Jaswal | OK |
| 62 | PB169761BL  | PB169761BL   | MB  | 09/24/25 16:26 |  | Jaswal | OK |
| 63 | CCV05       | CCV05        | CCV | 09/24/25 16:31 |  | Jaswal | OK |
| 64 | CCB05       | CCB05        | CCB | 09/24/25 16:35 |  | Jaswal | OK |
| 65 | Q3159-01    | VNJ-245      | SAM | 09/24/25 16:44 |  | Jaswal | OK |
| 66 | Q3159-03    | VNJ-255      | SAM | 09/24/25 16:48 |  | Jaswal | OK |
| 67 | Q3159-03DUP | VNJ-255DUP   | DUP | 09/24/25 16:52 |  | Jaswal | OK |
| 68 | Q3159-03L   | VNJ-255L     | SD  | 09/24/25 16:56 |  | Jaswal | OK |
| 69 | Q3159-03MS  | VNJ-255MS    | MS  | 09/24/25 17:00 |  | Jaswal | OK |
| 70 | Q3159-03MSD | VNJ-255MSD   | MSD | 09/24/25 17:04 |  | Jaswal | OK |
| 71 | Q3167-01    | RT4784       | SAM | 09/24/25 17:12 |  | Jaswal | OK |
| 72 | Q3159-03A   | VNJ-255A     | PS  | 09/24/25 17:25 |  | Jaswal | OK |
| 73 | PB169811BS  | PB169811BS   | LCS | 09/24/25 17:29 |  | Jaswal | OK |
| 74 | CCV06       | CCV06        | CCV | 09/24/25 17:33 |  | Jaswal | OK |
| 75 | CCB06       | CCB06        | CCB | 09/24/25 17:37 |  | Jaswal | OK |
| 76 | PB169811BL  | PB169811BL   | MB  | 09/24/25 17:41 |  | Jaswal | OK |
| 77 | Q3159-02    | VNJ-245      | SAM | 09/24/25 17:45 |  | Jaswal | OK |
| 78 | Q3159-04    | VNJ-255      | SAM | 09/24/25 17:50 |  | Jaswal | OK |

Instrument ID: P4

**Daily Analysis Runlog For Sequence/QC Batch ID # LB137295**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Janvi                                           | Review On    | 9/25/2025 3:54:42 PM |
| Supervise By     |                                                 | Supervise On |                      |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP87031,MP87038,MP87035,MP87034,MP87033,MP87032 |              |                      |
| ICV Standard     | MP87047,MP87038                                 |              |                      |
| CCV Standard     | MP87042                                         |              |                      |
| ICSA Standard    | MP87040,MP87048                                 |              |                      |
| CRI Standard     | MP87038                                         |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP87043,MP87044,MP87045,MP87046,MP87049         |              |                      |

|    |             |            |     |                |  |        |    |
|----|-------------|------------|-----|----------------|--|--------|----|
| 79 | Q3167-02    | RT4784     | SAM | 09/24/25 17:54 |  | Jaswal | OK |
| 80 | Q3167-02L   | RT4784L    | SD  | 09/24/25 18:03 |  | Jaswal | OK |
| 81 | Q3167-02DUP | RT4784DUP  | DUP | 09/24/25 18:07 |  | Jaswal | OK |
| 82 | Q3167-02MS  | RT4784MS   | MS  | 09/24/25 18:12 |  | Jaswal | OK |
| 83 | Q3167-02MSD | RT4784MSD  | MSD | 09/24/25 18:16 |  | Jaswal | OK |
| 84 | Q3167-02A   | RT4784A    | PS  | 09/24/25 18:20 |  | Jaswal | OK |
| 85 | PB169791TB  | PB169791TB | MB  | 09/24/25 18:24 |  | Jaswal | OK |
| 86 | CCV07       | CCV07      | CCV | 09/24/25 18:37 |  | Jaswal | OK |
| 87 | CCB07       | CCB07      | CCB | 09/24/25 18:42 |  | Jaswal | OK |
| 88 | PB169813BL  | PB169813BL | MB  | 09/24/25 18:46 |  | Jaswal | OK |
| 89 | PB169813BS  | PB169813BS | LCS | 09/24/25 18:59 |  | Jaswal | OK |
| 90 | CCV08       | CCV08      | CCV | 09/24/25 19:03 |  | Jaswal | OK |
| 91 | CCB08       | CCB08      | CCB | 09/24/25 19:07 |  | Jaswal | OK |