

Surrogate Summary

SDG No.: Q3181

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PO113719.D	PIBLK-PO113719.D	Tetrachloro-m-xyl	1	20	17.0	85		70 (60)	130 (140)
		Decachlorobiphen	1	20	17.3	86		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.0	85		70 (60)	130 (140)
		Decachlorobiphen	2	20	18.1	91		70 (60)	130 (140)
I.BLK-PO114081.D	PIBLK-PO114081.D	Tetrachloro-m-xyl	1	20	19.3	97		70 (60)	130 (140)
		Decachlorobiphen	1	20	20.9	104		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	18.1	90		70 (60)	130 (140)
		Decachlorobiphen	2	20	19.8	99		70 (60)	130 (140)
Q3181-02	WC-A3-01-C	Tetrachloro-m-xyl	1	20	16.9	84		30 (32)	150 (144)
		Decachlorobiphen	1	20	9.87	49		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	15.4	77		30 (32)	150 (144)
		Decachlorobiphen	2	20	9.50	48		30 (32)	150 (175)
Q3181-06	WC-A3-02-C	Tetrachloro-m-xyl	1	20	16.8	84		30 (32)	150 (144)
		Decachlorobiphen	1	20	13.1	65		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	16.8	84		30 (32)	150 (144)
		Decachlorobiphen	2	20	13.0	65		30 (32)	150 (175)
Q3181-10	WC-A3-03-C	Tetrachloro-m-xyl	1	20	16.5	82		30 (32)	150 (144)
		Decachlorobiphen	1	20	11.3	56		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	15.8	79		30 (32)	150 (144)
		Decachlorobiphen	2	20	11.1	55		30 (32)	150 (175)
I.BLK-PO114096.D	PIBLK-PO114096.D	Tetrachloro-m-xyl	1	20	18.5	93		70 (60)	130 (140)
		Decachlorobiphen	1	20	20.8	104		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.7	88		70 (60)	130 (140)
		Decachlorobiphen	2	20	19.6	98		70 (60)	130 (140)
I.BLK-PP075246.D	PIBLK-PP075246.D	Tetrachloro-m-xyl	1	20	19.9	100		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.4	97		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.9	90		70 (60)	130 (140)
		Decachlorobiphen	2	20	20.0	100		70 (60)	130 (140)
I.BLK-PP075356.D	PIBLK-PP075356.D	Tetrachloro-m-xyl	1	20	19.6	98		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.9	99		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	16.9	85		70 (60)	130 (140)
		Decachlorobiphen	2	20	20.9	104		70 (60)	130 (140)
PB169856BL	PB169856BL	Tetrachloro-m-xyl	1	20	22.9	114		30 (32)	150 (144)
		Decachlorobiphen	1	20	23.6	118		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	15.9	79		30 (32)	150 (144)
		Decachlorobiphen	2	20	24.9	125		30 (32)	150 (175)
PB169856BS	PB169856BS	Tetrachloro-m-xyl	1	20	15.1	76		30 (32)	150 (144)
		Decachlorobiphen	1	20	14.5	72		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	11.9	60		30 (32)	150 (144)
		Decachlorobiphen	2	20	15.2	76		30 (32)	150 (175)
Q3209-01MS	PL-02-092525MS	Tetrachloro-m-xyl	1	20	21.8	109		30 (32)	150 (144)

() = LABORATORY INHOUSE LIMIT

Surrogate Summary

SDG No.: Q3181

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
Q3209-01MS	PL-02-092525MS	Decachlorobiphen	1	20	16.4	82		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	12.5	63		30 (32)	150 (144)
		Decachlorobiphen	2	20	14.6	73		30 (32)	150 (175)
Q3209-01MSD	PL-02-092525MSD	Tetrachloro-m-xyl	1	20	22.6	113		30 (32)	150 (144)
		Decachlorobiphen	1	20	17.7	88		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	12.3	61		30 (32)	150 (144)
		Decachlorobiphen	2	20	16.3	82		30 (32)	150 (175)
I.BLK-PP075371.D	PIBLK-PP075371.D	Tetrachloro-m-xyl	1	20	19.3	96		70 (60)	130 (140)
		Decachlorobiphen	1	20	18.7	94		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	16.4	82		70 (60)	130 (140)
		Decachlorobiphen	2	20	18.1	91		70 (60)	130 (140)