

### Surrogate Summary

SDG No.: Q3181

Client: ENTACT

Analytical Method: 8151A

| Lab Sample ID    | Client ID        | Parameter | Column | Spike | Result | Recovery(%) | Qual | Limits(%) |           |
|------------------|------------------|-----------|--------|-------|--------|-------------|------|-----------|-----------|
|                  |                  |           |        |       |        |             |      | Low       | High      |
| I.BLK-PS031667.D | PIBLK-PS031667.D | 2,4-DCAA  | 1      | 500   | 445    | 89          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 525    | 105         |      | 70 (61)   | 130 (136) |
| I.BLK-PS031866.D | PIBLK-PS031866.D | 2,4-DCAA  | 1      | 500   | 524    | 105         |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 610    | 122         |      | 70 (61)   | 130 (136) |
| PB169875BL       | PB169875BL       | 2,4-DCAA  | 1      | 500   | 458    | 92          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 517    | 103         |      | 70 (61)   | 130 (136) |
| PB169875BS       | PB169875BS       | 2,4-DCAA  | 1      | 500   | 564    | 113         |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 589    | 118         |      | 70 (61)   | 130 (136) |
| PB169840TB       | PB169840TB       | 2,4-DCAA  | 1      | 500   | 415    | 83          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 493    | 99          |      | 70 (61)   | 130 (136) |
| Q3181-03         | WC-A3-01-C       | 2,4-DCAA  | 1      | 500   | 422    | 84          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 505    | 101         |      | 70 (61)   | 130 (136) |
| Q3181-03MS       | WC-A3-01-CMS     | 2,4-DCAA  | 1      | 500   | 470    | 94          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 581    | 116         |      | 70 (61)   | 130 (136) |
| Q3181-03MSD      | WC-A3-01-CMSD    | 2,4-DCAA  | 1      | 500   | 472    | 94          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 587    | 117         |      | 70 (61)   | 130 (136) |
| Q3181-07         | WC-A3-02-C       | 2,4-DCAA  | 1      | 500   | 441    | 88          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 495    | 99          |      | 70 (61)   | 130 (136) |
| Q3181-11         | WC-A3-03-C       | 2,4-DCAA  | 1      | 500   | 425    | 85          |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 471    | 94          |      | 70 (61)   | 130 (136) |
| I.BLK-PS031878.D | PIBLK-PS031878.D | 2,4-DCAA  | 1      | 500   | 520    | 104         |      | 70 (61)   | 130 (136) |
|                  |                  | 2,4-DCAA  | 2      | 500   | 583    | 117         |      | 70 (61)   | 130 (136) |