

Instrument ID: FID\_C

**Daily Analysis Runlog For Sequence/QC Batch ID # FC092325AL**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 9/23/2025 12:34:12 PM           |
| Supervise By                                                       | mohammad                                | Supervise On      | 9/25/2025 2:25:00 AM            |
| SubDirectory                                                       | FC092325AL                              | HP Acquire Method | HP Processing Method FC092325AL |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24768,PP24770,PP24771,PP24772,PP24773 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24771                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24769,PP24774                         |                   |                                 |

| Sr# | SampleId                    | Data File Name | Date-Time         | Operator | Status |
|-----|-----------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                       | FC069823.D     | 23 Sep 2025 09:03 | YP/AJ    | Ok     |
| 2   | I.BLK                       | FC069824.D     | 23 Sep 2025 09:47 | YP/AJ    | Ok     |
| 3   | 100 PPM ALIPHATIC HC STD1   | FC069825.D     | 23 Sep 2025 10:31 | YP/AJ    | Ok     |
| 4   | 50 PPM ALIPHATIC HC STD2    | FC069826.D     | 23 Sep 2025 11:16 | YP/AJ    | Ok     |
| 5   | 20 PPM ALIPHATIC HC STD3    | FC069827.D     | 23 Sep 2025 12:02 | YP/AJ    | Ok     |
| 6   | 10 PPM ALIPHATIC HC STD4    | FC069828.D     | 23 Sep 2025 12:47 | YP/AJ    | Ok     |
| 7   | 5 PPM ALIPHATIC HC STD5     | FC069829.D     | 23 Sep 2025 13:33 | YP/AJ    | Ok     |
| 8   | 20 PPM ALIPHATIC HC STD ICV | FC069830.D     | 23 Sep 2025 14:19 | YP/AJ    | Ok     |
| 9   | I.BLK                       | FC069831.D     | 23 Sep 2025 15:05 | YP/AJ    | Ok     |
| 10  | 20 PPM ALIPHATIC HC STD     | FC069832.D     | 23 Sep 2025 15:53 | YP/AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_C

**Daily Analysis Runlog For Sequence/QC Batch ID # FC092625AL**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 9/26/2025 10:12:39 AM           |
| Supervise By                                                       | mohammad                                | Supervise On      | 9/30/2025 1:28:00 AM            |
| SubDirectory                                                       | FC092625AL                              | HP Acquire Method | HP Processing Method FC092325AL |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24768,PP24770,PP24771,PP24772,PP24773 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24771                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24769,PP24774                         |                   |                                 |

| Sr# | SampleId                | Data File Name | Date-Time         | Operator | Status |
|-----|-------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                   | FC069851.D     | 26 Sep 2025 07:44 | YP/AJ    | Ok     |
| 2   | I.BLK                   | FC069852.D     | 26 Sep 2025 08:26 | YP/AJ    | Ok     |
| 3   | 20 PPM ALIPHATIC HC STD | FC069853.D     | 26 Sep 2025 09:09 | YP/AJ    | Ok     |
| 4   | PB169823BL              | FC069854.D     | 26 Sep 2025 10:37 | YP/AJ    | Ok     |
| 5   | PB169823BS              | FC069855.D     | 26 Sep 2025 13:39 | YP/AJ    | Ok     |
| 6   | PB169823BSD             | FC069856.D     | 26 Sep 2025 14:25 | YP/AJ    | Ok     |
| 7   | I.BLK                   | FC069857.D     | 26 Sep 2025 15:17 | YP/AJ    | Ok     |
| 8   | 20 PPM ALIPHATIC HC STD | FC069858.D     | 26 Sep 2025 16:04 | YP/AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD090425AR**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 9/4/2025 10:41:46 AM            |
| Supervise By                                                       | mohammad                                | Supervise On      | 9/8/2025 4:38:38 AM             |
| SubDirectory                                                       | FD090425AR                              | HP Acquire Method | HP Processing Method FD090425AR |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24500,PP24502,PP24503,PP24504,PP24505 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24503                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24501,PP24506                         |                   |                                 |

| Sr# | SampleId                   | Data File Name | Date-Time         | Operator | Status |
|-----|----------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                      | FD049716.D     | 04 Sep 2025 08:15 | YP/AJ    | Ok     |
| 2   | I.BLK                      | FD049717.D     | 04 Sep 2025 08:55 | YP/AJ    | Ok     |
| 3   | 100 PPM AROMATIC HC STD1   | FD049718.D     | 04 Sep 2025 09:36 | YP/AJ    | Ok,M   |
| 4   | 50 PPM AROMATIC HC STD2    | FD049719.D     | 04 Sep 2025 10:16 | YP/AJ    | Ok     |
| 5   | 20 PPM AROMATIC HC STD3    | FD049720.D     | 04 Sep 2025 10:56 | YP/AJ    | Ok     |
| 6   | 10 PPM AROMATIC HC STD4    | FD049721.D     | 04 Sep 2025 11:38 | YP/AJ    | Ok     |
| 7   | 5 PPM AROMATIC HC STD5     | FD049722.D     | 04 Sep 2025 12:19 | YP/AJ    | Ok     |
| 8   | 20 PPM AROMATIC HC STD ICV | FD049723.D     | 04 Sep 2025 13:00 | YP/AJ    | Ok     |
| 9   | I.BLK                      | FD049724.D     | 04 Sep 2025 14:24 | YP/AJ    | Ok     |
| 10  | 20 PPM AROMATIC HC STD     | FD049725.D     | 04 Sep 2025 15:08 | YP/AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD092625AR**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 9/26/2025 10:13:31 AM           |
| Supervise By                                                       | mohammad                                | Supervise On      | 9/30/2025 1:28:06 AM            |
| SubDirectory                                                       | FD092625AR                              | HP Acquire Method | HP Processing Method FD090425AR |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24500,PP24502,PP24503,PP24504,PP24505 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24503                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24501,PP24506                         |                   |                                 |

| Sr# | SampleId               | Data File Name | Date-Time         | Operator | Status |
|-----|------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                  | FD049763.D     | 26 Sep 2025 07:44 | YP/AJ    | Ok     |
| 2   | I.BLK                  | FD049764.D     | 26 Sep 2025 08:26 | YP/AJ    | Ok     |
| 3   | 20 PPM AROMATIC HC STD | FD049765.D     | 26 Sep 2025 09:09 | YP/AJ    | Ok     |
| 4   | Q3178-01               | FD049766.D     | 26 Sep 2025 09:52 | YP/AJ    | Ok     |
| 5   | Q3178-01D              | FD049767.D     | 26 Sep 2025 10:37 | YP/AJ    | Ok     |
| 6   | Q3178-02               | FD049768.D     | 26 Sep 2025 11:22 | YP/AJ    | Ok     |
| 7   | Q3178-03               | FD049769.D     | 26 Sep 2025 12:07 | YP/AJ    | Ok     |
| 8   | I.BLK                  | FD049770.D     | 26 Sep 2025 12:53 | YP/AJ    | Ok     |
| 9   | 20 PPM AROMATIC HC STD | FD049771.D     | 26 Sep 2025 13:39 | YP/AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD100125AR**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 10/1/2025 3:38:16 PM            |
| Supervise By                                                       | Sohil                                   | Supervise On      | 10/2/2025 1:04:55 PM            |
| SubDirectory                                                       | FD100125AR                              | HP Acquire Method | HP Processing Method FD090425AR |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24500,PP24502,PP24503,PP24504,PP24505 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24503                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24501,PP24506                         |                   |                                 |

| Sr# | SampleId               | Data File Name | Date-Time         | Operator | Status |
|-----|------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                  | FD049781.D     | 01 Oct 2025 08:17 | YP/AJ    | Ok     |
| 2   | I.BLK                  | FD049782.D     | 01 Oct 2025 08:58 | YP/AJ    | Ok     |
| 3   | 20 PPM AROMATIC HC STD | FD049783.D     | 01 Oct 2025 09:40 | YP/AJ    | Ok     |
| 4   | PB169823BS             | FD049784.D     | 01 Oct 2025 14:56 | YP/AJ    | Ok     |
| 5   | PB169823BL             | FD049785.D     | 01 Oct 2025 15:40 | YP/AJ    | Ok     |
| 6   | PB169823BSD            | FD049786.D     | 01 Oct 2025 16:25 | YP/AJ    | Ok     |
| 7   | I.BLK                  | FD049787.D     | 01 Oct 2025 17:10 | YP/AJ    | Ok     |
| 8   | 20 PPM AROMATIC HC STD | FD049788.D     | 01 Oct 2025 17:54 | YP/AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_F

**Daily Analysis Runlog For Sequence/QC Batch ID # FF092425AR**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 9/24/2025 12:57:21 PM           |
| Supervise By                                                       | mohammad                                | Supervise On      | 9/26/2025 4:25:38 AM            |
| SubDirectory                                                       | FF092425AR                              | HP Acquire Method | HP Processing Method FF092425AR |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24500,PP24502,PP24503,PP24504,PP24505 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24503                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24501,PP24506                         |                   |                                 |

| Sr# | SampleId                   | Data File Name | Date-Time         | Operator | Status |
|-----|----------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                      | FF016429.D     | 24 Sep 2025 13:00 | YP\AJ    | Ok     |
| 2   | I.BLK                      | FF016430.D     | 24 Sep 2025 13:30 | YP\AJ    | Ok     |
| 3   | 100 PPM AROMATIC HC STD1   | FF016431.D     | 24 Sep 2025 13:59 | YP\AJ    | Ok     |
| 4   | 50 PPM AROMATIC HC STD2    | FF016432.D     | 24 Sep 2025 14:29 | YP\AJ    | Ok     |
| 5   | 20 PPM AROMATIC HC STD3    | FF016433.D     | 24 Sep 2025 14:59 | YP\AJ    | Ok     |
| 6   | 10 PPM AROMATIC HC STD4    | FF016434.D     | 24 Sep 2025 15:28 | YP\AJ    | Ok     |
| 7   | 5 PPM AROMATIC HC STD5     | FF016435.D     | 24 Sep 2025 15:58 | YP\AJ    | Ok     |
| 8   | 20 PPM AROMATIC HC STD ICV | FF016436.D     | 24 Sep 2025 16:28 | YP\AJ    | Ok     |
| 9   | I.BLK                      | FF016437.D     | 24 Sep 2025 17:27 | YP\AJ    | Ok     |
| 10  | 20 PPM AROMATIC HC STD     | FF016438.D     | 24 Sep 2025 17:57 | YP\AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_F

**Daily Analysis Runlog For Sequence/QC Batch ID # FF092525AR**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 9/25/2025 2:53:23 PM            |
| Supervise By                                                       | Sohil                                   | Supervise On      | 9/26/2025 3:57:30 PM            |
| SubDirectory                                                       | FF092525AR                              | HP Acquire Method | HP Processing Method FF092425AR |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24500,PP24502,PP24503,PP24504,PP24505 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24503                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24501,PP24506                         |                   |                                 |

| Sr# | SampleId               | Data File Name | Date-Time         | Operator | Status |
|-----|------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                  | FF016449.D     | 25 Sep 2025 13:49 | YP\AJ    | Ok     |
| 2   | I.BLK                  | FF016450.D     | 25 Sep 2025 14:19 | YP\AJ    | Ok     |
| 3   | 20 PPM AROMATIC HC STD | FF016451.D     | 25 Sep 2025 14:49 | YP\AJ    | Ok     |
| 4   | Q3192-01               | FF016452.D     | 25 Sep 2025 15:50 | YP\AJ    | Ok     |
| 5   | Q3192-02               | FF016453.D     | 25 Sep 2025 16:20 | YP\AJ    | Ok     |
| 6   | Q3176-01               | FF016454.D     | 25 Sep 2025 16:52 | YP\AJ    | Ok     |
| 7   | Q3178-01               | FF016455.D     | 25 Sep 2025 17:22 | YP\AJ    | ReRun  |
| 8   | Q3178-01D              | FF016456.D     | 25 Sep 2025 17:52 | YP\AJ    | ReRun  |
| 9   | Q3178-01MS             | FF016457.D     | 25 Sep 2025 18:22 | YP\AJ    | Not Ok |
| 10  | Q3178-01MSD            | FF016458.D     | 25 Sep 2025 18:52 | YP\AJ    | Not Ok |
| 11  | PB169823BL             | FF016459.D     | 25 Sep 2025 19:22 | YP\AJ    | Not Ok |
| 12  | PB169823BS             | FF016460.D     | 25 Sep 2025 19:52 | YP\AJ    | Not Ok |
| 13  | I.BLK                  | FF016461.D     | 26 Sep 2025 07:41 | YP\AJ    | Ok     |
| 14  | 20 PPM AROMATIC HC STD | FF016462.D     | 26 Sep 2025 08:42 | YP\AJ    | Ok     |
| 15  | Q3176-02               | FF016463.D     | 26 Sep 2025 09:11 | YP\AJ    | Ok     |
| 16  | Q3176-03               | FF016464.D     | 26 Sep 2025 09:41 | YP\AJ    | Ok     |
| 17  | I.BLK                  | FF016465.D     | 26 Sep 2025 10:10 | YP\AJ    | Ok     |
| 18  | 20 PPM AROMATIC HC STD | FF016466.D     | 26 Sep 2025 10:40 | YP\AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_G

**Daily Analysis Runlog For Sequence/QC Batch ID # FG092425AL**

|                                                                    |                                         |                   |                                 |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                          | yogesh                                  | Review On         | 9/24/2025 1:00:45 PM            |
| Supervise By                                                       | mohammad                                | Supervise On      | 9/26/2025 4:25:56 AM            |
| SubDirectory                                                       | FG092425AL                              | HP Acquire Method | HP Processing Method FG092425AL |
| <b>STD. NAME</b>                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP24768,PP24770,PP24771,PP24772,PP24773 |                   |                                 |
| CCC<br>Internal Standard/PEM                                       | PP24771                                 |                   |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24769,PP24774                         |                   |                                 |

| Sr# | SampleId                    | Data File Name | Date-Time         | Operator | Status |
|-----|-----------------------------|----------------|-------------------|----------|--------|
| 1   | MECL2                       | FG016625.D     | 24 Sep 2025 13:00 | YP\AJ    | Ok     |
| 2   | I.BLK                       | FG016626.D     | 24 Sep 2025 13:30 | YP\AJ    | Ok     |
| 3   | 100 PPM ALIPHATIC HC STD1   | FG016627.D     | 24 Sep 2025 13:59 | YP\AJ    | Ok     |
| 4   | 50 PPM ALIPHATIC HC STD2    | FG016628.D     | 24 Sep 2025 14:29 | YP\AJ    | Ok     |
| 5   | 20 PPM ALIPHATIC HC STD3    | FG016629.D     | 24 Sep 2025 14:59 | YP\AJ    | Ok     |
| 6   | 10 PPM ALIPHATIC HC STD4    | FG016630.D     | 24 Sep 2025 15:28 | YP\AJ    | Ok     |
| 7   | 5 PPM ALIPHATIC HC STD5     | FG016631.D     | 24 Sep 2025 15:58 | YP\AJ    | Ok     |
| 8   | 20 PPM ALIPHATIC HC STD ICV | FG016632.D     | 24 Sep 2025 16:28 | YP\AJ    | Ok     |
| 9   | I.BLK                       | FG016633.D     | 24 Sep 2025 17:27 | YP\AJ    | Ok     |
| 10  | 20 PPM ALIPHATIC HC STD     | FG016634.D     | 24 Sep 2025 17:57 | YP\AJ    | Ok     |

M : Manual Integration



Instrument ID: FID\_G

**Daily Analysis Runlog For Sequence/QC Batch ID # FG092525AL**

|                                                                                                    |                                         |                   |                                 |
|----------------------------------------------------------------------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                                                                                          | yogesh                                  | Review On         | 9/25/2025 2:52:00 PM            |
| Supervise By                                                                                       | Sohil                                   | Supervise On      | 10/2/2025 1:07:29 PM            |
| SubDirectory                                                                                       | FG092525AL                              | HP Acquire Method | HP Processing Method FG092425AL |
| <b>STD. NAME</b>                                                                                   | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds                                                            | PP24768,PP24770,PP24771,PP24772,PP24773 |                   |                                 |
| CCC<br>Internal Standard/PEM<br>ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP24771<br><br>PP24769,PP24774          |                   |                                 |

| Sr# | SampleId                | Data File Name | Date-Time         | Operator | Status   |
|-----|-------------------------|----------------|-------------------|----------|----------|
| 1   | MECL2                   | FG016645.D     | 25 Sep 2025 14:19 | YP\AJ    | Ok       |
| 2   | I.BLK                   | FG016646.D     | 25 Sep 2025 14:49 | YP\AJ    | Ok       |
| 3   | 20 PPM ALIPHATIC HCSTD  | FG016647.D     | 25 Sep 2025 15:19 | YP\AJ    | Ok       |
| 4   | Q3192-01                | FG016648.D     | 25 Sep 2025 15:50 | YP\AJ    | Dilution |
| 5   | Q3192-02                | FG016649.D     | 25 Sep 2025 16:20 | YP\AJ    | Dilution |
| 6   | Q3192-01DL              | FG016650.D     | 25 Sep 2025 16:52 | YP\AJ    | Ok       |
| 7   | Q3192-02DL              | FG016651.D     | 25 Sep 2025 17:22 | YP\AJ    | Not Ok   |
| 8   | Q3192-02DL              | FG016652.D     | 25 Sep 2025 17:52 | YP\AJ    | Ok       |
| 9   | Q3176-01                | FG016653.D     | 25 Sep 2025 18:22 | YP\AJ    | Ok       |
| 10  | Q3176-02                | FG016654.D     | 25 Sep 2025 18:52 | YP\AJ    | Ok       |
| 11  | Q3176-03                | FG016655.D     | 25 Sep 2025 19:22 | YP\AJ    | Ok       |
| 12  | Q3178-01                | FG016656.D     | 25 Sep 2025 19:52 | YP\AJ    | Ok       |
| 13  | I.BLK                   | FG016657.D     | 26 Sep 2025 07:41 | YP\AJ    | Ok       |
| 14  | 20 PPM ALIPHATIC HC STD | FG016658.D     | 26 Sep 2025 08:42 | YP\AJ    | Ok,M     |
| 15  | Q3178-01D               | FG016659.D     | 26 Sep 2025 09:11 | YP\AJ    | Ok       |
| 16  | Q3178-02                | FG016660.D     | 26 Sep 2025 09:41 | YP\AJ    | Ok,M     |
| 17  | Q3178-03                | FG016661.D     | 26 Sep 2025 10:10 | YP\AJ    | Ok,M     |
| 18  | I.BLK                   | FG016662.D     | 26 Sep 2025 10:40 | YP\AJ    | Ok       |
| 19  | 20 PPM ALIPHATIC HC STD | FG016663.D     | 26 Sep 2025 11:10 | YP\AJ    | Ok,M     |

M : Manual Integration