

Surrogate Summary

SDG No.: **Q3238**

Client: **Yannuzzi Group, Inc.**

Analytical Method: **8082A**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PP075246.D	PIBLK-PP075246.D	Tetrachloro-m-xyl	1	20	19.9	100		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.4	97		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.9	90		70 (60)	130 (140)
		Decachlorobiphen	2	20	20.0	100		70 (60)	130 (140)
I.BLK-PP075438.D	PIBLK-PP075438.D	Tetrachloro-m-xyl	1	20	20.5	102		70 (60)	130 (140)
		Decachlorobiphen	1	20	22.4	112		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	17.6	88		70 (60)	130 (140)
		Decachlorobiphen	2	20	24.2	121		70 (60)	130 (140)
PB169907BL	PB169907BL	Tetrachloro-m-xyl	1	20	22.2	111		30 (32)	150 (144)
		Decachlorobiphen	1	20	25.5	127		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	21.6	108		30 (32)	150 (144)
		Decachlorobiphen	2	20	26.6	133		30 (32)	150 (175)
PB169907BS	PB169907BS	Tetrachloro-m-xyl	1	20	15.8	79		30 (32)	150 (144)
		Decachlorobiphen	1	20	17.1	85		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	8.57	43		30 (32)	150 (144)
		Decachlorobiphen	2	20	17.1	86		30 (32)	150 (175)
Q3238-01	RCA	Tetrachloro-m-xyl	1	20	20.6	103		30 (32)	150 (144)
		Decachlorobiphen	1	20	21.9	110		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	20.5	102		30 (32)	150 (144)
		Decachlorobiphen	2	20	23.2	116		30 (32)	150 (175)
Q3240-01MS	OK-02-092925MS	Tetrachloro-m-xyl	1	20	21.7	108		30 (32)	150 (144)
		Decachlorobiphen	1	20	20.8	104		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	20.8	104		30 (32)	150 (144)
		Decachlorobiphen	2	20	19.0	95		30 (32)	150 (175)
Q3240-01MSD	OK-02-092925MSD	Tetrachloro-m-xyl	1	20	22.1	111		30 (32)	150 (144)
		Decachlorobiphen	1	20	20.9	105		30 (32)	150 (175)
		Tetrachloro-m-xyl	2	20	23.9	120		30 (32)	150 (144)
		Decachlorobiphen	2	20	18.8	94		30 (32)	150 (175)
I.BLK-PP075453.D	PIBLK-PP075453.D	Tetrachloro-m-xyl	1	20	17.7	88		70 (60)	130 (140)
		Decachlorobiphen	1	20	19.6	98		70 (60)	130 (140)
		Tetrachloro-m-xyl	2	20	16.3	81		70 (60)	130 (140)
		Decachlorobiphen	2	20	18.8	94		70 (60)	130 (140)