

Surrogate Summary

SDG No.: Q3257

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS031667.D	PIBLK-PS031667.D	2,4-DCAA	1	500	445	89		61	136
		2,4-DCAA	2	500	525	105		61	136
I.BLK-PS031912.D	PIBLK-PS031912.D	2,4-DCAA	1	500	468	94		61	136
		2,4-DCAA	2	500	540	108		61	136
PB169931BL	PB169931BL	2,4-DCAA	1	500	426	85		10	141
		2,4-DCAA	2	500	499	100		10	141
PB169931BS	PB169931BS	2,4-DCAA	1	500	561	112		10	141
		2,4-DCAA	2	500	531	106		10	141
I.BLK-PS031917.D	PIBLK-PS031917.D	2,4-DCAA	1	500	461	92		61	136
		2,4-DCAA	2	500	523	105		61	136
I.BLK-PS031920.D	PIBLK-PS031920.D	2,4-DCAA	1	500	440	88		61	136
		2,4-DCAA	2	500	512	102		61	136
Q3257-01	ENV1-6FT	2,4-DCAA	1	500	214	43		10	141
		2,4-DCAA	2	500	238	48		10	141
Q3257-02	ENV2-6FT	2,4-DCAA	1	500	420	84		10	141
		2,4-DCAA	2	500	365	73		10	141
I.BLK-PS031928.D	PIBLK-PS031928.D	2,4-DCAA	1	500	464	93		61	136
		2,4-DCAA	2	500	537	107		61	136
I.BLK-PS031940.D	PIBLK-PS031940.D	2,4-DCAA	1	500	460	92		61	136
		2,4-DCAA	2	500	530	106		61	136
Q3257-02MS	ENV2-6FTMS	2,4-DCAA	1	500	445	89		10	141
		2,4-DCAA	2	500	432	86		10	141
Q3257-02MSD	ENV2-6FTMSD	2,4-DCAA	1	500	443	89		10	141
		2,4-DCAA	2	500	427	85		10	141
I.BLK-PS031944.D	PIBLK-PS031944.D	2,4-DCAA	1	500	478	96		61	136
		2,4-DCAA	2	500	544	109		61	136