

Surrogate Summary

SDG No.: Q3266

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Recovery(%)	Qual	Limits(%)	
								Low	High
I.BLK-PS031667.D	PIBLK-PS031667.D	2,4-DCAA	1	500	445	89		61	136
		2,4-DCAA	2	500	525	105		61	136
I.BLK-PS031954.D	PIBLK-PS031954.D	2,4-DCAA	1	500	448	90		61	136
		2,4-DCAA	2	500	537	107		61	136
PB169991BL	PB169991BL	2,4-DCAA	1	500	444	89		10	141
		2,4-DCAA	2	500	486	97		10	141
PB169991BS	PB169991BS	2,4-DCAA	1	500	619	124		10	141
		2,4-DCAA	2	500	558	112		10	141
Q3266-01	ENV-3-6FT	2,4-DCAA	1	500	292	58		10	141
		2,4-DCAA	2	500	324	65		10	141
Q3266-01MS	ENV-3-6FTMS	2,4-DCAA	1	500	453	91		10	141
		2,4-DCAA	2	500	495	99		10	141
Q3266-01MSD	ENV-3-6FTMSD	2,4-DCAA	1	500	532	106		10	141
		2,4-DCAA	2	500	502	100		10	141
Q3266-02	ENV-4-6FT	2,4-DCAA	1	500	386	77		10	141
		2,4-DCAA	2	500	167	33		10	141
I.BLK-PS031967.D	PIBLK-PS031967.D	2,4-DCAA	1	500	488	98		61	136
		2,4-DCAA	2	500	547	109		61	136